



PERIPHERAL INTERVENTIONS

2024 PRODUCT CATALOG



INDEX

Customer Order Information 3

C-Code Cross Reference Guide Notes 6

Alphabetical Index by Product Category 8

Product Catalog 15



Customer Ordering Information

Terms and Conditions:

For Customer convenience, the Boston Scientific ordering system can accept the UPN (Universal Product Number) e.g., H74939053052070 found on the product labeling as well as Boston Scientific’s historical order numbers e.g., 39053-05207. Please use either one of these numbers when placing an order via phone, fax, Electronic Data Interchange (EDI) or online. Customer Service Representatives will be happy to help you with any questions. We are eager to help make the ordering process as simple and efficient as possible.

WARNING REGARDING SINGLE USE DEVICES: Products marked as “single-use” are for single-use only. Do not reuse, reprocess or re-sterilize single-use products. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Boston Scientific assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties, express or implied, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose, with respect to any such instrument.

NOTE

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device.

For more information, please contact your Boston Scientific sales representative.

Products shown are subject to modification without prior notice

2024 Peripheral Interventions Product Catalog



Peripheral Intervention Products

Place your order directly with Boston Scientific Customer Service by mail, telephone, fax, Electronic Data Interchange (EDI) or online.

Boston Scientific Corporation

Peripheral Intervention

500 Commander Shea Blvd.
Quincy, MA 02171-1518

To Place an Order

1.888.272.1001 Customer Service Phone
1.888.272.9444 Customer Service Fax
bscemailorders@bsci.com

Cap and Risk Fax Orders

Fax: 1.888.809.7120
prior

Billing Inquiries

Phone: 1.888.272.9442

Hours

Customer Service is available from 8:30 am EST to 8:00 pm EST weekdays excluding national holidays. All orders need to be placed by 7:00 pm. Our toll-free fax number is available at all times. In addition, Technical and Emergency Service is available 24 hours a day, 7 days a week.

Terms

- Terms are Net 30 Days from date of invoice.
- For peripheral products: All Shipments are F.O.B. Quincy, Massachusetts.
- Prices are subject to change without notice.

Standing Orders

Boston Scientific accepts standing orders on all products. This convenience provides you with scheduled shipments using a single purchase order. Contact the Customer Service DataBase Department at 1.888.325.1060 for details.

Product Exchanges

Product exchanges will no longer be processed with a no charge bill. Instead, customer service will process a return for credit and billing for the replacement product.

Returned Goods

Authorization to return merchandise is required for all returns and may be withheld at the sole discretion of Boston Scientific. Contact the Customer Service Returns Department for a Return Goods Authorization (RGA) number.

Customer Service Returns Department: **1-800-862-1284**

The request for authorization must include:

- Reason for the return
- Order number, lot number, quantity
- Purchase order number, Boston Scientific invoice number, date of invoice

Boston Scientific Corporation
Customer Service Returns Department
500 Commander Shea Blvd.
Quincy, MA 02171

Returns Requirements and Policies

- A copy of the packing slip or invoice must accompany all merchandise returned to Boston Scientific.
- The Returned Goods Authorization (RGA) number must be clearly marked on the shipping container. Merchandise returned without an RGA number will not be accepted or considered returned for purposes of the warranty.
- Timeframe for returns is 90 days from date of shipment and product must have 6 months of shelf life remaining to be authorized for returns unless the product shipped out with less than 6 months. Restocking fees may apply.
- Products manufactured as Specials or Custom products are not returnable.
- Material that is opened, damaged, with customer labels applied on packaging, with pen marks, expired or tampered with is not returnable for credit.
- Capital Equipment and Accessories processed in error by Boston Scientific will receive full credit if returned within 30 days of invoice date in saleable condition. Cannot be returned once customer has taken title to the equipment.
- Even Exchanges of customer owned inventory require one of the following approvals: Regional Manager, Divisional Manager approval or Sales Rep cost center.
- Customer absorbs shipping costs, unless BSC or transportation carrier error.
- RGA requests stay open for 45 days. If product is not received in that timeframe a new request must be initiated.

SPECIFICATION CHANGES: NOTICE

As we are continually improving our products, the equipment supplied may not agree in all details with the description and specification shown in the catalog and/or price list. All prices and specifications are subject to change without notice. Federal (USA) law restricts these devices to sale by or on the order of a physician. Indications, contraindications, and instructions for use, as well as product warranty information, can be found in the product labeling supplied with each device.

Boston Scientific Corporation makes no warranties, express or implied, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose, except as specifically set forth in the product labeling supplied with each device.

TRADEMARK NOTICE

CPT® is a registered trademark of the American Medical Association

SUPER SHEATH IS A TRADEMARK OF MEDIKIT CO., LTD.

C-Code Cross-Reference Guide to Boston Scientific Peripheral Interventions Products

C-Codes:

C-Codes should be included on all Medicare claims where procedures are billed under the Outpatient Prospective Payment System (OPPS) and a C-Code is assigned. They are used to track the types of devices used during applicable outpatient procedures and to ensure that both procedures and devices are appropriately billed. Separate device reimbursement is not usually driven by C-Code use, but CMS uses the cost data supplied on the claim to set appropriate procedure payment. Some devices can qualify for additional hospital reimbursement through Medicare transitional pass-through payment.

Key Points of the Guide:

- The Guide cross-references Boston Scientific products at the UPN/order number level to CMS associated Category C-Codes.
- Medicare does not have a C-Code category for all BSC PI products. The letters 'N/A' will appear to denote products without a Medicare assigned C-Code.

Changes to Centers for Medicare and Medicaid Services C-Code Policy

Effective January 1, 2005, Medicare identified a number of device-dependent Ambulatory Payment Classifications (APC) that would require C-Code reporting for outpatient procedures to pass C-Code edits. This requirement was enforced beginning April 1, 2005 when CMS began to review claims at the CPT®/HCPCS code level. In the CY 2015 final rule, CMS finalized a policy and implemented claims processing edits that require any of the device codes used in the previous device-to-procedure edits to be present on the claim whenever a procedure code assigned to any of the former device dependent APCs are reported on the claim (79 FR 66795). Despite the expiration of pass-through payment status for device category codes, hospitals are still required to report the device category C-codes on claims when such devices are used in conjunction with procedures billed and paid under the OPPS. For additional information: <https://www.cms.gov/Medicare/Medicare-Fee-for-Service-Payment/HospitalOutpatientPPS/Downloads/Compleat-list-DeviceCats-OPPS.pdf>

IMPORTANT INFORMATION

Health economic and reimbursement information provided by Boston Scientific Corporation is gathered from third-party sources and is subject to change without notice as a result of complex and frequently changing laws, regulations, rules, and policies. This information is presented for illustrative purposes only and does not constitute reimbursement or legal advice. Boston Scientific encourages providers to submit accurate and appropriate claims for services. **It is always the provider's sole responsibility to determine medical necessity, the proper site for delivery of any services, and to submit appropriate codes, charges, and modifiers for services rendered.** It is also always the provider's responsibility to understand and comply with Medicare national coverage determinations (NCD), Medicare local coverage determinations (LCD), and any other coverage requirements established by relevant payers which can be updated frequently. Boston Scientific recommends that you consult with your payers, reimbursement specialists, and/or legal counsel regarding coding, coverage, and reimbursement matters.

2024 Peripheral Interventions Product Catalog



Boston Scientific does not promote the use of its products outside its FDA-approved label. Payer policies will vary and should be verified before treatment for limitations on diagnosis, coding, or site of service requirements. All trademarks are the property of their respective owners.

This coding information may include codes for procedures for which Boston Scientific currently offers no cleared or approved products. In those instances, such codes have been included solely in the interest of providing users with comprehensive coding information and are not intended to promote the use of any Boston Scientific products for which they are not cleared or approved. The Health Care Provider (HCP) is solely responsible for selecting the site of service and treatment modalities appropriate for the patient based on medically appropriate needs of that patient and the independent medical judgment of the HCP.

CPT ® Copyright 2020 American Medical Association. All rights reserved. CPT is a registered trademark of the American Medical Association. Applicable FARS/DFARS Restrictions Apply to Government Use. Fee schedules, relative value units, conversion factors and/or related components are not assigned by the AMA, are not part of CPT, and the AMA is not recommending their use. The AMA does not directly or indirectly practice medicine or dispense medical services. The AMA assumes no liability for data contained or not contained herein. All trademarks are the property of their respective owners.

In response to feedback, we have created a C-Code finder that crosswalks Boston Scientific products and their corresponding C-Codes. The finder can be found at this website: <http://www.bostonscientific.com/en-US/reimbursement/ccode-finder.html>

(Current as of September 2024. Information subject to change without notice)

		Page
ABLATION	Radiofrequency Ablation	
	LeVeen™ Needle Electrode	15
	LeVeen SuperSlim™ Needle Electrode	15
	LeVeen CoAccess™ Electrode System	15
	Soloist™ Electrode	15
	Cryoablation Needles	
	Cryoablation CX Needles	15
	Cryoablation Classic Needles	16
	Cryoablation MRI Needles	16
	Cryoablation Systems	
	Visual-ICE™ Cryoablation System	16
	Visual-ICE™ MRI Cryoablation System	16
ICEfx™ Cryoablation System	16	
Cryoablation Accessories	16	
ACCESS DEVICES AND ACCESSORIES	Introducers & Accessories	
	AccuStick™ II with RO Marker Introducer Systems	17
	FloSwitch™ HP Device	17
	Gateway™ Plus Y-Adapter	17
	Percufix™ Catheter Cuff Kit	17
	Connecting Tube for Drainage Bag	17
	ESP™ Radiation Reduction Examination Gloves	17
	Inflation Device	
	Encore™ 26 Inflation Device	18
	Encore™ 40 Inflation Device	18
	Introducer Sheaths	
	Super Sheath™ Introducer	18
Catapult™ Guide Sheath	20	
ANGIOGRAPHY	Tip Shape Guide	21
	Imager™ II Select Catheters	
	Imager™ II - 5 F Select Catheters	22

	Page
ANGIOPLASTY BALLOON DILATATION CATHETERS	0.014" Low Profile
	Coyote™ Over-the-Wire Balloon Dilatation Catheter 23
	Coyote™ Monorail Balloon Dilatation Catheter 24
	Coyote™ ES Over-the-Wire Balloon Dilatation Catheter 26
	Coyote™ ES Monorail™ Balloon Dilatation Catheter 26
	0.018" Low Profile
	Sterling™ Over-the-Wire Balloon Dilatation Catheter 27
	Sterling™ Monorail™ Balloon Dilatation Catheter 30
	Sterling™ SL Over-the-Wire Balloon Dilatation Catheter 32
	Sterling™ SL Monorail™ Balloon Dilatation Catheter 32
	0.035" High Pressure
	Mustang™ Balloon Dilatation Catheter 33
	Charger™ 0.035" Balloon Dilatation Catheter 38
	0.035" Ultra-High Pressure
	Athletis™ Ultra-High Pressure Balloon 40
	0.035" Large Vessel
	XXL™ Balloon Dilatation Catheter 42
	Occlusion Balloons
	Equalizer™ Occlusion Balloon Catheter 42
	Specialty Balloons
	Peripheral Cutting Balloon™ Microsurgical Dilatation Device 43
ATHERECTOMY	Peripheral Rotablator™ Rotational Atherectomy System 44
	Peripheral ROTAPRO™ Rotational Atherectomy System 44
	Jetstream™ Atherectomy System 45
DRAINAGE	Abscess Drainage
	All-Purpose Drainage Catheters (Flexima™ Material) 46
	vanSonnenberg™ Sump Catheters 46
	All-Purpose Drainage Catheters (Percuflex™ Material) 47
	vanSonnenberg™ Chest Drainage Catheters 47

		Page
DRAINAGE (Continued)	Biliary Drainage	
	Biliary Drainage Catheters	47
	Biliary Drainage Catheters with RO Marker	47
	Urinary Drainage	
	Flexima™ Material (Nephrostomy Catheters and Kits)	48
	Flexima™ Material (Ureteral Stents and Kits)	48
	Percuflex® Material (Catheters)	49
	Percuflex® Material (Stents)	49
	Percuflex® Material (Kits)	50
DRUG-ELUTING	Stents	
	Eluvia™ Drug-Eluting Vascular Stent System	51
	Balloons	
	Ranger™ Drug-Coated Balloon	51
EMBOLIC BEADS & PARTICLES	Microspheres	
	LC Bead M1™	53
	LC Bead LUMI™	53
	LC Bead™	53
	Particles	
	Contour™ - PVA Particles	53
GUIDE CATHETERS	Mach 1™ Peripheral Guide Catheter	54
	Rubicon™ Support Catheter	55
	Rubicon Control Support Catheter	55

GUIDEWIRES	Specialty Guidewires	
	Amplatz Super Stiff™ Guidewires	57
	Meier Guidewire	57
	Platinum Plus™ Guidewires	57
	SAVION™ 0.014 Guidewires	58
	Thruway™ 0.014" Guidewires	58
	Thruway™ 0.018" Guidewires	58
	V-14™ ControlWire™ Guidewires	59
	V-18™ ControlWire™ Guidewires	59
	Victory™ 14 Guidewire	59
	Victory™ 18 Guidewire	59
	ZIPwire® Guide Wire - J-Tip	60
	ZIPwire® Guide Wire - Long Taper Standard Shaft	60
	ZIPwire® Guide Wire - Long Taper Stiff Shaft	60
	ZIPwire® Guide Wire - Small Vessel	60
	Standard Guidewires	
	Starter™ Guidewires (PTFE Coated)	60
	Starter™ Guidewires (PTFE/Heparin Coated)	62
	ZIPwire™ Hydrophilic Guidewires	
	ZIPwire™ Guide Wire - Standard Shaft	62
	ZIPwire™ Guide Wire - Stiff Shaft	63
	<i>*For Radiology Guidewires, please see Embolization section.</i>	
IO & EMBOLIZATION	Guidewires	
	Fathom™ Peripheral Embolization Guidewires	64
	Transend™ Peripheral Embolization Guidewires	64
	Microcatheters	
	TruSelect™ Microcatheter	64
	Direxion™ Torqueable Microcatheters	65
	Direxion™ Fathom™ - 16 System Pre-Loaded Torqueable Microcatheters	65
	Direxion™ Transend™ - 14 System Pre-Loaded Torqueable Microcatheters	65
	Direxion™ HI-FLO™ Torqueable Microcatheters	65
	Direxion™ HI-FLO Fathom™ - 16 System Pre-Loaded Torqueable Microcatheters	65
	Direxion™ HI-FLO Transend™ - 18 System Pre-Loaded Torqueable Microcatheters	66
	Renegade™ Fiber Braided Microcatheters with Hydro Pass™ Hydrophilic Coating	66

		Page
IO & EMBOLIZATION (Continued)	Microcatheters	
	Renegade™ STC-18 Microcatheters	66
	Renegade™ HI-FLO™ Microcatheters	66
	Microcatheter Kits/Systems	
	Renegade® HI-FLO Microcatheter Kits with Transend™-18 Guidewires	67
	Renegade® HI-FLO™ Fathom® Pre-Loaded System with Fathom®-16 Guidewires	67
	Occlusion Balloons	
	Standard Occlusion Balloon Catheters	67
	Berenstein Occlusion Balloon Catheters	67
	Conformable Embolics	
	Obsidio™	67
	Fully Retractable Coils	
	EMBOLD™ Fibered Detachable Coil System	68
	EMBOLD™ Soft Detachable Coil System	68
	EMBOLD™ Packing Detachable Coil System	69
	Detachable Coils	
	Interlock™ - 35 Fibered IDC Occlusion System	69
	Interlock™ - 18 Fibered IDC Occlusion System	70
	IDC™ Interlocking Detachable Coil	70
	Pushable Coils	
	VortX™ Diamond - 18	70
	VortX™ - 18 Coil	71
	Straight - 18 Coil	71
	Figure 8 - 18 Coil	71
	Multi-Loop -18 Coil	71
	Complex Helical - 18 Coil	71
	Coil Pusher -16 Coil	71
	2D Helical - 35 Coil	71
	VortX™ - 35 Coil	72
RADIOEMBOLIZATION (Y-90)	Y-90 Glass Microspheres TheraSphere™	73
DOSIMETRY	Y-90 Personalized Dosimetry Simplicit90Y™ Personalized Dosimetry	74

	Page
STENTS (NONVASCULAR)	Metallic Stents Express™ LD Biliary Over-the-Wire Premounted Stent System 75 Express™ SD Renal Monorail™ Premounted Stent System 75 WALLSTENT™ RP Endoprosthesis 76 WALLSTENT™ Endoprosthesis with Unistep™ Plus Delivery System 76 WallFlex™ Biliary Transhepatic Stent System 77
STENTS (VASCULAR)	Metallic Stents Carotid WALLSTENT™ Monorail™ Endoprosthesis 78 Innova™ Vascular Self-Expanding Stent System 78 Express™ LD Iliac Over-the-Wire Premounted Stent System 79 Express™ SD Renal Monorail™ Premounted Stent System 80 WALLSTENT™ RP Endoprosthesis 80 Epic™ Self-Expanding Vascular Stent System 81
STENTS (VENOUS)	VENOUS WALLSTENT™ 83
THROMBECTOMY	Thrombectomy System AngioJet® Ultra Thrombectomy Console 84 Thrombectomy Catheters SmartClaw® Thrombectomy Catheter ZelanteDVT™ Thrombectomy Set Ultra 84 Solent® Proxi 84 Solent® Omni 84 Solent® Dista 84 AVX® 84 Accessories Ultra Power Pulse® Delivery Kit 84
THROMBOLYSIS	EKOS™ Endovascular System Control Unit 4.0 85 EKOS™ Endovascular System 85 EKOS+™ Endovascular System 85



	Page
ULTRASOUND	Peripheral Imaging Catheters and Accessories 86
	Avvigo™ + Multi-Modality Guidance System (Integrated Version) 86
	Avvigo™ + Multi-Modality Guidance System (Cart Version) 86
	Avvigo™ + Guidance System Accessories 86
VARITHENA	Varithena 87

RADIOFREQUENCY ABLATION

LeVeen™ Needle Electrode

UPN / Order Code	Catalog No.	GTIN	Diameter (cm)	Length (cm)	Gauge	C-Code	Price
M001262160	26-216	08714729454915	5	15	13	N/A	\$2,300
M001262170	26-217	08714729767619	5	25	13	N/A	\$2,300
M001262130	26-213	08714729295594	4	15	14	N/A	\$2,070
M001262310	26-231	08714729120735	4	25	14	N/A	\$2,070
M001262020	26-202	08714729300601	3.5	12	15	N/A	\$1,835
M001262030	26-203	08714729300618	3.5	15	15	N/A	\$1,835
M001262150	26-215	08714729349594	3.5	25	15	N/A	\$1,835
M001262040	26-204	08714729300625	3	12	15	N/A	\$1,610
M001262050	26-205	08714729300632	3	15	15	N/A	\$1,610

LeVeen SuperSlim™ Needle Electrode

UPN / Order Code	Catalog No.	GTIN	Diameter (cm)	Length (cm)	Gauge	C-Code	Price
M001262290	26-229	08714729441724	3	25	17	N/A	\$1,580
M001262280	26-228	08714729441717	3	15	17	N/A	\$1,580
M001262270	26-227	08714729441700	2	25	17	N/A	\$1,445
M001262260	26-226	08714729441694	2	15	17	N/A	\$1,445

LeVeen CoAccess™ Electrode System

UPN / Order Code	Catalog No.	GTIN	Diameter (cm)	Length (cm)	Gauge	C-Code	Price
M001262220	26-222	08714729454922	3	15	15	N/A	\$1,720
M001262230	26-223	08714729295617	3.5	15	15	N/A	\$1,950
M001262240	26-224	08714729454939	4	15	15	N/A	\$2,180
M001262250	26-225	08714729349570			16.5 (ID)	CoAccess Introducer Set N/A	\$145

Soloist™ Electrode

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Gauge	C-Code	Price
M001262500	26-250	08714729122487	0.9	18	16.5 (ID)	N/A	\$1,375

CRYOABLATION NEEDLES

Cryoablation CX Needles

Product Name	UPN	GTIN	Color Handle	Length (cm)	Gauge	C-Code*	Price
IcePearl 2.1 CX 90° Needle	FPRPR3601	07290104830990	White	17.5	14	C2618	\$2,550
IcePearl 2.1 CX L 90° Needle	FPRPR3617	07290104831072	White	23	14	C2618	\$2,550
IcePearl 2.1 CX Straight Needle	FPRPR3603	07290104831010	White	17.5	14	C2618	\$2,550
IceFORCE 2.1 CX 90° Needle	FPRPR3602	07290104831003	Gray	17.5	14	C2618	\$2,550
IceFORCE 2.1 CX L 90° Needle	FPRPR3618	07290104831089	Gray	23	14	C2618	\$2,550
IceFORCE 2.1 CX Straight Needle	FPRPR3604	07290104831027	Gray	17.5	14	C2618	\$2,550
IceSphere 1.5 CX 90° Needle	FPRPR3573	07290104831133	Yellow	17.5	17	C2618	\$2,550
IceRod 1.5 CX 90° Needle	FPRPR3533	07290104830587	Red	17.5	17	C2618	\$2,550
IceSeed 1.5 CX 90° Needle	H7493967433170	00191506032692	Black	17.5	17	C2618	\$2,550
IceSeed 1.5 CX S 90° Needle	H7493967233100	00191506032678	Black	10	17	C2618	\$2,550

*C-Code: Must be billed per needle used

CRYOABLATION NEEDLES

Cryoablation Classic Needles

Product Name	UPN	GTIN	Color Handle	Length (cm)	Gauge	C-Code*	Price
IceRod 1.5 i-Thaw Straight Needle	FPRPR4009	07290104830204	Red	17.5	17	C2618	\$2,445
IceSeed 1.5 Straight Needle	FPRPR3201	07290104830181	Black	17.5	17	C2618	\$2,135
IceSeed 1.5 90° Needle	FPRPR3202	07290104830235	Black	17.5	17	C2618	\$2,135
IceSphere 1.5 90° Short Needle	FPRPR3561	07290104830938	Yellow	10	17	C2618	\$2,135

*C-Code: Must be billed per needle used

Cryoablation MRI Needles

Product Name	UPN	GTIN	Color Handle	Length (cm)	Gauge	C-Code*	Price
IceSeed 1.5 Straight MRI Needle	FPRPR3192	07290104830266	Black	17.5	17	C2618	\$2,850
IceSeed 1.5 90° MRI Needle	FPRPR3194	07290104830280	Black	17.5	17	C2618	\$2,850
IceRod 1.5 Straight MRI Needle	FPRPR3193	07290104830273	Red	17.5	17	C2618	\$2,850
IceRod 1.5 90° MRI Needle	FPRPR3195	07290104830297	Red	17.5	17	C2618	\$2,850

*C-Code: Must be billed per needle used

CRYOABLATION SYSTEMS

Visual-ICE Cryoablation System

UPN	GTIN	Gas	Needles Supported	Number of Needles	C-Code	Price
FPRCH6000-02	00859019006417	Argon and Helium	All CX and Classic Needles	Up to 20	N/A	\$115,995

Visual-ICE MRI Cryoablation System

UPN	GTIN	Gas	Needles Supported	Number of Needles	C-Code	Price
H7493961070000	859019006554	Argon and Helium	MRI Compatible Needles	Up to 16	N/A	\$310,995

ICEfx Cryoablation System

UPN	GTIN	Gas	Needles Supported	Number of Needles	C-Code	Price
FPRCH8000-02	00859019006660	Argon	CX Needles, IceSeed Classic Needle*		N/A	\$89,995

*Passive Thaw must be used

CRYOABLATION ACCESSORIES

Product Name	UPN	GTIN	Description	C-Code	Price
MRI Mobile Connection Panel	H7493970071000	00859019006332	Mobile Connection Panel to support MRI needle	N/A	\$40,500
ICEfx Cryoablation System Cart	FPRCH8010	00859019006387	Cart to place ICEfx console on, sold separately from console	N/A	\$9,995
Multi-Point Thermal Sensor	FPRPR4010	07290104830303	One (1) Multi-Point Thermal Sensor 1.5 mm, individually packaged	N/A	\$1,195
EZ CONNECT2 Dual Cylinder Adapter, psi	FPRCH4004	00859019006196	One (1) Dual Cylinder Adapter Set in film pouch: One (1) 5 m Teflon™ lined high pressure gas supply line connected to a four-way adapter with an argon pressure gauge (psi) with blue rubber sleeve; One (1) 1.22 m Teflon™ lined high pressure gas supply line with quick connect fittings	N/A	\$4,995

ACCESS DEVICES AND ACCESSORIES

Introducers & Accessories

AccuStick™ II with RO Marker Introducer Systems

UPN / Order Code	Catalog No.	GTIN	Description	Units Per Box	C-Code	Price
M001207020	20-702	08714729157588	Accustick II with 0.038" J Tip Wire and 0.018" Stainless Steel Wire*	1	C1894	\$140
M001207030	20-703	08714729157601	Accustick II with 0.038" Straight Tip Wire and 0.018" Stainless Steel Wire*	1	C1894	\$140
M001207040	20-704	08714729194514	Accustick II with 0.018" Stainless Steel Wire*	1	C1894	\$115
M001207050	20-705	08714729157625	Accustick II without Guidewires*	1	C1894	\$90
M001207100	20-710	08714729201137	Accustick II with 0.038" J Tip Wire and 0.018" Nitinol Wire*	1	C1894	\$145
M001207111	20-711	08714729201144	Nitinol 60cm 0.018" Guidewire with 6cm Floppy Platinum Tip	5	C1894	\$305
M001201081	20-108	08714729115533	Stainless Steel 60cm 0.018" Guidewire with 7.5cm Floppy Platinum Tip	5	C1894	\$230
M001206041	20-604	08714729133766	21G Diagnostic Needle x 15cm	10	C1894	\$160

*Each Accustick II Includes: 1 Accustick II Introducer (Coaxial 6F Sheath with RO Marker/4F 20cm Dilator) with Locking Stiffening Cannula; 1 21G 15cm Diagnostic Needle with Stylet

FloSwitch™ HP Device

Flo Control Device

UPN / Order Code	Catalog No.	GTIN	Description	Units Per Box	C-Code	Price
M001442001	44-200	08714729115700	FloSwitch	12	N/A	\$90
M001442011	44-201	08714729049746	FloSwitch	24	N/A	\$170

Gateway™ Plus Y-Adapter

UPN / Order Code	Catalog No.	GTIN	Description	Units Per Box	C-Code	Price
M001153223	15-322	08714729190493	Gateway Plus Y-Adapter	10	N/A	\$410

Percufix™ Catheter Cuff Kit

All

UPN / Order Code	Catalog No.	GTIN	Description	Units Per Box	C-Code	Price
M001207001	20-700	08714729115526	PCC	10	N/A	\$215

Connecting Tube for Drainage Bag

UPN / Order Code	Catalog No.	GTIN	Description	Units Per Box	C-Code	Price
M001229101	22-910	08714729115564	Ctube	5	N/A	\$100

ESP™ Radiation Reduction Examination Gloves

UPN / Order Code	Catalog No.	GTIN	Description (Size)	C-Code	Price
M001201101	20-110	08714729175605	7	N/A	\$350
M001201111	20-111	08714729175612	7 1/2	N/A	\$350
M001201121	20-112	08714729175629	8	N/A	\$350
M001201131	20-113	08714729175636	8 1/2	N/A	\$350
M001201141	20-114	08714729175643	9	N/A	\$350

Packaged 3/box

Inflation Device

Encore™ 26 Inflation Device

UPN / Order Code	Catalog No.	GTIN	Description	Units Per Box	C-Code	Price
M001151050	15-105	08714729183624	Encore Inflation Device	1	N/A	\$97
M001151062	15-106	08714729137542	Encore Inflation Device	5	N/A	\$476

Encore™ 40 Inflation Device

UPN / Order Code	Catalog No.	GTIN	Description	Units Per Box	C-Code	Price
M001394472540		856877002445	Encore 40 Inflation Device	1	N/A	\$130

Introducer Sheaths

Super Sheath™ Introducer

4 F - 9 F Introducer Sheath and Guidewire

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Sheath Length (cm)	Guidewire Size (in)	Units Per Box	C-Code	Price
M00115710B1	15-710B	14543527181560	4	11	0.035	10	C1894	\$266
M00115711B1	15-711B	14543527181577	5	11	0.035	10	C1894	\$266
M00115712B1	15-712B	14543527181584	6	11	0.035	10	C1894	\$266
M00115713B1	15-713B	14543527181591	7	11	0.035	10	C1894	\$266
M00115714B1	15-714B	14543527181607	8	11	0.035	10	C1894	\$266
M00115715B1	15-715B	14543527181614	9	11	0.035	10	C1894	\$266

10 F - 14 F Introducer Sheath and Guidewire

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Sheath Length (cm)	Guidewire Size (in)	Units Per Box	C-Code	Price
M00115716B1	15-716B	14543527185193	10	11	0.035	10	C1894	\$266
M00115717B1	15-717B	14543527185209	11	11	0.035	10	C1894	\$266
M00115718B1	15-718B	14543527185216	12	11	0.035	10	C1894	\$266
M00115741B1	15-741B	14543527185223	14	11	0.035	10	C1894	\$266

4 F - 9 F Introducer Sheath (Guidewire Not Included)

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Sheath Length (cm)	Max Guidewire Size (in)	Units Per Box	C-Code	Price
M00115720B1	15-720B	14543527181621	4	11	0.038	10	C1894	\$255
M00115721B1	15-721B	14543527181638	5	11	0.038	10	C1894	\$255
M00115722B1	15-722B	14543527181645	6	11	0.038	10	C1894	\$255
M00115723B1	15-723B	14543527181652	7	11	0.038	10	C1894	\$255
M00115724B1	15-724B	14543527181669	8	11	0.038	10	C1894	\$255

Super Sheath™ Introducer (cont.)

4 F - 9 F Introducer Sheath (Guidewire Not Included)

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Sheath Length (cm)	Max Guidewire Size (in)	Units Per Box	C-Code	Price
M00115725B1	15-725B	14543527181676	9	11	0.038	10	C1894	\$255
M00115730B1	15-730B	14543527181683	5	25	0.038	10	C1894	\$473
M00115731B1	15-731B	14543527181690	6	25	0.038	10	C1894	\$473
M00115732B1	15-732B	14543527181706	7	25	0.038	10	C1894	\$473
M00115733B1	15-733B	14543527181713	8	25	0.038	10	C1894	\$473
M00115734B1	15-734B	14543527181720	9	25	0.038	10	C1894	\$473

10 F - 14 F Introducer Sheath (Guidewire Not Included)

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Sheath Length (cm)	Max Guidewire Size (in)	Units Per Box	C-Code	Price
M00115726B1	15-726B	14543527185230	10	11	0.038	10	C1894	\$255
M00115727B1	15-727B	14543527185247	11	11	0.038	10	C1894	\$255
M00115728B1	15-728B	14543527185254	12	11	0.038	10	C1894	\$255
M00115740B1	15-740B	14543527185261	14	11	0.038	10	C1894	\$255
M00115735B1	15-735B	14543527185278	10	25	0.038	10	C1894	\$473
M00115736B1	15-736B	14543527185285	11	25	0.038	10	C1894	\$473
M00115737B1	15-737B	14543527185292	12	25	0.038	10	C1894	\$473
M00115739B1	15-739B	14543527185308	14	25	0.038	10	C1894	\$473

4 F - 9 F R/O Sheath (Guidewire Not Included)

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Sheath Length (cm)	Max Guidewire Size (in)	Units Per Box	C-Code	Price
M00115956B1	15-956B	14543527181737	5	7	0.038	10	C1894	\$386
M00115957B1	15-957B	14543527181744	6	7	0.038	10	C1894	\$386
M00115958B1	15-958B	14543527181751	7	7	0.038	10	C1894	\$386
M00115959B1	15-959B	14543527181768	8	7	0.038	10	C1894	\$386
M00115955B1	15-955B	14543527181775	9	7	0.038	10	C1895	\$386
M00115960B1	15-960B	14543527181782	4	11	0.035	10	C1896	\$412
M00115961B1	15-961B	14543527181799	5	11	0.038	10	C1894	\$412
M00115962B1	15-962B	14543527181805	6	11	0.038	10	C1894	\$412
M00115963B1	15-963B	14543527181812	7	11	0.038	10	C1894	\$412
M00115964B1	15-964B	14543527181829	8	11	0.038	10	C1894	\$412
M00115969B1	15-969B	14543527181836	9	11	0.038	10	C1894	\$412
M00115970B1	15-970B	14543527181843	4	25	0.035	10	C1894	\$519
M00115971B1	15-971B	14543527181850	5	25	0.038	10	C1894	\$519
M00115965B1	15-965B	14543527181867	6	25	0.038	10	C1894	\$519
M00115966B1	15-966B	14543527181874	7	25	0.038	10	C1894	\$519
M00115967B1	15-967B	14543527181881	8	25	0.038	10	C1894	\$519
M00115968B1	15-968B	14543527181898	9	25	0.038	10	C1894	\$519

Catapult™ Guide Sheath

4 F - 8 F Introducer Sheath (Guidewire Not Included)

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Sheath Length (cm)	Max Guidewire Size (in)	Units Per Box	C-Code	Price
US-34015-F-ST-H	n/a	24251244501970	4	15	0.018 & 0.035	5	C1894	\$883
US-34045-F-ST-H	n/a	24251244501987	4	45	0.035	5	C1894	\$1,325
US-34060-F-ST-H	n/a	24251244501994	4	60	0.035	5	C1894	\$1,619
US-34090-F-ST-H	n/a	24251244502007	4	90	0.035	5	C1894	\$1,699
US-35015-F-ST-H	n/a	24251244502014	5	15	0.018 & 0.035	5	C1894	\$883
US-35045-F-ST-H	n/a	24251244502021	5	45	0.035	5	C1894	\$1,325
US-35045-F-MP-H	n/a	24251244502038	5	45	0.035	5	C1894	\$1,325
US-35045-F-RDC-H	n/a	24251244502045	5	45	0.035	5	C1894	\$1,325
US-35045-F-HS-H	n/a	24251244502052	5	45	0.035	5	C1894	\$1,325
US-35060-F-ST-H	n/a	24251244502069	5	60	0.035	5	C1894	\$1,619
US-35090-F-ST-H	n/a	24251244502076	5	90	0.035	5	C1894	\$1,699
US-36015-F-ST-H	n/a	24251244502083	6	15	0.018 & 0.035	5	C1894	\$883
US-36045-F-ST-H	n/a	24251244502090	6	45	0.035	5	C1894	\$1,325
US-36045-F-MP-H	n/a	24251244502106	6	45	0.035	5	C1894	\$1,325
US-36045-F-LIMA-H	n/a	24251244502113	6	45	0.035	5	C1894	\$1,325
US-36045-F-RDC-H	n/a	24251244502120	6	45	0.035	5	C1894	\$1,325
US-36045-F-HS-H	n/a	24251244502137	6	45	0.035	5	C1894	\$1,325
US-36060-F-ST-H	n/a	24251244502144	6	60	0.035	5	C1894	\$1,619
US-36090-F-ST-H	n/a	24251244502151	6	90	0.035	5	C1894	\$1,699
US-36090-F-MP-H	n/a	24251244502168	6	90	0.035	5	C1894	\$1,699
US-36090-F-RDC-H	n/a	24251244502175	6	90	0.035	5	C1894	\$1,699
US-36090-F-HS-H	n/a	24251244502182	6	90	0.035	5	C1894	\$1,699
US-37015-F-ST-H	n/a	24251244502199	7	15	0.018 & 0.035	5	C1894	\$883
US-37045-F-ST-H	n/a	24251244502205	7	45	0.035	5	C1894	\$1,325
US-37045-F-MP-H	n/a	24251244502212	7	45	0.035	5	C1894	\$1,325
US-37045-F-LIMA-H	n/a	24251244502229	7	45	0.035	5	C1894	\$1,325
US-37045-F-RDC-H	n/a	24251244502236	7	45	0.035	5	C1894	\$1,325
US-37045-F-HS-H	n/a	24251244502243	7	45	0.035	5	C1894	\$1,325
US-37060-F-ST-H	n/a	24251244502250	7	60	0.035	5	C1894	\$1,619
US-37090-F-ST-H	n/a	24251244502267	7	90	0.035	5	C1894	\$1,699
US-37090-F-MP-H	n/a	24251244502274	7	90	0.035	5	C1894	\$1,699
US-37090-F-HS-H	n/a	24251244502281	7	90	0.035	5	C1894	\$1,699
US-38015-F-ST-H	n/a	24251244502298	8	15	0.018 & 0.035	5	C1894	\$883
US-38045-F-ST-H	n/a	24251244502304	8	45	0.035	5	C1894	\$1,325
US-38060-F-ST-H	n/a	24251244502311	8	60	0.035	5	C1894	\$1,619
US-38090-F-ST-H	n/a	24251244502328	8	90	0.035	5	C1894	\$1,699

Cerebral Curves

Headhunter

Commonly used for carotid, innominate and subclavian access




Simmons


Facilitates access of common carotid, subclavian, renal and contralateral common iliac



Other Cerebral Curves




ST



BARN

Commonly used for guide wire redirect





BERN

Facilitates access of the great vessels

Visceral Curves


Cobra

Commonly used to access abdominal aortic vessels and internal iliacs



Renal

Best suited for steeply angled renal arteries



Chuang

Longer tips facilitate distal access to the celiac and mesenteric systems



Other Visceral Curves



MIK

Assists with access to perpendicular ostia



CONTRA2

Commonly used for contralateral access



SHK 1.0

Commonly used for contralateral access

Imager™ II Select Catheters

Imager™ II - 5 F Select Catheters

Braided Torque Cerebral Selective Catheters

UPN / Order Code	Catalog No.	GTINs	French Size	Tip Shape	Length (cm)	Guidewire Compatibility (in)	Maximum Flow Rate (mL/sec)	C-Code	Price
M001314001	31-400	08714729354826	5	ST	65	0.038	20	N/A	\$155
M001314011	31-401	08714729354833	5	ST	100	0.038	18	N/A	\$155
M001314021	31-402	08714729354840	5	ST	65	0.035	20	N/A	\$155
M001314031	31-403	08714729354857	5	ST	100	0.035	18	N/A	\$155
M001314041	31-404	08714729354864	5	BARNH	100	0.038	18	N/A	\$155
M001314051	31-405	08714729354871	5	BERN	100	0.038	18	N/A	\$155
M001314061	31-406	08714729354888	5	BERN	100	0.035	18	N/A	\$155
M001314081	31-408	08714729354901	5	BERN	40	0.038	20	N/A	\$155
M001314101	31-410	08714729354925	5	BERN	65	0.038	20	N/A	\$155
M001314141	31-414	08714729354963	5	H1	100	0.038	12	N/A	\$155
M001314151	31-415	08714729354970	5	H1	100	0.035	12	N/A	\$155
M001314301	31-430	08714729355120	5	SIM1	65	0.038	20	N/A	\$155
M001314321	31-432	08714729355144	5	SIM1	100	0.038	20	N/A	\$155
M001314341	31-434	08714729355168	5	SIM2	100	0.038	20	N/A	\$155
M001314351	31-435	08714729355175	5	SIM2	100	0.035	20	N/A	\$155
M001314361	31-436	08714729355182	5	SIM3	100	0.038	20	N/A	\$155

Packed 5 per box

Braided Torque Visceral Selective Catheters

UPN / Order Code	Catalog No.	GTINs	French Size	Tip Shape	Length (cm)	Guidewire Compatibility (in)	Maximum Flow Rate (mL/sec)	C-Code	Price
M001314521	31-452	08714729355342	5	CHGB	65	0.038	10	N/A	\$155
M001314581	31-458	08714729355403	5	CONTRA 2	65	0.038	10	N/A	\$155
M001314621	31-462	08714729355441	5	C1	65	0.038	20	N/A	\$155
M001314661	31-466	08714729355489	5	C2	65	0.038	20	N/A	\$155
M001314671	31-467	08714729355496	5	C2	65	0.035	20	N/A	\$155
M001314761	31-476	08714729355588	5	MIK	65	0.038	10	N/A	\$155
M001314771	31-477	08714729355595	5	MIK	65	0.035	10	N/A	\$155
M001314821	31-482	08714729355649	5	SHK1.0	65	0.038	10	N/A	\$155
M001314861	31-486	08714729355687	5	RC2	65	0.038	10	N/A	\$155
M001314881	31-488	08714729355700	5	RDC	65	0.038	20	N/A	\$155

Packed 5 per box

† Distal segment has 2 side holes

Coyote™ Over-the-Wire Balloon Dilatation Catheter

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939186150490	39186-15049	08714729797111	1.5	40	4	6	14	C1725	\$690
H74939186200690	39186-20069	08714729797128	2	60	4	6	14	C1725	\$690
H74939186200890	39186-20089	08714729796435	2	80	4	6	14	C1725	\$690
H74939186201090	39186-20109	08714729797142	2	100	4	6	14	C1725	\$690
H74939186201290	39186-20129	08714729796442	2	120	4	6	14	C1725	\$795
H74939186201590	39186-20159	08714729797166	2	150	4	6	14	C1725	\$795
H74939186202290	39186-20229	08714729797173	2	220	4	6	14	C1725	\$835
H74939186250690	39186-25069	08714729796480	2.5	60	4	6	14	C1725	\$690
H74939186250890	39186-25089	08714729796497	2.5	80	4	6	14	C1725	\$690
H74939186251090	39186-25109	08714729796510	2.5	100	4	6	14	C1725	\$690
H74939186251290	39186-25129	08714729797197	2.5	120	4	6	14	C1725	\$795
H74939186251590	39186-25159	08714729797203	2.5	150	4	6	14	C1725	\$795
H74939186252290	39186-25229	08714729797210	2.5	220	4	6	14	C1725	\$835
H74939186300690	39186-30069	08714729797234	3	60	4	6	14	C1725	\$690
H74939186300890	39186-30089	08714729796558	3	80	4	6	14	C1725	\$690
H74939186301090	39186-30109	08714729797241	3	100	4	6	14	C1725	\$690
H74939186301290	39186-30129	08714729796572	3	120	4	6	14	C1725	\$795
H74939186301590	39186-30159	08714729797258	3	150	4	6	14	C1725	\$795
H74939186302290	39186-30229	08714729797272	3	220	4	6	14	C1725	\$835
H74939186350690	39186-35069	08714729797296	3.5	60	4	6	14	C1725	\$690
H74939186350890	39186-35089	08714729796602	3.5	80	4	6	14	C1725	\$690
H74939186351090	39186-35109	08714729796619	3.5	100	4	6	14	C1725	\$690
H74939186351290	39186-35129	08714729796626	3.5	120	4	6	14	C1725	\$795
H74939186351590	39186-35159	08714729797319	3.5	150	4	6	14	C1725	\$795
H74939186352290	39186-35229	08714729797333	3.5	220	4	6	14	C1725	\$835
H74939186400690	39186-40069	08714729797340	4	60	4	6	14	C1725	\$690
H74939186400890	39186-40089	08714729797357	4	80	4	6	14	C1725	\$690
H74939186401090	39186-40109	08714729796657	4	100	4	6	14	C1725	\$690
H74939186401290	39186-40129	08714729797388	4	120	4	6	14	C1725	\$795
H74939186401590	39186-40159	08714729797395	4	150	4	6	14	C1725	\$795
H74939186402290	39186-40229	08714729797401	4	220	4	6	14	C1725	\$835

150 cm Usable Length

H74939186150410	39186-15041	08714729796411	1.5	40	4	6	14	C1725	\$690
H74939186200610	39186-20061	08714729796732	2	60	4	6	14	C1725	\$690
H74939186200810	39186-20081	08714729796428	2	80	4	6	14	C1725	\$690
H74939186201010	39186-20101	08714729797135	2	100	4	6	14	C1725	\$690
H74939186201210	39186-20121	08714729797159	2	120	4	6	14	C1725	\$795
H74939186201510	39186-20151	08714729796459	2	150	4	6	14	C1725	\$795
H74939186202210	39186-20221	08714729796466	2	220	4	6	14	C1725	\$835
H74939186250610	39186-25061	08714729796473	2.5	60	4	6	14	C1725	\$690
H74939186250810	39186-25081	08714729797180	2.5	80	4	6	14	C1725	\$690
H74939186251010	39186-25101	08714729796503	2.5	100	4	6	14	C1725	\$690
H74939186251210	39186-25121	08714729796527	2.5	120	4	6	14	C1725	\$795
H74939186251510	39186-25151	08714729796534	2.5	150	4	6	14	C1725	\$795

Coyote™ Over-the-Wire Balloon Dilatation Catheter (cont.)

150 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939186252210	39186-25221	08714729796718	2.5	220	4	6	14	C1725	\$835
H74939186300610	39186-30061	08714729797227	3	60	4	6	14	C1725	\$690
H74939186300810	39186-30081	08714729796541	3	80	4	6	14	C1725	\$690
H74939186301010	39186-30101	08714729796565	3	100	4	6	14	C1725	\$690
H74939186301210	39186-30121	08714729796725	3	120	4	6	14	C1725	\$795
H74939186301510	39186-30151	08714729796589	3	150	4	6	14	C1725	\$795
H74939186302210	39186-30221	08714729797265	3	220	4	6	14	C1725	\$835
H74939186350610	39186-35061	08714729797289	3.5	60	4	6	14	C1725	\$690
H74939186350810	39186-35081	08714729796596	3.5	80	4	6	14	C1725	\$690
H74939186351010	39186-35101	08714729797302	3.5	100	4	6	14	C1725	\$690
H74939186351210	39186-35121	08714729796756	3.5	120	4	6	14	C1725	\$795
H74939186351510	39186-35151	08714729796749	3.5	150	4	6	14	C1725	\$795
H74939186352210	39186-35221	08714729797326	3.5	220	4	6	14	C1725	\$835
H74939186400610	39186-40061	08714729796633	4	60	4	6	14	C1725	\$690
H74939186400810	39186-40081	08714729796640	4	80	4	6	14	C1725	\$690
H74939186401010	39186-40101	08714729797364	4	100	4	6	14	C1725	\$690
H74939186401210	39186-40121	08714729797371	4	120	4	6	14	C1725	\$795
H74939186401510	39186-40151	08714729796664	4	150	4	6	14	C1725	\$795
H74939186402210	39186-40221	08714729796671	4	220	4	6	14	C1725	\$835

Coyote™ Monorail Balloon Dilatation Catheter

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939185150490	39185-15049	08714729796770	1.5	40	4	6	14	C1725	\$690
H74939185200690	39185-20069	08714729796152	2	60	4	6	14	C1725	\$690
H74939185200890	39185-20089	08714729796176	2	80	4	6	14	C1725	\$690
H74939185201090	39185-20109	08714729796794	2	100	4	6	14	C1725	\$690
H74939185201290	39185-20129	08714729796190	2	120	4	6	14	C1725	\$795
H74939185201590	39185-20159	08714729796206	2	150	4	6	14	C1725	\$795
H74939185202290	39185-20229	08714729796213	2	220	4	6	14	C1725	\$835
H74939185250690	39185-25069	08714729796831	2.5	60	4	6	14	C1725	\$690
H74939185250890	39185-25089	08714729796855	2.5	80	4	6	14	C1725	\$690
H74939185251090	39185-25109	08714729796220	2.5	100	4	6	14	C1725	\$690
H74939185251290	39185-25129	08714729796237	2.5	120	4	6	14	C1725	\$795
H74939185251590	39185-25159	08714729796244	2.5	150	4	6	14	C1725	\$795
H74939185252290	39185-25229	08714729796251	2.5	220	4	6	14	C1725	\$835
H74939185300690	39185-30069	08714729796909	3	60	4	6	14	C1725	\$690
H74939185300890	39185-30089	08714729796275	3	80	4	6	14	C1725	\$690
H74939185301090	39185-30109	08714729796923	3	100	4	6	14	C1725	\$690
H74939185301290	39185-30129	08714729796299	3	120	4	6	14	C1725	\$795
H74939185301590	39185-30159	08714729796954	3	150	4	6	14	C1725	\$795
H74939185302290	39185-30229	08714729796978	3	220	4	6	14	C1725	\$835
H74939185350690	39185-35069	08714729796305	3.5	60	4	6	14	C1725	\$690
H74939185350890	39185-35089	08714729796329	3.5	80	4	6	14	C1725	\$690
H74939185351090	39185-35109	08714729797005	3.5	100	4	6	14	C1725	\$690
H74939185351290	39185-35129	08714729797012	3.5	120	4	6	14	C1725	\$795

Coyote™ Monorail Balloon Dilatation Catheter (cont.)

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939185351590	39185-35159	08714729796343	3.5	150	4	6	14	C1725	\$795
H74939185352290	39185-35229	08714729797043	3.5	220	4	6	14	C1725	\$835
H74939185400690	39185-40069	08714729796367	4	60	4	6	14	C1725	\$690
H74939185400890	39185-40089	08714729796374	4	80	4	6	14	C1725	\$690
H74939185401090	39185-40109	08714729797074	4	100	4	6	14	C1725	\$690
H74939185401290	39185-40129	08714729796381	4	120	4	6	14	C1725	\$795
H74939185401590	39185-40159	08714729796398	4	150	4	6	14	C1725	\$795
H74939185402290	39185-40229	08714729797104	4	220	4	6	14	C1725	\$835

150 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939185150410	39185-15041	08714729796763	1.5	40	4	6	14	C1725	\$690
H74939185200610	39185-20061	08714729796145	2	60	4	6	14	C1725	\$690
H74939185200810	39185-20081	08714729796169	2	80	4	6	14	C1725	\$690
H74939185201010	39185-20101	08714729796787	2	100	4	6	14	C1725	\$690
H74939185201210	39185-20121	08714729796183	2	120	4	6	14	C1725	\$795
H74939185201510	39185-20151	08714729796800	2	150	4	6	14	C1725	\$795
H74939185202210	39185-20221	08714729796817	2	220	4	6	14	C1725	\$835
H74939185250610	39185-25061	08714729796824	2.5	60	4	6	14	C1725	\$690
H74939185250810	39185-25081	08714729796848	2.5	80	4	6	14	C1725	\$690
H74939185251010	39185-25101	08714729796862	2.5	100	4	6	14	C1725	\$690
H74939185251210	39185-25121	08714729796879	2.5	120	4	6	14	C1725	\$795
H74939185251510	39185-25151	08714729796886	2.5	150	4	6	14	C1725	\$795
H74939185252210	39185-25221	08714729796893	2.5	220	4	6	14	C1725	\$835
H74939185300610	39185-30061	08714729796268	3	60	4	6	14	C1725	\$690
H74939185300810	39185-30081	08714729796916	3	80	4	6	14	C1725	\$690
H74939185301010	39185-30101	08714729796282	3	100	4	6	14	C1725	\$690
H74939185301210	39185-30121	08714729796930	3	120	4	6	14	C1725	\$795
H74939185301510	39185-30151	08714729796947	3	150	4	6	14	C1725	\$795
H74939185302210	39185-30221	08714729796961	3	220	4	6	14	C1725	\$835
H74939185350610	39185-35061	08714729796985	3.5	60	4	6	14	C1725	\$690
H74939185350810	39185-35081	08714729796312	3.5	80	4	6	14	C1725	\$690
H74939185351010	39185-35101	08714729796992	3.5	100	4	6	14	C1725	\$690
H74939185351210	39185-35121	08714729796336	3.5	120	4	6	14	C1725	\$795
H74939185351510	39185-35151	08714729797029	3.5	150	4	6	14	C1725	\$795
H74939185352210	39185-35221	08714729797036	3.5	220	4	6	14	C1725	\$835
H74939185400610	39185-40061	08714729796350	4	60	4	6	14	C1725	\$690
H74939185400810	39185-40081	08714729797050	4	80	4	6	14	C1725	\$690
H74939185401010	39185-40101	08714729797067	4	100	4	6	14	C1725	\$690
H74939185401210	39185-40121	08714729797081	4	120	4	6	14	C1725	\$795
H74939185401510	39185-40151	08714729797098	4	150	4	6	14	C1725	\$795
H74939185402210	39185-40221	08714729796404	4	220	4	6	14	C1725	\$835

Coyote™ ES Over-the-Wire Balloon Dilatation Catheter

14 ATM Rated Burst Pressure

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Introducer Sheath (F)	Guide Catheter (F)	C-Code	Price
H74939134152010	39134-15201	08714729767329	1.5	20	142	4	6	C1725	\$670
H74939134202010	39134-20201	08714729767336	2	20	142	4	6	C1725	\$670
H74939134203010	39134-20301	08714729767343	2	30	143	4	6	C1725	\$670
H74939134204010	39134-20401	08714729767350	2	40	144	4	6	C1725	\$670
H74939134252010	39134-25201	08714729767367	2.5	20	143	4	6	C1725	\$670
H74939134253010	39134-25301	8714729767374	2.5	30	144	4	6	C1725	\$670
H74939134254010	39134-25401	8714729767381	2.5	40	145	4	6	C1725	\$670
H74939134302010	39134-30201	08714729767398	3	20	143	4	6	C1725	\$670
H74939134303010	39134-30301	08714729767404	3	30	144	4	6	C1725	\$670
H74939134304010	39134-30401	08714729767411	3	40	145	4	6	C1725	\$670

12 ATM Rated Burst Pressure

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Introducer Sheath (F)	Guide Catheter (F)	C-Code	Price
H74939134352010	39134-35201	08714729767428	3.5	20	143	4	6	C1725	\$670
H74939134353010	39134-35301	08714729767435	3.5	30	144	4	6	C1725	\$670
H74939134354010	39134-35401	08714729767442	3.5	40	145	4	6	C1725	\$670
H74939134402010	39134-40201	08714729767459	4	20	143	4	6	C1725	\$670
H74939134403010	39134-40301	08714729767466	4	30	144	4	6	C1725	\$670
H74939134404010	39134-40401	08714729767473	4	40	145	4	6	C1725	\$670

Coyote™ ES Monorail™ Balloon Dilatation Catheter

14 ATM Rated Burst Pressure

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Introducer Sheath (F)	Guide Catheter (F)	C-Code	Price
H74939135152010	39135-15201	08714729767169	1.5	20	143	4	5	C1725	\$695
H74939135202010	39135-20201	08714729767176	2	20	143	4	5	C1725	\$695
H74939135203010	39135-20301	08714729767183	2	30	144	4	5	C1725	\$695
H74939135204010	39135-20401	08714729767190	2	40	145	4	5	C1725	\$695
H74939135252010	39135-25201	08714729767206	2.5	20	144	4	5	C1725	\$695
H74939135253010	39135-25301	08714729767213	2.5	30	145	4	5	C1725	\$695
H74939135254010	39135-25401	08714729767220	2.5	40	146	4	5	C1725	\$695
H74939135302010	39135-30201	08714729767237	3	20	144	4	5	C1725	\$695
H74939135303010	39135-30301	08714729767244	3	30	145	4	5	C1725	\$695
H74939135304010	39135-30401	08714729767251	3	40	146	4	5	C1725	\$695

12 ATM Rated Burst Pressure

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Introducer Sheath (F)	Guide Catheter (F)	C-Code	Price
H74939135352010	39135-35201	08714729767268	3.5	20	144	4	5	C1725	\$695
H74939135353010	39135-35301	08714729767275	3.5	30	145	4	5	C1725	\$695
H74939135354010	39135-35401	08714729767282	3.5	40	146	4	5	C1725	\$695
H74939135402010	39135-40201	08714729767299	4	20	144	4	5	C1725	\$695
H74939135403010	39135-40301	08714729767305	4	30	145	4	5	C1725	\$695
H74939135404010	39135-40401	08714729767312	4	40	146	4	5	C1725	\$695

Sterling™ Over-the-Wire Balloon Dilatation Catheter

40 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939032402040	39032-40204	08714729123538	4	20	4	6	14	C1725	\$670
H74939032404040	39032-40404	08714729123620	4	40	4	6	14	C1725	\$670
H74939032502040	39032-50204	08714729122845	5	20	4	6	14	C1725	\$670
H74939032504040	39032-50404	08714729122913	5	40	4	6	14	C1725	\$670
H74939032602040	39032-60204	08714729123743	6	20	4	6	14	C1725	\$670
H74939032604040	39032-60404	08714729123781	6	40	4	6	14	C1725	\$670
H74939032702040	39032-70204	08714729124016	7	20	4	6	14	C1725	\$670
H74939032704040	39032-70404	08714729124092	7	40	4	6	14	C1725	\$670
H74939032802040	39032-80204	08714729720638	8	20	4	6	12	C1725	\$670
H74939032804040	39032-80404	08714729720683	8	40	4	6	12	C1725	\$670

80 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939032402080	39032-40208	08714729123545	4	20	4	6	14	C1725	\$670
H74939032403080	39032-40308	08714729123569	4	30	4	6	14	C1725	\$670
H74939032404080	39032-40408	08714729123637	4	40	4	6	14	C1725	\$670
H74939032406080	39032-40608	08714729123651	4	60	4	6	14	C1725	\$670
H74939032408080	39032-40808	08714729123675	4	80	4	6	14	C1725	\$670
H74939032401080	39032-40108	08714729123514	4	100	4	6	14	C1725	\$670
H74939032502080	39032-50208	08714729122852	5	20	4	6	14	C1725	\$670
H74939032503080	39032-50308	08714729122883	5	30	4	6	14	C1725	\$670
H74939032504080	39032-50408	08714729122920	5	40	4	6	14	C1725	\$670
H74939032506080	39032-50608	08714729122944	5	60	4	6	14	C1725	\$670
H74939032508080	39032-50808	08714729122968	5	80	4	6	14	C1725	\$670
H74939032501080	39032-50108	08714729122821	5	100	4	6	14	C1725	\$670
H74939032602080	39032-60208	08714729123750	6	20	4	6	14	C1725	\$670
H74939032603080	39032-60308	08714729123774	6	30	4	6	14	C1725	\$670
H74939032604080	39032-60408	08714729123798	6	40	4	6	14	C1725	\$670
H74939032606080	39032-60608	08714729123811	6	60	4	6	14	C1725	\$670
H74939032608080	39032-60808	08714729123965	6	80	5	7	14	C1725	\$670
H74939032601080	39032-60108	08714729123699	6	100	5	7	14	C1725	\$670
H74939032702080	39032-70208	08714729124047	7	20	4	6	14	C1725	\$670
H74939032703080	39032-70308	08714729124061	7	30	4	6	14	C1725	\$670
H74939032704080	39032-70408	08714729124108	7	40	4	6	14	C1725	\$670
H74939032706080	39032-70608	08714729124122	7	60	5	7	14	C1725	\$670
H74939032708080	39032-70808	08714729124146	7	80	6	8	14	C1725	\$670
H74939032701080	39032-70108	08714729123989	7	100	6	8	14	C1725	\$670
H74939032802080	39032-80208	08714729720645	8	20	4	6	12	C1725	\$670
H74939032803080	39032-80308	08714729720669	8	30	4	6	12	C1725	\$670
H74939032804080	39032-80408	08714729720690	8	40	4	6	12	C1725	\$670
H74939032806080	39032-80608	08714729720713	8	60	5	7	12	C1725	\$670
H74939032808080	39032-80808	08714729720737	8	80	6	8	12	C1725	\$670
H74939032902080	39032-90208	08714729720751	9	20	5	7	10	C1725	\$670
H74939032903080	39032-90308	08714729720775	9	30	5	7	10	C1725	\$670
H74939032904080	39032-90408	08714729720799	9	40	5	7	10	C1725	\$670
H74939032906080	39032-90608	08714729720812	9	60	5	7	10	C1725	\$670
H74939032908080	39032-90808	08714729720836	9	80	6	8	10	C1725	\$670
H74939032102080	39032-10208	08714729123415	10	20	5	7	10	C1725	\$670
H74939032103080	39032-10308	08714729123439	10	30	5	7	10	C1725	\$670
H74939032104080	39032-10408	08714729123453	10	40	5	7	10	C1725	\$670
H74939032106080	39032-10608	08714729123477	10	60	6	8	10	C1725	\$670
H74939032108080	39032-10808	08714729123491	10	80	6	8	10	C1725	\$670

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.018" Low Profile

Sterling™ Over-the-Wire Balloon Dilatation Catheter

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939032200290	39032-20029	08714729859901	2	20	4	6	14	C1725	\$670
H74939032200390	39032-20039	08714729859925	2	30	4	6	14	C1725	\$670
H74939032200490	39032-20049	08714729859949	2	40	4	6	14	C1725	\$670
H74939032200690	39032-20069	08714729859963	2	60	4	6	14	C1725	\$670
H74939032250290	39032-25029	08714729860006	2.5	20	4	6	14	C1725	\$670
H74939032250390	39032-25039	08714729860020	2.5	30	4	6	14	C1725	\$670
H74939032250490	39032-25049	08714729860044	2.5	40	4	6	14	C1725	\$670
H74939032250690	39032-25069	08714729860068	2.5	60	4	6	14	C1725	\$670
H74939032300290	39032-30029	08714729860105	3	20	4	6	14	C1725	\$670
H74939032300390	39032-30039	08714729860129	3	30	4	6	14	C1725	\$670
H74939032300490	39032-30049	08714729860143	3	40	4	6	14	C1725	\$670
H74939032300690	39032-30069	08714729860167	3	60	4	6	14	C1725	\$670
H74939032350290	39032-35029	08714729860204	3.5	20	4	6	14	C1725	\$670
H74939032350390	39032-35039	08714729860228	3.5	30	4	6	14	C1725	\$670
H74939032350490	39032-35049	08714729860242	3.5	40	4	6	14	C1725	\$670
H74939032350690	39032-35069	08714729860266	3.5	60	4	6	14	C1725	\$670
H74939032202290	39032-20229	08714729859987	2	220	4	6	14	C1725	\$800
H74939032252290	39032-25229	08714729860082	2.5	220	4	6	14	C1725	\$800
H74939032302290	39032-30229	08714729860181	3	220	4	6	14	C1725	\$800
H74939032352290	39032-35229	08714729860280	3.5	220	4	6	14	C1725	\$800
H74939032402290	39032-40229	08714729860327	4	220	4	6	14	C1725	\$800
H74939032501290	39032-50129	08714729845355	5	120	5	7	14	C1725	\$700
H74939032501590	39032-50159	08714729845379	5	150	5	7	14	C1725	\$700
H74939032502090	39032-50209	08714729845386	5	200	5	7	14	C1725	\$800
H74939032502290	39032-50229	08714729845416	5	220	5	7	14	C1725	\$800
H74939032601290	39032-60129	08714729845430	6	120	5	7	14	C1725	\$700
H74939032601590	39032-60159	08714729845454	6	150	5	7	14	C1725	\$700
H74939032602090	39032-60209	08714729845461	6	200	5	7	14	C1725	\$800
H74939032602290	39032-60229	08714729845492	6	220	5	7	14	C1725	\$800
H74939032701290	39032-70129	08714729845515	7	120	5	7	14	C1725	\$700
H74939032701590	39032-70159	08714729845539	7	150	5	7	14	C1725	\$700
H74939032702090	39032-70209	08714729845546	7	200	5	7	14	C1725	\$800
H74939032702290	39032-70229	08714729845577	7	220	5	7	14	C1725	\$800

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939032402010	39032-40201	08714729123521	4	20	4	6	14	C1725	\$670
H74939032403010	39032-40301	08714729123552	4	30	4	6	14	C1725	\$670
H74939032404010	39032-40401	08714729123576	4	40	4	6	14	C1725	\$670
H74939032406010	39032-40601	08714729123644	4	60	4	6	14	C1725	\$670
H74939032408010	39032-40801	08714729123668	4	80	4	6	14	C1725	\$670
H74939032401010	39032-40101	08714729123507	4	100	4	6	14	C1725	\$670
H74939032502010	39032-50201	08714729122838	5	20	4	6	14	C1725	\$670
H74939032503010	39032-50301	08714729122876	5	30	4	6	14	C1725	\$670
H74939032504010	39032-50401	08714729122906	5	40	4	6	14	C1725	\$670
H74939032506010	39032-50601	08714729122937	5	60	4	6	14	C1725	\$670
H74939032508010	39032-50801	08714729122951	5	80	4	6	14	C1725	\$670
H74939032501010	39032-50101	08714729122814	5	100	4	6	14	C1725	\$670
H74939032602010	39032-60201	08714729123729	6	20	4	6	14	C1725	\$670
H74939032603010	39032-60301	08714729123767	6	30	4	6	14	C1725	\$670
H74939032604010	39032-60401	08714729122784	6	40	4	6	14	C1725	\$670

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.018" Low Profile

Sterling™ Over-the-Wire Balloon Dilatation Catheter (cont.)

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939032606010	39032-60601	08714729123804	6	60	4	6	14	C1725	\$670
H74939032608010	39032-60801	08714729123828	6	80	5	7	14	C1725	\$670
H74939032601010	39032-60101	08714729123682	6	100	5	7	14	C1725	\$670
H74939032702010	39032-70201	08714729123996	7	20	4	6	14	C1725	\$670
H74939032703010	39032-70301	08714729124054	7	30	4	6	14	C1725	\$670
H74939032704010	39032-70401	08714729124078	7	40	4	6	14	C1725	\$670
H74939032706010	39032-70601	08714729124115	7	60	5	7	14	C1725	\$670
H74939032708010	39032-70801	08714729124139	7	80	6	8	14	C1725	\$670
H74939032701010	39032-70101	08714729123972	7	100	6	8	14	C1725	\$670
H74939032802010	39032-80201	08714729720621	8	20	4	6	12	C1725	\$670
H74939032803010	39032-80301	08714729720652	8	30	4	6	12	C1725	\$670
H74939032804010	39032-80401	08714729720676	8	40	4	6	12	C1725	\$670
H74939032806010	39032-80601	08714729720706	8	60	5	7	12	C1725	\$670
H74939032808010	39032-80801	08714729720720	8	80	6	8	12	C1725	\$670
H74939032902010	39032-90201	08714729720744	9	20	5	7	10	C1725	\$670
H74939032903010	39032-90301	08714729720768	9	30	5	7	10	C1725	\$670
H74939032904010	39032-90401	08714729720782	9	40	5	7	10	C1725	\$670
H74939032906010	39032-90601	08714729720805	9	60	5	7	10	C1725	\$670
H74939032908010	39032-90801	08714729720829	9	80	6	8	10	C1725	\$670
H74939032102010	39032-10201	08714729123408	10	20	5	7	10	C1725	\$670
H74939032103010	39032-10301	08714729123422	10	30	5	7	10	C1725	\$670
H74939032104010	39032-10401	08714729123446	10	40	5	7	10	C1725	\$670
H74939032106010	39032-10601	08714729123460	10	60	6	8	10	C1725	\$670
H74939032108010	39032-10801	08714729123484	10	80	6	8	10	C1725	\$670

150 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939032200210	39032-20021	08714729859895	2	20	4	6	14	C1725	\$670
H74939032200310	39032-20031	08714729859918	2	30	4	6	14	C1725	\$670
H74939032200410	39032-20041	08714729859932	2	40	4	6	14	C1725	\$670
H74939032200610	39032-20061	08714729859956	2	60	4	6	14	C1725	\$670
H74939032250210	39032-25021	08714729859994	2.5	20	4	6	14	C1725	\$670
H74939032250310	39032-25031	08714729860013	2.5	30	4	6	14	C1725	\$670
H74939032250410	39032-25041	08714729860037	2.5	40	4	6	14	C1725	\$670
H74939032250610	39032-25061	08714729860051	2.5	60	4	6	14	C1725	\$670
H74939032300210	39032-30021	08714729860099	3	20	4	6	14	C1725	\$670
H74939032300310	39032-30031	08714729860112	3	30	4	6	14	C1725	\$670
H74939032300410	39032-30041	08714729860136	3	40	4	6	14	C1725	\$670
H74939032300610	39032-30061	08714729860150	3	60	4	6	14	C1725	\$670
H74939032350210	39032-35021	08714729860198	3.5	20	4	6	14	C1725	\$670
H74939032350310	39032-35031	08714729860211	3.5	30	4	6	14	C1725	\$670
H74939032350410	39032-35041	08714729860235	3.5	40	4	6	14	C1725	\$670
H74939032350610	39032-35061	08714729860259	3.5	60	4	6	14	C1725	\$670
H74939032202210	39032-20221	08714729859970	2	220	4	6	14	C1725	\$800
H74939032252210	39032-25221	08714729860075	2.5	220	4	6	14	C1725	\$800
H74939032302210	39032-30221	08714729860174	3	220	4	6	14	C1725	\$800
H74939032352210	39032-35221	08714729860273	3.5	220	4	6	14	C1725	\$800
H74939032402210	39032-40221	08714729860310	4	220	4	6	14	C1725	\$800
H74939032501210	39032-50121	08714729845348	5	120	5	7	14	C1725	\$700
H74939032501510	39032-50151	08714729845362	5	150	5	7	14	C1725	\$700
H74939032502110	39032-50211	08714729845393	5	200	5	7	14	C1725	\$800
H74939032502210	39032-50221	08714729845409	5	220	5	7	14	C1725	\$800

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.018" Low Profile

Sterling™ Over-the-Wire Balloon Dilatation Catheter (cont.)

150 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939032601210	39032-60121	08714729845423	6	120	5	7	14	C1725	\$700
H74939032601510	39032-60151	08714729845447	6	150	5	7	14	C1725	\$700
H74939032602110	39032-60211	08714729845478	6	200	5		14	C1725	\$800
H74939032602210	39032-60221	08714729845485	6	220	5	7	14	C1725	\$800
H74939032701210	39032-70121	08714729845508	7	120	5	7	14	C1725	\$700
H74939032701510	39032-70151	08714729845522	7	150	5	7	14	C1725	\$700
H74939032702110	39032-70211	08714729845553	7	200	5	7	14	C1725	\$800
H74939032702210	39032-70221	08714729845560	7	220	5	7	14	C1725	\$800

Sterling™ Monorail™ Balloon Dilatation Catheter

80 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939031401080	39031-40108	08714729756934	4	10	4	6	14	C1725	\$695
H74939031401580	39031-40158	08714729756941	4	15	4	6	14	C1725	\$695
H74939031402080	39031-40208	08714729756965	4	20	4	6	14	C1725	\$695
H74939031403080	39031-40308	08714729756989	4	30	4	6	14	C1725	\$695
H74939031404080	39031-40408	08714729757009	4	40	4	6	14	C1725	\$695
H74939031451580	39031-45158	08714729757023	4.5	15	4	6	14	C1725	\$695
H74939031452080	39031-45208	08714729757047	4.5	20	4	6	14	C1725	\$695
H74939031453080	39031-45308	08714729757061	4.5	30	4	6	14	C1725	\$695
H74939031454080	39031-45408	08714729757085	4.5	40	4	6	14	C1725	\$695
H74939031501080	39031-50108	08714729757092	5	10	4	6	14	C1725	\$695
H74939031501580	39031-50158	08714729757108	5	15	4	6	14	C1725	\$695
H74939031502080	39031-50208	08714729757115	5	20	4	6	14	C1725	\$695
H74939031503080	39031-50308	08714729757139	5	30	4	6	14	C1725	\$695
H74939031504080	39031-50408	08714729757153	5	40	4	6	14	C1725	\$695
H74939031551580	39031-55158	08714729757177	5.5	15	4	6	14	C1725	\$695
H74939031552080	39031-55208	08714729757184	5.5	20	4	6	14	C1725	\$695
H74939031553080	39031-55308	08714729757207	5.5	30	4	6	14	C1725	\$695
H74939031554080	39031-55408	08714729757221	5.5	40	4	6	14	C1725	\$695
H74939031601080	39031-60108	08714729757238	6	10	4	6	14	C1725	\$695
H74939031601580	39031-60158	08714729757245	6	15	4	6	14	C1725	\$695
H74939031602080	39031-60208	08714729757269	6	20	4	6	14	C1725	\$695
H74939031603080	39031-60308	08714729757283	6	30	4	6	14	C1725	\$695
H74939031604080	39031-60408	08714729757306	6	40	4	6	14	C1725	\$695
H74939031651580	39031-65158	08714729757320	6.5	15	4	6	14	C1725	\$695
H74939031652080	39031-65208	08714729757344	6.5	20	4	6	14	C1725	\$695
H74939031653080	39031-65308	08714729757368	6.5	30	4	6	14	C1725	\$695
H74939031654080	39031-65408	08714729757382	6.5	40	4	6	14	C1725	\$695
H74939031701080	39031-70108	08714729757399	7	10	4	6	14	C1725	\$695
H74939031701580	39031-70158	08714729757405	7	15	4	6	14	C1725	\$695
H74939031702080	39031-70208	08714729757429	7	20	4	6	14	C1725	\$695
H74939031703080	39031-70308	08714729757436	7	30	4	6	14	C1725	\$695
H74939031704080	39031-70408	08714729757450	7	40	4	6	14	C1725	\$695
H74939031801080	39031-80108	08714729757474	8	10	4	6	12	C1725	\$695
H74939031801580	39031-80158	08714729757481	8	15	4	6	12	C1725	\$695
H74939031802080	39031-80208	08714729757504	8	20	4	6	12	C1725	\$695
H74939031803080	39031-80308	08714729757528	8	30	4	6	12	C1725	\$695
H74939031804080	39031-80408	08714729757542	8	40	4	6	12	C1725	\$695

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.018" Low Profile

Sterling™ Monorail™ Balloon Dilatation Catheter

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939031302010	39031-30201	08714729755586	3	20	4	6	14	C1725	\$695
H74939031303010	39031-30301	08714729756873	3	30	4	6	14	C1725	\$695
H74939031304010	39031-30401	08714729756880	3	40	4	6	14	C1725	\$695
H74939031306010	39031-30601	08714729756897	3	60	4	6	14	C1725	\$695
H74939031352010	39031-35201	08714729756903	3.5	20	4	6	14	C1725	\$695
H74939031353010	39031-35301	08714729756910	3.5	30	4	6	14	C1725	\$695
H74939031354010	39031-35401	08714729756927	3.5	40	4	6	14	C1725	\$695
H74939031402010	39031-40201	08714729756958	4	20	4	6	14	C1725	\$695
H74939031403010	39031-40301	08714729756972	4	30	4	6	14	C1725	\$695
H74939031404010	39031-40401	08714729756996	4	40	4	6	14	C1725	\$695
H74939031406010	39031-40601	08714729757016	4	60	4	6	14	C1725	\$695
H74939031452010	39031-45201	08714729757030	4.5	20	4	6	14	C1725	\$695
H74939031453010	39031-45301	08714729757054	4.5	30	4	6	14	C1725	\$695
H74939031454010	39031-45401	08714729757078	4.5	40	4	6	14	C1725	\$695
H74939031502010	39031-50201	08714729756194	5	20	4	6	14	C1725	\$695
H74939031503010	39031-50301	08714729757122	5	30	4	6	14	C1725	\$695
H74939031504010	39031-50401	08714729757146	5	40	4	6	14	C1725	\$695
H74939031506010	39031-50601	08714729757160	5	60	4	6	14	C1725	\$695
H74939031552010	39031-55201	08714729755579	5.5	20	4	6	14	C1725	\$695
H74939031553010	39031-55301	08714729757191	5.5	30	4	6	14	C1725	\$695
H74939031554010	39031-55401	08714729757214	5.5	40	4	6	14	C1725	\$695
H74939031602010	39031-60201	08714729757252	6	20	4	6	14	C1725	\$695
H74939031603010	39031-60301	08714729757276	6	30	4	6	14	C1725	\$695
H74939031604010	39031-60401	08714729757290	6	40	4	6	14	C1725	\$695
H74939031606010	39031-60601	08714729757313	6	60	4	6	14	C1725	\$695
H74939031652010	39031-65201	08714729757337	6.5	20	4	6	14	C1725	\$695
H74939031653010	39031-65301	08714729757351	6.5	30	4	6	14	C1725	\$695
H74939031654010	39031-65401	08714729757375	6.5	40	4	6	14	C1725	\$695
H74939031702010	39031-70201	08714729757412	7	20	4	6	14	C1725	\$695
H74939031703010	39031-70301	08714729757566	7	30	4	6	14	C1725	\$695
H74939031704010	39031-70401	08714729757443	7	40	4	6	14	C1725	\$695
H74939031706010	39031-70601	08714729757467	7	60	5	7	14	C1725	\$695
H74939031802010	39031-80201	08714729757498	8	20	4	6	12	C1725	\$695
H74939031803010	39031-80301	08714729757511	8	30	4	6	12	C1725	\$695
H74939031804010	39031-80401	08714729757535	8	40	4	6	12	C1725	\$695
H74939031806010	39031-80601	08714729757559	8	60	5	7	12	C1725	\$695

150 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Guide Catheter (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939031202010	39031-20201	08714729859048	2	20	4	6	14	C1725	\$695
H74939031204010	39031-20401	08714729859086	2	40	4	6	14	C1725	\$695
H74939031206010	39031-20601	08714729859109	2	60	4	6	14	C1725	\$695
H74939031202210	39031-20221	08714729859062	2	220	4	6	14	C1725	\$800
H74939031252010	39031-25201	08714729859123	2.5	20	4	6	14	C1725	\$695
H74939031254010	39031-25401	08714729859161	2.5	40	4	6	14	C1725	\$695
H74939031256010	39031-25601	08714729859185	2.5	60	4	6	14	C1725	\$695
H74939031252210	39031-25221	08714729859147	2.5	220	4	6	14	C1725	\$800
H74939031302210	39031-30221	08714729859208	3	220	4	6	14	C1725	\$800
H74939031356010	39031-35601	08714729859246	3.5	60	4	6	14	C1725	\$695
H74939031352210	39031-35221	08714729859222	3.5	220	4	6	14	C1725	\$800
H74939031402210	39031-40221	08714729859291	4	220	4	6	14	C1725	\$800

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.018" Low Profile

Sterling™ SL Over-the-Wire Balloon Dilatation Catheter

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939148208090	39148-20809	08714729782421	2	80	4	14	C1725	\$650
H74939148201090	39148-20109	08714729782391	2	100	4	14	C1725	\$650
H74939148201290	39148-20129	08714729782407	2	120	4	14	C1725	\$745
H74939148201590	39148-20159	08714729782414	2	150	4	14	C1725	\$745
H74939148258090	39148-25809	08714729782469	2.5	80	4	14	C1725	\$650
H74939148251090	39148-25109	08714729782438	2.5	100	4	14	C1725	\$650
H74939148251290	39148-25129	08714729782445	2.5	120	4	14	C1725	\$745
H74939148251590	39148-25159	08714729782452	2.5	150	4	14	C1725	\$745
H74939148308090	39148-30809	08714729782506	3	80	4	14	C1725	\$650
H74939148301090	39148-30109	08714729782476	3	100	4	14	C1725	\$650
H74939148301290	39148-30129	08714729782483	3	120	4	14	C1725	\$745
H74939148301590	39148-30159	08714729782490	3	150	4	14	C1725	\$745
H74939148358090	39148-35809	08714729782544	3.5	80	4	14	C1725	\$650
H74939148351090	39148-35109	08714729782513	3.5	100	4	14	C1725	\$650
H74939148351290	39148-35129	08714729782520	3.5	120	4	14	C1725	\$745
H74939148351590	39148-35159	08714729782537	3.5	150	4	14	C1725	\$745
H74939148401290	39148-40129	08714729782568	4	120	4	14	C1725	\$745
H74939148401590	39148-40159	08714729782575	4	150	4	14	C1725	\$745

150 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939148208010	39148-20801	08714729782247	2	80	4	14	C1725	\$650
H74939148201010	39148-20101	08714729782216	2	100	4	14	C1725	\$650
H74939148201210	39148-20121	08714729782223	2	120	4	14	C1725	\$745
H74939148201510	39148-20151	08714729782230	2	150	4	14	C1725	\$745
H74939148258010	39148-25801	08714729782285	2.5	80	4	14	C1725	\$650
H74939148251010	39148-25101	08714729782254	2.5	100	4	14	C1725	\$650
H74939148251210	39148-25121	08714729782261	2.5	120	4	14	C1725	\$745
H74939148251510	39148-25151	08714729782278	2.5	150	4	14	C1725	\$745
H74939148308010	39148-30801	08714729782322	3	80	4	14	C1725	\$650
H74939148301010	39148-30101	08714729782292	3	100	4	14	C1725	\$650
H74939148301210	39148-30121	08714729782308	3	120	4	14	C1725	\$745
H74939148301510	39148-30151	08714729782315	3	150	4	14	C1725	\$745
H74939148358010	39148-35801	08714729782360	3.5	80	4	14	C1725	\$650
H74939148351010	39148-35101	08714729782339	3.5	100	4	14	C1725	\$650
H74939148351210	39148-35121	08714729782346	3.5	120	4	14	C1725	\$745
H74939148351510	39148-35151	08714729782353	3.5	150	4	14	C1725	\$745
H74939148408010	39148-40801	08714729782582	4	80	4	14	C1725	\$650
H74939148401010	39148-40101	08714729782551	4	100	4	14	C1725	\$650
H74939148401210	39148-40121	08714729782377	4	120	4	14	C1725	\$745
H74939148401510	39148-40151	08714729782384	4	150	4	14	C1725	\$745

Sterling™ SL Monorail™ Balloon Dilatation Catheter

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939147208090	39147-20809	08714729782827	2	80	4	14	C1725	\$675
H74939147201090	39147-20109	08714729782797	2	100	4	14	C1725	\$675
H74939147201290	39147-20129	08714729782803	2	120	4	14	C1725	\$780
H74939147201590	39147-20159	08714729782810	2	150	4	14	C1725	\$780

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.018" Low Profile

Sterling™ SL Monorail™ Balloon Dilatation Catheter (cont.)

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939147258090	39147-25809	08714729782865	2.5	80	4	14	C1725	\$675
H74939147251090	39147-25109	08714729782834	2.5	100	4	14	C1725	\$675
H74939147251290	39147-25129	08714729782841	2.5	120	4	14	C1725	\$780
H74939147251590	39147-25159	08714729782858	2.5	150	4	14	C1725	\$780
H74939147308090	39147-30809	08714729782902	3	80	4	14	C1725	\$675
H74939147301090	39147-30109	08714729782872	3	100	4	14	C1725	\$675
H74939147301290	39147-30129	08714729782889	3	120	4	14	C1725	\$780
H74939147301590	39147-30159	08714729782896	3	150	4	14	C1725	\$780
H74939147358090	39147-35809	08714729782940	3.5	80	4	14	C1725	\$675
H74939147351090	39147-35109	08714729782919	3.5	100	4	14	C1725	\$675
H74939147351290	39147-35129	08714729782926	3.5	120	4	14	C1725	\$780
H74939147351590	39147-35159	08714729782933	3.5	150	4	14	C1725	\$780
H74939147408090	39147-40809	08714729782988	4	80	4	14	C1725	\$675
H74939147401090	39147-40109	08714729782957	4	100	4	14	C1725	\$675
H74939147401290	39147-40129	08714729782964	4	120	4	14	C1725	\$780
H74939147401590	39147-40159	08714729782971	4	150	4	14	C1725	\$780

150 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939147208010	39147-20801	08714729782629	2	80	4	14	C1725	\$675
H74939147201010	39147-20101	08714729782599	2	100	4	14	C1725	\$675
H74939147201210	39147-20121	08714729782605	2	120	4	14	C1725	\$780
H74939147201510	39147-20151	08714729782612	2	150	4	14	C1725	\$780
H74939147258010	39147-25801	08714729782667	2.5	80	4	14	C1725	\$675
H74939147251010	39147-25101	08714729782636	2.5	100	4	14	C1725	\$675
H74939147251210	39147-25121	08714729782643	2.5	120	4	14	C1725	\$780
H74939147251510	39147-25151	08714729782650	2.5	150	4	14	C1725	\$780
H74939147308010	39147-30801	08714729782704	3	80	4	14	C1725	\$675
H74939147301010	39147-30101	08714729782674	3	100	4	14	C1725	\$675
H74939147301210	39147-30121	08714729782681	3	120	4	14	C1725	\$780
H74939147301510	39147-30151	08714729782698	3	150	4	14	C1725	\$780
H74939147358010	39147-35801	08714729782742	3.5	80	4	14	C1725	\$675
H74939147351010	39147-35101	08714729782711	3.5	100	4	14	C1725	\$675
H74939147351210	39147-35121	08714729782728	3.5	120	4	14	C1725	\$780
H74939147351510	39147-35151	08714729782735	3.5	150	4	14	C1725	\$780
H74939147408010	39147-40801	08714729782780	4	80	4	14	C1725	\$675
H74939147401010	39147-40101	08714729782759	4	100	4	14	C1725	\$675
H74939147401210	39147-40121	08714729782766	4	120	4	14	C1725	\$780
H74939147401510	39147-40151	08714729782773	4	150	4	14	C1725	\$780

0.035" High Pressure

Mustang™ Balloon Dilatation Catheter

40 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939171030240	39171-03024	08714729794639	3	20	5	24	C1725	\$460
H74939171030340	39171-03034	08714729795025	3	30	5	24	C1725	\$460
H74939171030440	39171-03044	08714729794646	3	40	5	24	C1725	\$460
H74939171030640	39171-03064	08714729795032	3	60	5	24	C1725	\$460

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" High Pressure

Mustang™ Balloon Dilatation Catheter (cont.)

40 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939171030840	39171-03084	08714729795049	3	80	5	24	C1725	\$460
H74939171031040	39171-03104	08714729795308	3	100	5	24	C1725	\$460
H74939171040240	39171-04024	08714729794172	4	20	5	24	C1725	\$460
H74939171040340	39171-04034	08714729795063	4	30	5	24	C1725	\$460
H74939171040440	39171-04044	08714729794103	4	40	5	24	C1725	\$460
H74939171040640	39171-04064	08714729795070	4	60	5	24	C1725	\$460
H74939171040840	39171-04084	08714729795087	4	80	5	24	C1725	\$460
H74939171041040	39171-04104	08714729795315	4	100	5	24	C1725	\$460
H74939171050240	39171-05024	08714729794066	5	20	5	24	C1725	\$460
H74939171050340	39171-05034	08714729795094	5	30	5	24	C1725	\$460
H74939171050440	39171-05044	08714729793731	5	40	5	24	C1725	\$460
H74939171050640	39171-05064	08714729795100	5	60	5	24	C1725	\$460
H74939171050840	39171-05084	08714729795117	5	80	5	24	C1725	\$460
H74939171051040	39171-05104	08714729795124	5	100	5	24	C1725	\$460
H74939171051240	39171-05124	08714729794882	5	120	5	24	C1725	\$560
H74939171060240	39171-06024	08714729793908	6	20	5	24	C1725	\$460
H74939171060340	39171-06034	08714729795131	6	30	5	24	C1725	\$460
H74939171060440	39171-06044	08714729793564	6	40	5	24	C1725	\$460
H74939171060640	39171-06064	08714729794783	6	60	5	24	C1725	\$460
H74939171060840	39171-06084	08714729795148	6	80	5	24	C1725	\$460
H74939171061040	39171-06104	08714729794899	6	100	5	24	C1725	\$460
H74939171061240	39171-06124	08714729794905	6	120	5	22	C1725	\$560
H74939171070240	39171-07024	08714729794370	7	20	5	20	C1725	\$460
H74939171070340	39171-07034	08714729795155	7	30	5	20	C1725	\$460
H74939171070440	39171-07044	08714729793588	7	40	5	20	C1725	\$460
H74939171070640	39171-07064	08714729794790	7	60	5	20	C1725	\$460
H74939171070840	39171-07084	08714729794806	7	80	5	20	C1725	\$460
H74939171071040	39171-07104	08714729794813	7	100	5	20	C1725	\$460
H74939171080240	39171-08024	08714729794660	8	20	6	20	C1725	\$460
H74939171080340	39171-08034	08714729794400	8	30	6	20	C1725	\$460
H74939171080440	39171-08044	08714729793632	8	40	6	20	C1725	\$460
H74939171080640	39171-08064	08714729794820	8	60	6	20	C1725	\$460
H74939171080840	39171-08084	08714729794837	8	80	6	20	C1725	\$460
H74939171081040	39171-08104	08714729794844	8	100	6	20	C1725	\$460
H74939171090240	39171-09024	08714729794912	9	20	6	18	C1725	\$460
H74939171090340	39171-09034	08714729795179	9	30	6	18	C1725	\$460
H74939171090440	39171-09044	08714729794622	9	40	6	18	C1725	\$460
H74939171090640	39171-09064	08714729795193	9	60	6	18	C1725	\$460
H74939171090840	39171-09084	08714729795216	9	80	6	18	C1725	\$460
H74939171100240	39171-10024	08714729795322	10	20	6	14	C1725	\$460
H74939171100340	39171-10034	08714729795339	10	30	6	14	C1725	\$460
H74939171100440	39171-10044	08714729794226	10	40	6	14	C1725	\$460
H74939171100640	39171-10064	08714729795247	10	60	6	14	C1725	\$460
H74939171100840	39171-10084	08714729794721	10	80	6	14	C1725	\$460
H74939171120240	39171-12024	08714729794929	12	20	7	14	C1725	\$460
H74939171120340	39171-12034	08714729795346	12	30	7	14	C1725	\$460
H74939171120440	39171-12044	08714729794592	12	40	7	14	C1725	\$460
H74939171120640	39171-12064	08714729795278	12	60	7	14	C1725	\$460
H74939171120840	39171-12084	08714729795292	12	80	7	14	C1725	\$460

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" High Pressure

Mustang™ Balloon Dilatation Catheter (cont.)

75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939171030270	39171-03027	08714729793991	3	20	5	24	C1725	\$460
H74939171030370	39171-03037	08714729794981	3	30	5	24	C1725	\$460
H74939171030470	39171-03047	08714729793946	3	40	5	24	C1725	\$460
H74939171030670	39171-03067	08714729794943	3	60	5	24	C1725	\$460
H74939171030870	39171-03087	08714729794950	3	80	5	24	C1725	\$460
H74939171031070	39171-03107	08714729794967	3	100	5	24	C1725	\$460
H74939171031270	39171-03127	08714729794554	3	120	5	24	C1725	\$560
H74939171031570	39171-03157	08714729794561	3	150	5	24	C1725	\$560
H74939171031870	39171-03187	08714729794547	3	180	5	24	C1725	\$610
H74939171032070	39171-03207	08714729794530	3	200	5	24	C1725	\$610
H74939171040270	39171-04027	08714729793496	4	20	5	24	C1725	\$460
H74939171040370	39171-04037	08714729794752	4	30	5	24	C1725	\$460
H74939171040470	39171-04047	08714729793403	4	40	5	24	C1725	\$460
H74939171040670	39171-04067	08714729794349	4	60	5	24	C1725	\$460
H74939171040870	39171-04087	08714729794363	4	80	5	24	C1725	\$460
H74939171041070	39171-04107	08714729794240	4	100	5	24	C1725	\$460
H74939171041270	39171-04127	08714729794158	4	120	5	24	C1725	\$560
H74939171041570	39171-04157	08714729794462	4	150	5	24	C1725	\$560
H74939171041870	39171-04187	08714729794479	4	180	5	24	C1725	\$610
H74939171042070	39171-04207	08714729794455	4	200	5	24	C1725	\$610
H74939171050270	39171-05027	08714729793427	5	20	5	24	C1725	\$460
H74939171050370	39171-05037	08714729794707	5	30	5	24	C1725	\$460
H74939171050470	39171-05047	08714729793366	5	40	5	24	C1725	\$460
H74939171050670	39171-05067	08714729793977	5	60	5	24	C1725	\$460
H74939171050870	39171-05087	08714729793960	5	80	5	24	C1725	\$460
H74939171051070	39171-05107	08714729793892	5	100	5	24	C1725	\$460
H74939171051270	39171-05127	08714729793915	5	120	5	24	C1725	\$560
H74939171051570	39171-05157	08714729794042	5	150	5	24	C1725	\$560
H74939171051870	39171-05187	08714729794486	5	180	5	24	C1725	\$610
H74939171052070	39171-05207	08714729794202	5	200	5	24	C1725	\$610
H74939171060270	39171-06027	08714729793397	6	20	5	24	C1725	\$460
H74939171060370	39171-06037	08714729794318	6	30	5	24	C1725	\$460
H74939171060470	39171-06047	08714729793335	6	40	5	24	C1725	\$460
H74939171060670	39171-06067	08714729793762	6	60	5	24	C1725	\$460
H74939171060870	39171-06087	08714729793854	6	80	5	24	C1725	\$460
H74939171061070	39171-06107	08714729793625	6	100	5	24	C1725	\$460
H74939171061270	39171-06127	08714729794127	6	120	5	22	C1725	\$560
H74939171061570	39171-06157	08714729793984	6	150	5	22	C1725	\$560
H74939171061870	39171-06187	08714729794493	6	180	5	22	C1725	\$610
H74939171062070	39171-06207	08714729794509	6	200	5	22	C1725	\$610
H74939171070270	39171-07027	08714729793502	7	20	5	20	C1725	\$460
H74939171070370	39171-07037	08714729794271	7	30	5	20	C1725	\$460
H74939171070470	39171-07047	08714729793342	7	40	5	20	C1725	\$460
H74939171070670	39171-07067	08714729793861	7	60	5	20	C1725	\$460
H74939171070870	39171-07087	08714729794264	7	80	5	20	C1725	\$460
H74939171071070	39171-07107	08714729793779	7	100	5	20	C1725	\$460
H74939171071270	39171-07127	08714729794424	7	120	5	18	C1725	\$560
H74939171071570	39171-07157	08714729794578	7	150	5	18	C1725	\$560
H74939171071870	39171-07187	08714729794608	7	180	5	18	C1725	\$610

Mustang™ Balloon Dilatation Catheter (cont)

75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939171072070	39171-07207	08714729794189	7	200	5	18	C1725	\$610
H74939171080270	39171-08027	08714729793519	8	20	6	20	C1725	\$460
H74939171080370	39171-08037	08714729793809	8	30	6	20	C1725	\$460
H74939171080470	39171-08047	08714729793328	8	40	6	20	C1725	\$460
H74939171080670	39171-08067	08714729793885	8	60	6	20	C1725	\$460
H74939171080870	39171-08087	08714729793656	8	80	6	20	C1725	\$460
H74939171081070	39171-08107	08714729793700	8	100	6	20	C1725	\$460
H74939171081270	39171-08127	08714729794431	8	120	6	18	C1725	\$560
H74939171081570	39171-08157	08714729794585	8	150	6	18	C1725	\$560
H74939171081870	39171-08187	08714729794615	8	180	6	18	C1725	\$610
H74939171082070	39171-08207	08714729794196	8	200	6	18	C1725	\$610
H74939171090270	39171-09027	08714729793755	9	20	6	18	C1725	\$460
H74939171090370	39171-09037	08714729794851	9	30	6	18	C1725	\$460
H74939171090470	39171-09047	08714729793489	9	40	6	18	C1725	\$460
H74939171090670	39171-09067	08714729794691	9	60	6	18	C1725	\$460
H74939171090870	39171-09087	08714729794868	9	80	6	18	C1725	\$460
H74939171100270	39171-10027	08714729793601	10	20	6	14	C1725	\$460
H74939171100370	39171-10037	08714729794394	10	30	6	14	C1725	\$460
H74939171100470	39171-10047	08714729793359	10	40	6	14	C1725	\$460
H74939171100670	39171-10067	08714729794004	10	60	6	14	C1725	\$460
H74939171100870	39171-10087	08714729794332	10	80	6	14	C1725	\$460
H74939171120270	39171-12027	08714729793922	12	20	7	14	C1725	\$460
H74939171120370	39171-12037	08714729794875	12	30	7	14	C1725	\$460
H74939171120470	39171-12047	08714729793441	12	40	7	14	C1725	\$460
H74939171120670	39171-12067	08714729794738	12	60	7	14	C1725	\$460
H74939171120870	39171-12087	08714729794745	12	80	7	14	C1725	\$460

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939171030210	39171-03021	08714729794141	3	20	5	24	C1725	\$460
H74939171030310	39171-03031	08714729794974	3	30	5	24	C1725	\$460
H74939171030410	39171-03041	08714729794035	3	40	5	24	C1725	\$460
H74939171030610	39171-03061	08714729794998	3	60	5	24	C1725	\$460
H74939171030810	39171-03081	08714729794936	3	80	5	24	C1725	\$460
H74939171031010	39171-03101	08714729794295	3	100	5	24	C1725	\$460
H74939171031210	39171-03121	08714729793816	3	120	5	24	C1725	\$560
H74939171031510	39171-03151	08714729794516	3	150	5	24	C1725	\$560
H74939171031810	39171-03181	08714729794523	3	180	5	24	C1725	\$610
H74939171032010	39171-03201	08714729794387	3	200	5	24	C1725	\$610
H74939171040210	39171-04021	08714729793724	4	20	5	24	C1725	\$460
H74939171040310	39171-04031	08714729795056	4	30	5	24	C1725	\$460
H74939171040410	39171-04041	08714729793557	4	40	5	24	C1725	\$460
H74939171040610	39171-04061	08714729794257	4	60	5	24	C1725	\$460
H74939171040810	39171-04081	08714729794028	4	80	5	24	C1725	\$460
H74939171041010	39171-04101	08714729793717	4	100	5	24	C1725	\$460
H74939171041210	39171-04121	08714729793670	4	120	5	24	C1725	\$560

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" High Pressure

Mustang™ Balloon Dilatation Catheter (cont)

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939171041510	39171-04151	08714729793694	4	150	5	24	C1725	\$560
H74939171041810	39171-04181	08714729794653	4	180	5	24	C1725	\$610
H74939171042010	39171-04201	08714729793595	4	200	5	24	C1725	\$610
H74939171050210	39171-05021	08714729793618	5	20	5	24	C1725	\$460
H74939171050310	39171-05031	08714729794769	5	30	5	24	C1725	\$460
H74939171050410	39171-05041	08714729793465	5	40	5	24	C1725	\$460
H74939171050610	39171-05061	08714729793878	5	60	5	24	C1725	\$460
H74939171050810	39171-05081	08714729793687	5	80	5	24	C1725	\$460
H74939171051010	39171-05101	08714729793533	5	100	5	24	C1725	\$460
H74939171051210	39171-05121	08714729793410	5	120	5	24	C1725	\$560
H74939171051510	39171-05151	08714729793434	5	150	5	24	C1725	\$560
H74939171051810	39171-05181	08714729794219	5	180	5	24	C1725	\$610
H74939171052010	39171-05201	08714729793373	5	200	5	24	C1725	\$610
H74939171060210	39171-06021	08714729793571	6	20	5	24	C1725	\$460
H74939171060310	39171-06031	08714729794776	6	30	5	24	C1725	\$460
H74939171060410	39171-06041	08714729793472	6	40	5	24	C1725	\$460
H74939171060610	39171-06061	08714729793939	6	60	5	24	C1725	\$460
H74939171060810	39171-06081	08714729793793	6	80	5	24	C1725	\$460
H74939171061010	39171-06101	08714729793540	6	100	5	24	C1725	\$460
H74939171061210	39171-06121	08714729793526	6	120	5	22	C1725	\$560
H74939171061510	39171-06151	08714729793458	6	150	5	22	C1725	\$560
H74939171061810	39171-06181	08714729794059	6	180	5	22	C1725	\$610
H74939171062010	39171-06201	08714729793380	6	200	5	22	C1725	\$610
H74939171070210	39171-07021	08714729793830	7	20	5	20	C1725	\$460
H74939171070310	39171-07031	08714729795001	7	30	5	20	C1725	\$460
H74939171070410	39171-07041	08714729793649	7	40	5	20	C1725	\$460
H74939171070610	39171-07061	08714729794325	7	60	5	20	C1725	\$460
H74939171070810	39171-07081	08714729794288	7	80	5	20	C1725	\$460
H74939171071010	39171-07101	08714729794011	7	100	5	20	C1725	\$460
H74939171071210	39171-07121	08714729794097	7	120	5	18	C1725	\$560
H74939171071510	39171-07151	08714729793786	7	150	5	18	C1725	\$560
H74939171071810	39171-07181	08714729794677	7	180	5	18	C1725	\$610
H74939171072010	39171-07201	08714729793748	7	200	5	18	C1725	\$610
H74939171080210	39171-08021	08714729793847	8	20	6	20	C1725	\$460
H74939171080310	39171-08031	08714729794417	8	30	6	20	C1725	\$460
H74939171080410	39171-08041	08714729793663	8	40	6	20	C1725	\$460
H74939171080610	39171-08061	08714729794714	8	60	6	20	C1725	\$460
H74939171080810	39171-08081	08714729794301	8	80	6	20	C1725	\$460
H74939171081010	39171-08101	08714729794233	8	100	6	20	C1725	\$460
H74939171081210	39171-08121	08714729794448	8	120	6	18	C1725	\$560
H74939171081510	39171-08151	08714729794110	8	150	6	18	C1725	\$560
H74939171081810	39171-08181	08714729794684	8	180	6	18	C1725	\$610
H74939171082010	39171-08201	08714729794080	8	200	6	18	C1725	\$610
H74939171090210	39171-09021	08714729794165	9	20	6	18	C1725	\$460
H74939171090310	39171-09031	08714729795162	9	30	6	18	C1725	\$460
H74939171090410	39171-09041	08714729794134	9	40	6	18	C1725	\$460
H74939171090610	39171-09061	08714729795186	9	60	6	18	C1725	\$460
H74939171090810	39171-09081	08714729795209	9	80	6	18	C1725	\$460

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" High Pressure

Mustang™ Balloon Dilatation Catheter (cont)

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939171100210	39171-10021	08714729794073	10	20	6	14	C1725	\$460
H74939171100310	39171-10031	08714729795223	10	30	6	14	C1725	\$460
H74939171100410	39171-10041	08714729793823	10	40	6	14	C1725	\$460
H74939171100610	39171-10061	08714729795230	10	60	6	14	C1725	\$460
H74939171100810	39171-10081	08714729795254	10	80	6	14	C1725	\$460
H74939171120210	39171-12021	08714729794356	12	20	7	14	C1725	\$460
H74939171120310	39171-12031	08714729795261	12	30	7	14	C1725	\$460
H74939171120410	39171-12041	08714729793953	12	40	7	14	C1725	\$460
H74939171120610	39171-12061	08714729795018	12	60	7	14	C1725	\$460
H74939171120810	39171-12081	08714729795285	12	80	7	14	C1725	\$460

Charger™ Balloon Dilatation Catheter

75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939206030270	39206-03027	08714729803942	3	20	5	20	C1725	\$425
H74939206030370	39206-03037	08714729803966	3	30	5	20	C1725	\$425
H74939206030470	39206-03047	08714729803980	3	40	5	20	C1725	\$425
H74939206030670	39206-03067	08714729803997	3	60	5	20	C1725	\$425
H74939206040270	39206-04027	08714729804079	4	20	5	20	C1725	\$425
H74939206040370	39206-04037	08714729804093	4	30	5	20	C1725	\$425
H74939206040470	39206-04047	08714729804116	4	40	5	20	C1725	\$425
H74939206040670	39206-04067	08714729804123	4	60	5	20	C1725	\$425
H74939206040870	39206-04087	08714729804147	4	80	5	20	C1725	\$425
H74939206041070	39206-04107	08714729804024	4	100	5	20	C1725	\$425
H74939206050270	39206-05027	08714729804215	5	20	5	20	C1725	\$425
H74939206050370	39206-05037	08714729804239	5	30	5	20	C1725	\$425
H74939206050470	39206-05047	08714729804253	5	40	5	20	C1725	\$425
H74939206050670	39206-05067	08714729804260	5	60	5	20	C1725	\$425
H74939206050870	39206-05087	08714729804284	5	80	5	20	C1725	\$425
H74939206051070	39206-05107	08714729804161	5	100	5	20	C1725	\$425
H74939206060270	39206-06027	08714729804352	6	20	5	18	C1725	\$425
H74939206060370	39206-06037	08714729804376	6	30	5	18	C1725	\$425
H74939206060470	39206-06047	08714729804390	6	40	5	18	C1725	\$425
H74939206060670	39206-06067	08714729804406	6	60	5	18	C1725	\$425
H74939206060870	39206-06087	08714729804420	6	80	5	18	C1725	\$425
H74939206061070	39206-06107	08714729804307	6	100	5	16	C1725	\$425
H74939206070270	39206-07027	08714729804499	7	20	5	18	C1725	\$425
H74939206070370	39206-07037	08714729804512	7	30	5	18	C1725	\$425
H74939206070470	39206-07047	08714729804536	7	40	5	18	C1725	\$425
H74939206070670	39206-07067	08714729804543	7	60	5	18	C1725	\$425
H74939206070870	39206-07087	08714729804567	7	80	5	18	C1725	\$425
H74939206071070	39206-07107	08714729804444	7	100	6	16	C1725	\$425
H74939206080270	39206-08027	08714729804635	8	20	6	18	C1725	\$425
H74939206080370	39206-08037	08714729804659	8	30	6	18	C1725	\$425
H74939206080470	39206-08047	08714729804673	8	40	6	18	C1725	\$425
H74939206080670	39206-08067	08714729804680	8	60	6	18	C1725	\$425
H74939206080870	39206-08087	08714729804703	8	80	6	18	C1725	\$425
H74939206081070	39206-08107	08714729804581	8	100	6	16	C1725	\$425
H74939206090270	39206-09027	08714729804727	9	20	6	16	C1725	\$425

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" High Pressure

Charger™ Balloon Dilatation Catheter (cont.)

75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939206090370	39206-09037	08714729804734	9	30	6	16	C1725	\$425
H74939206090470	39206-09047	08714729804758	9	40	6	16	C1725	\$425
H74939206090670	39206-09067	08714729804765	9	60	6	16	C1725	\$425
H74939206090870	39206-09087	08714729804772	9	80	6	16	C1725	\$425
H74939206100270	39206-10027	08714729803768	10	20	6	14	C1725	\$425
H74939206100370	39206-10037	08714729803775	10	30	6	14	C1725	\$425
H74939206100470	39206-10047	08714729803799	10	40	6	14	C1725	\$425
H74939206100670	39206-10067	08714729803805	10	60	6	14	C1725	\$425
H74939206100870	39206-10087	08714729803812	10	80	6	14	C1725	\$425
H74939206120270	39206-12027	08714729803836	12	20	7	12	C1725	\$425
H74939206120370	39206-12037	08714729803843	12	30	7	12	C1725	\$425
H74939206120470	39206-12047	08714729803867	12	40	7	12	C1725	\$425
H74939206120670	39206-12067	08714729803874	12	60	7	12	C1725	\$425
H74939206120870	39206-12087	08714729803881	12	80	7	12	C1725	\$425

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939206030210	39206-03021	08714729803935	3	20	5	20	C1725	\$425
H74939206030410	39206-03041	08714729803973	3	40	5	20	C1725	\$425
H74939206030810	39206-03081	08714729804000	3	80	5	20	C1725	\$425
H74939206031010	39206-03101	08714729803898	3	100	5	20	C1725	\$425
H74939206031210	39206-03121	08714729803904	3	120	5	20	C1725	\$520
H74939206031510	39206-03151	08714729803911	3	150	5	20	C1725	\$520
H74939206031810	39206-03181	08714729803928	3	180	5	20	C1725	\$570
H74939206032010	39206-03201	08714729803959	3	200	5	20	C1725	\$570
H74939206040210	39206-04021	08714729804062	4	20	5	20	C1725	\$425
H74939206040410	39206-04041	08714729804109	4	40	5	20	C1725	\$425
H74939206040610	39206-04061	08714729825128	4	60	5	20	C1725	\$425
H74939206040810	39206-04081	08714729804130	4	80	5	20	C1725	\$425
H74939206041010	39206-04101	08714729804017	4	100	5	20	C1725	\$425
H74939206041210	39206-04121	08714729804031	4	120	5	20	C1725	\$520
H74939206041510	39206-04151	08714729804048	4	150	5	20	C1725	\$520
H74939206041810	39206-04181	08714729804055	4	180	5	20	C1725	\$570
H74939206042010	39206-04201	08714729804086	4	200	5	20	C1725	\$570
H74939206050210	39206-05021	08714729804208	5	20	5	20	C1725	\$425
H74939206050410	39206-05041	08714729804246	5	40	5	20	C1725	\$425
H74939206050610	39206-05061	08714729825241	5	60	5	20	C1725	\$425
H74939206050810	39206-05081	08714729804277	5	80	5	20	C1725	\$425
H74939206051010	39206-05101	08714729804154	5	100	5	20	C1725	\$425
H74939206051210	39206-05121	08714729804178	5	120	5	20	C1725	\$520
H74939206051510	39206-05151	08714729804185	5	150	5	20	C1725	\$520
H74939206051810	39206-05181	08714729804192	5	180	5	20	C1725	\$570
H74939206052010	39206-05201	08714729804222	5	200	5	20	C1725	\$570
H74939206060210	39206-06021	08714729804345	6	20	5	18	C1725	\$425
H74939206060410	39206-06041	08714729804383	6	40	5	18	C1725	\$425
H74939206060610	39206-06061	08714729825371	6	60	5	18	C1725	\$425
H74939206060810	39206-06081	08714729804413	6	80	5	18	C1725	\$425
H74939206061010	39206-06101	08714729804291	6	100	5	16	C1725	\$425
H74939206061210	39206-06121	08714729804314	6	120	5	16	C1725	\$520
H74939206061510	39206-06151	08714729804321	6	150	5	16	C1725	\$520
H74939206061810	39206-06181	08714729804338	6	180	6	16	C1725	\$570

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" High Pressure

Charger™ Balloon Dilatation Catheter (cont.)

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939206062010	39206-06201	08714729804369	6	200	6	16	C1725	\$570
H74939206070210	39206-07021	08714729804482	7	20	5	18	C1725	\$425
H74939206070410	39206-07041	08714729804529	7	40	5	18	C1725	\$425
H74939206070610	39206-07061	08714729825500	7	60	5	18	C1725	\$425
H74939206070810	39206-07081	08714729804550	7	80	5	18	C1725	\$425
H74939206071010	39206-07101	08714729804437	7	100	6	16	C1725	\$425
H74939206071210	39206-07121	08714729804451	7	120	6	16	C1725	\$520
H74939206071510	39206-07151	08714729804468	7	150	6	16	C1725	\$520
H74939206071810	39206-07181	08714729804475	7	180	6	16	C1725	\$570
H74939206072010	39206-07201	08714729804505	7	200	6	16	C1725	\$570
H74939206080210	39206-08021	08714729804628	8	20	6	18	C1725	\$425
H74939206080410	39206-08041	08714729804666	8	40	6	18	C1725	\$425
H74939206080810	39206-08081	08714729804697	8	80	6	18	C1725	\$425
H74939206081010	39206-08101	08714729804574	8	100	6	16	C1725	\$425
H74939206081210	39206-08121	08714729804598	8	120	6	16	C1725	\$520
H74939206081510	39206-08151	08714729804604	8	150	6	16	C1725	\$520
H74939206081810	39206-08181	08714729804611	8	180	6	16	C1725	\$570
H74939206082010	39206-08201	08714729804642	8	200	6	16	C1725	\$570
H74939206090210	39206-09021	08714729804710	9	20	6	16	C1725	\$425
H74939206090410	39206-09041	08714729804741	9	40	6	16	C1725	\$425
H74939206100210	39206-10021	08714729803751	10	20	6	14	C1725	\$425
H74939206100410	39206-10041	08714729803782	10	40	6	14	C1725	\$425
H74939206120210	39206-12021	08714729803829	12	20	7	12	C1725	\$425
H74939206120410	39206-12041	08714729803850	12	40	7	12	C1725	\$425

0.035" Ultra-High Pressure

Athletis™ Ultra-High Pressure Balloon

UPN / Order Code	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Length (cm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939347040250	08714729973539	4	20	50	6	40	C1725	\$580
H74939347040450	08714729973553	4	40	50	6	40	C1725	\$580
H74939347040650	08714729973577	4	60	50	6	40	C1725	\$580
H74939347040850	08714729973591	4	80	50	6	40	C1725	\$580
H74939347041050	08714729973614	4	100	50	6	40	C1725	\$580
H74939347050250	08714729973638	5	20	50	6	40	C1725	\$580
H74939347050450	08714729973669	5	40	50	6	40	C1725	\$580
H74939347050650	08714729973683	5	60	50	6	40	C1725	\$580
H74939347050850	08714729973706	5	80	50	6	40	C1725	\$580
H74939347051050	08714729973720	5	100	50	6	40	C1725	\$580
H74939347060250	08714729973751	6	20	50	6	40	C1725	\$580
H74939347060450	08714729973782	6	40	50	6	40	C1725	\$580
H74939347060650	08714729973805	6	60	50	6	40	C1725	\$580
H74939347060850	08714729973829	6	80	50	6	40	C1725	\$580
H74939347061050	08714729973843	6	100	50	6	40	C1725	\$580
H74939347070250	08714729973874	7	20	50	6	34	C1725	\$580
H74939347070450	08714729973904	7	40	50	6	34	C1725	\$580
H74939347070650	08714729973928	7	60	50	6	34	C1725	\$580
H74939347070850	08714729973942	7	80	50	6	34	C1725	\$580
H74939347071050	08714729973966	7	100	50	6	34	C1725	\$580
H74939347080250	08714729973997	8	20	50	6	31	C1725	\$580
H74939347080350	08714729974017	8	30	50	6	31	C1725	\$580
H74939347080450	08714729974048	8	40	50	6	31	C1725	\$580

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" Ultra-High Pressure

Athletis™ Ultra-High Pressure Balloon (cont.)

UPN / Order Code	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Catheter Length (cm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939347080650	08714729974062	8	60	50	6	31	C1725	\$580
H74939347080850	08714729974086	8	80	50	6	31	C1725	\$580
H74939347081050	08714729974109	8	100	50	6	31	C1725	\$580
H74939347090250	08714729974123	9	20	50	7	30	C1725	\$580
H74939347090450	08714729974154	9	40	50	7	30	C1725	\$580
H74939347090850	08714729974178	9	80	50	7	30	C1725	\$580
H74939347040270	08714729973546	4	20	75	6	40	C1725	\$580
H74939347040470	08714729973560	4	40	75	6	40	C1725	\$580
H74939347040670	08714729973584	4	60	75	6	40	C1725	\$580
H74939347040870	08714729973607	4	80	75	6	40	C1725	\$580
H74939347041070	08714729973621	4	100	75	6	40	C1725	\$580
H74939347050270	08714729973645	5	20	75	6	40	C1725	\$580
H74939347050470	08714729973676	5	40	75	6	40	C1725	\$580
H74939347050670	08714729973690	5	60	75	6	40	C1725	\$580
H74939347050870	08714729973713	5	80	75	6	40	C1725	\$580
H74939347051070	08714729973737	5	100	75	6	40	C1725	\$580
H74939347060270	08714729973768	6	20	75	6	40	C1725	\$580
H74939347060470	08714729973799	6	40	75	6	40	C1725	\$580
H74939347060670	08714729973812	6	60	75	6	40	C1725	\$580
H74939347060870	08714729973836	6	80	75	6	40	C1725	\$580
H74939347061070	08714729973850	6	100	75	6	40	C1725	\$580
H74939347070270	08714729973881	7	20	75	6	34	C1725	\$580
H74939347070470	08714729973911	7	40	75	6	34	C1725	\$580
H74939347070670	08714729973935	7	60	75	6	34	C1725	\$580
H74939347070870	08714729973959	7	80	75	6	34	C1725	\$580
H74939347071070	08714729973973	7	100	75	6	34	C1725	\$580
H74939347080270	08714729974000	8	20	75	6	31	C1725	\$580
H74939347080370	08714729974024	8	30	75	6	31	C1725	\$580
H74939347080470	08714729974055	8	40	75	6	31	C1725	\$580
H74939347080670	08714729974079	8	60	75	6	31	C1725	\$580
H74939347080870	08714729974093	8	80	75	6	31	C1725	\$580
H74939347081070	08714729974116	8	100	75	6	31	C1725	\$580
H74939347090270	08714729974130	9	20	75	7	30	C1725	\$580
H74939347090470	08714729974161	9	40	75	7	30	C1725	\$580
H74939347090870	08714729974185	9	80	75	7	30	C1725	\$580
H74939347100270	08714729974192	10	20	75	7	30	C1725	\$580
H74939347100470	08714729974215	10	40	75	7	30	C1725	\$580
H74939347100870	08714729974222	10	80	75	7	30	C1725	\$580
H74939347120270	08714729974246	12	20	75	8	20	C1725	\$580
H74939347120470	08714729974260	12	40	75	8	20	C1725	\$580
H74939347120670	08714729974284	12	60	75	8	20	C1725	\$580
H74939347050410	08714729973652	5	40	135	6	40	C1725	\$580
H74939347060210	08714729973744	6	20	135	6	40	C1725	\$580
H74939347060410	08714729973775	6	40	135	6	40	C1725	\$580
H74939347070210	08714729973867	7	20	135	6	34	C1725	\$580
H74939347070410	08714729973898	7	40	135	6	34	C1725	\$580
H74939347080210	08714729973980	8	20	135	6	31	C1725	\$580
H74939347080410	08714729974031	8	40	135	6	31	C1725	\$580
H74939347090410	08714729974147	9	40	135	7	30	C1725	\$580
H74939347100410	08714729974208	10	40	135	7	30	C1725	\$580
H74939347120210	08714729974239	12	20	135	8	20	C1725	\$580
H74939347120410	08714729974253	12	40	135	8	20	C1725	\$580
H74939347120610	08714729974277	12	60	135	8	20	C1725	\$580

ANGIOPLASTY BALLOON DILATATION CATHETERS

0.035" Large Vessel

XXL™ Balloon Dilatation Catheter

75 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
M001145080	14-508	08714729162865	12	20	7	8	C1725	\$525
M001145100	14-510	08714729096146	12	40	7	8	C1725	\$525
M001145300	14-530	08714729162919	12	60	7	8	C1725	\$525
M001145120	14-512	08714729162926	14	20	7	8	C1725	\$525
M001145140	14-514	08714729190462	14	40	7	8	C1725	\$525
M001145320	14-532	08714729190486	14	60	7	8	C1725	\$525
M001145500	14-550	08714729163039	16	20	N/A	5	C1725	\$580
M001145520	14-552	08714729163060	16	40	N/A	5	C1725	\$580
M001145540	14-554	08714729163114	16	60	N/A	5	C1725	\$580
M001145550	14-555	08714729163121	18	20	N/A	5	C1725	\$580
M001145570	14-557	08714729163169	18	40	N/A	5	C1725	\$580
M001145590	14-559	08714729163190	18	60	N/A	5	C1725	\$580

120 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
M001145090	14-509	08714729190455	12	20	7	8	C1725	\$525
M001145110	14-511	08714729162902	12	40	7	8	C1725	\$525
M001145130	14-513	08714729162964	14	20	7	8	C1725	\$525
M001145150	14-515	08714729190479	14	40	7	8	C1725	\$525
M001145510	14-551	08714729163053	16	20	N/A	5	C1725	\$580
M001145530	14-553	08714729163084	16	40	N/A	5	C1725	\$580
M001145560	14-556	08714729163145	18	20	N/A	5	C1725	\$580
M001145580	14-558	08714729163176	18	40	N/A	5	C1725	\$580

Occlusion Balloons

Equalizer™ Occlusion Balloon Catheter

65 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Inflated Balloon Diameter (mm)	Recommended Inflation Volume mL (cc)	Introducer Sheath (F)	C-Code	Price
M001171050	17-105	08714729402350	20	4.8	14	C2628	\$360
M001171070	17-107	08714729402374	27	10	14	C2628	\$360
M001171090	17-109	08714729402398	33	19	16	C2628	\$360
M001171110	17-111	08714729402411	40	34.5	N/A	C2628	\$360

100 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Inflated Balloon Diameter (mm)	Recommended Inflation Volume mL (cc)	Introducer Sheath (F)	C-Code	Price
M001171060	17-106	08714729402367	20	4.8	14	C2628	\$360
M001171080	17-108	08714729402381	27	10	14	C2628	\$360
M001171100	17-110	08714729402404	33	19	16	C2628	\$360
M001171120	17-112	08714729402428	40	34.5	N/A	C2628	\$360

Peripheral Cutting Balloon™ Microsurgical Dilatation Device

50 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	Recommended Guidewire (in)	C-Code	Price
M001PCB5020500	PCB502050	08714729813705	5	20	6	10	0.018	C1725	\$1,325
M001PCB6020500	PCB602050	08714729813736	6	20	6	10	0.018	C1725	\$1,325
M001PCB7020500	PCB702050	08714729813767	7	20	7	10	0.018	C1725	\$1,325
M001PCB8020500	PCB802050	08714729813798	8	20	7	10	0.018	C1725	\$1,325

90 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	Recommended Guidewire (in)	C-Code	Price
M001PCB5020900	PCB502090	08714729813712	5	20	6	10	0.018	C1725	\$1,325
M001PCB6020900	PCB602090	08714729813743	6	20	6	10	0.018	C1725	\$1,325
M001PCB7020900	PCB702090	08714729813774	7	20	7	10	0.018	C1725	\$1,325
M001PCB8020900	PCB802090	08714729813804	8	20	7	10	0.018	C1725	\$1,325

135 cm Usable Length

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	Recommended Guidewire (in)	C-Code	Price
M001PCB50201350	PCB5020135	08714729813699	5	20	6	10	0.018	C1725	\$1,325
M001PCB60201350	PCB6020135	08714729813729	6	20	6	10	0.018	C1725	\$1,325
M001PCB70201350	PCB7020135	08714729813750	7	20	7	10	0.018	C1725	\$1,325
M001PCB80201350	PCB8020135	08714729813781	8	20	7	10	0.018	C1725	\$1,325

Peripheral Rotablator™ Rotational Atherectomy System

Peripheral Rotalink™ Plus Burr Catheter

135cm Catheter Length, 0.009" OTW, 4.5F (0.058") Protective Sheath

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Description	C-Code	Price
M0011401250	140-125	08714729857808	5	1.25mm Peripheral RotaLink Plus	C1724	\$3,306
M0011401500	140-150	08714729838906	5	1.50mm Peripheral RotaLink Plus	C1724	\$3,306
M0011401750	140-175	08714729838913	6	1.75mm Peripheral RotaLink Plus	C1724	\$3,306
M0011402000	140-200	08714729838920	7	2.00mm Peripheral RotaLink Plus	C1724	\$3,306
M0011402250	140-225	08714729838944	7	2.25mm Peripheral RotaLink Plus	C1724	\$3,306
M0011402500	140-250	08714729838968	8	2.50mm Peripheral RotaLink Plus	C1724	\$3,306

Peripheral ROTAPRO™ Burr Catheter

135cm Catheter Length, 0.009" OTW, 4.5F (0.058") Protective Sheath

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Description	C-Code	Price
H749395141250	n/a	00191506007119	5	1.25mm Peripheral ROTAPRO Burr Catheter	C1724	\$3,495
H749395141500	n/a	00191506007126	5	1.50mm Peripheral ROTAPRO Burr Catheter	C1724	\$3,495
H749395141750	n/a	00191506007133	6	1.75mm Peripheral ROTAPRO Burr Catheter	C1724	\$3,495
H749395142000	n/a	00191506007140	7	2.00mm Peripheral ROTAPRO Burr Catheter	C1724	\$3,495
H749395142250	n/a	00191506007164	7	2.25mm Peripheral ROTAPRO Burr Catheter	C1724	\$3,495
H749395142500	n/a	00191506007188	8	2.50mm Peripheral ROTAPRO Burr Catheter	C1724	\$3,495

Peripheral RotaWire™ Guide Wire and WireClip™ Torquer

UPN / Order Code	Catalog No.	GTIN	Description	Quantity	C-Code	Price
M0011373321	137-3321	08714729838890	Extra Support Peripheral RotaWire	5	C1769	\$1,379

ROTAPRO™ System Console

UPN / Order Code	Catalog No.	GTIN	Description	C-Code	Price
H74939309010	n/a	08714729975076	ROTAPRO Console Kit Assy UL/CSA	N/A	\$20,000
H80215901011	15901-01	08714729263371	HTI Regulator Kit - US	N/A	\$570

Rotablator Console

UPN / Order Code	Catalog No.	GTIN	Description	C-Code	Price
H802220200391	22020-039	08714729353300	Rotablator Console Kit (includes Console, foot pedal and 20ft. Air Supply Hose)	N/A	\$19,995
H802216000031	21600-003	08714729076964	20ft. Air Supply Hose - Replacement	N/A	\$431
H802224360021	22436-002	08714729353317	Dynaglide™ Foot Pedal - Replacement	N/A	\$1,133

Accessories

UPN / Order Code	Catalog No.	GTIN	Description	C-Code	Price
H74939369010	n/a	0874729974307	ROTAGO™ Cart and Infusion Stand	N/A	\$7,995
M0013948400172	n/a	00191506005580	Peripheral RotaGlide™ Lubricant (Box of 10)	N/A	\$1,240
H802221960032	22196-003	08714729195535	WireClip Torquer (Box of 5)	N/A	\$155

Jetstream™ Console

UPN / Order Code	Catalog No.	GTIN	Description	C-Code	Price
050599-001	PVCN100	08714729890430	Jetstream Atherectomy System	NA	\$29,995

Jetstream™ XC Atherectomy Catheters

UPN / Order Code	Catalog No.	GTIN	Catheter Length (cm)	Min. Introducer Size (F)	Max. Guidewire Diameter (in)	Tip Diameter (mm)	Target Therapy Speed (rpm)	C-Code	Price
112266-001	PV41340	08714729889922	120	7	0.014	2.4 / 3.4	70K	C1724	\$5,890
112264-001	PV31300	08714729889892	135	7	0.014	2.1 / 3.0	70K	C1724	\$5,890

Jetstream™ SC Atherectomy Catheters

UPN / Order Code	Catalog No.	GTIN	Catheter Length (cm)	Min. Introducer Size (F)	Max. Guidewire Diameter (in)	Tip Diameter (mm)	Target Therapy Speed (rpm)	C-Code	Price
112262-001	PV3118F	08714729889861	145	7	0.014	1.85	73K	C1724	\$5,890
112260-001	PV3116F	08714729889830	145	7	0.014	1.6	73K	C1724	\$5,890

Accessories

Thruway™ 0.014 Guidewire, 300 cm Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492971	49-297	08714729717188	0.014	300	Short	Straight	1	C1769	\$210

All-Purpose Drainage Catheters

Flexima™ Material Catheters

UPN / Order Code	Catalog No.	GTIN	APD	Size (F)	Length (cm)	Durometer K (Kit)	C-Code	Price
M001271240	27-124	08714729323044	APD Loop	10	25	Firm	C1729	\$140
M001271250	27-125	08714729271369	APD	10	20	Firm	C1729	\$140
M001271260	27-126	08714729271376	APD	6	20	Firm	C1729	\$140
M001271270	27-127	08714729271383	APD	8	20	Firm	C1729	\$140
M001271280	27-128	08714729323051	APD	6	20	Reg	C1729	\$140
M001271290	27-129	08714729323068	APD	8	20	Reg	C1729	\$140
M001271300	27-130	08714729323075	APD	10	20	Reg	C1729	\$140
M001271310	27-131	08714729323082	APD Loop	6	25	Firm	C1729	\$140
M001271320	27-132	08714729191759	APD Loop	8	25	Firm	C1729	\$140
M001271330	27-133	08714729323099	APD Loop	6	25	Reg	C1729	\$140
M001271340	27-134	08714729323105	APD Loop	8	25	Reg	C1729	\$140
M001271350	27-135	08714729323112	APD Loop	10	25	Reg	C1729	\$140
M001271360	27-136	08714729323129	APD	12	25	Reg	C1729	\$140
M001271370	27-137	08714729323136	APD	14	25	Reg	C1729	\$140
M001271380	27-138	08714729323143	APD Loop	12	25	Reg	C1729	\$140
M001271390	27-139	08714729323150	APD Loop	14	25	Reg	C1729	\$140
M001271460	27-146	08714729897170	Quickstick Loop	8	25	Reg	C1729	\$120
M001271470	27-147	08714729821250	Quickstick Loop	10	25	Reg	C1729	\$120

All catheters packaged with a metal stiffening cannula and trocar; in addition, all catheters except 6F contain a flexible stiffening cannula.

Flexima™ Material Kits

UPN / Order Code	Catalog No.	GTIN	APD	Size (F)	Length (cm)	Durometer K (Kit)	C-Code	Price
M001271960	27-196	08714729323471	APD Loop	8	25	Reg-K	C1729	\$255
M001271970	27-197	08714729323488	APD Loop	10	25	Reg-K	C1729	\$255
M001271980	27-198	08714729323495	APD Loop	8	25	Firm-K	C1729	\$255
M001271990	27-199	08714729323501	APD Loop	10	25	Firm-K	C1729	\$255

All catheters packaged with a trocar as well as both metal and flexible stiffening cannulas.

Kits contain Accustick™ II Introducer Set with 0.038" and 0.018" wires, Percufix™ catheter cuff kit, dilators, and connecting tube with stop cock.

vanSonnenberg™ Sump Catheters

Flexima™ Material Catheters

UPN / Order Code	Catalog No.	GTIN	vanSonnenberg Sump	Size (F)	Length (cm)	Durometer K (Kit)	C-Code	Price
M001273000	27-300	08714729137733	VSC	12	30	Reg	C1729	\$195
M001273010	27-301	08714729138105	VSC	14	30	Reg	C1729	\$195
M001273030	27-303	08714729138112	VSC Loop	12	35	Reg	C1729	\$195
M001273040	27-304	08714729138129	VSC Loop	14	35	Reg	C1729	\$195

All catheters packaged with a metal stiffening cannula, trocar and microbial air filter.

Flexima™ Material Kits

UPN / Order Code	Catalog No.	GTIN	vanSonnenberg Sump	Size (F)	Length (cm)	Durometer K (Kit)	C-Code	Price
M001273060	27-306	08714729137658	VSDC	12	30	Reg-K	C1729	\$315
M001273070	27-307	08714729137665	VSDC	14	30	Reg-K	C1729	\$315
M001273090	27-309	08714729137740	VSDC Loop	12	35	Reg-K	C1729	\$315
M001273100	27-310	08714729149118	VSDC Loop	14	35	Reg-K	C1729	\$315

All-Purpose Drainage Catheters

Percuflex™ Material Glidex™ Coated Catheters

UPN / Order Code	Catalog No.	GTIN	APD	Size (F)	Length (cm)	C-Code	Price
M001245071	24-507	08714729115243	APDG Loop	8	25	C1729	\$775
M001245081	24-508	08714729115250	APDG Loop	10	25	C1729	\$775

All catheters packaged with both metal stiffening cannula and trocar.

vanSonnenberg™ Chest Drainage Catheters

Percuflex™ Material Glidex™ Coated Catheters

UPN / Order Code	Catalog No.	GTIN	VCD	Size (F)	Length (cm)	C-Code	Price
M001245130	24-513	08714729271352	VCDG Loop	12	25	C1729	\$155
M001245140	24-514	08714729048541	VCDG Loop	14	25	C1729	\$155

All catheters packaged with metal stiffening cannula

Biliary Drainage

Biliary Drainage Catheters

Flexima™ Material Catheters

UPN / Order Code	Catalog No.	GTIN	Biliary	Size (F)	Length (cm) K (Kit)	C-Code	Price
M001271540	27-154	08714729323167	Biliary	8	35 Firm	C1729	\$145
M001271550	27-155	08714729323174	Biliary	10	35 Firm	C1729	\$145
M001271560	27-156	08714729323181	Biliary	8	35 Reg	C1729	\$145
M001271570	27-157	08714729323198	Biliary	10	35 Reg	C1729	\$145
M001271580	27-158	08714729323204	Biliary	12	35 Reg	C1729	\$145
M001271590	27-159	08714729323211	Biliary	14	35 Reg	C1729	\$145

All catheters packaged with a luer cap with septum as well as both metal and flexible stiffening cannulas.

Flexima™ Material Kits

UPN / Order Code	Catalog No.	GTIN	Biliary	Size (F)	Length (cm) K (Kit)	C-Code	Price
M001271630	27-163	08714729323259	Biliary	8	35 Firm-K	C1729	\$260
M001271640	27-164	08714729323266	Biliary	10	35 Firm-K	C1729	\$260
M001271650	27-165	08714729323273	Biliary	8	35 Reg-K	C1729	\$260
M001271660	27-166	08714729323280	Biliary	10	35 Reg-K	C1729	\$260

All catheters packaged with a luer cap with septum as well as both metal and flexible stiffening cannulas.

Biliary Drainage Catheters with RO Marker

Flexima™ Material Catheters

UPN / Order Code	Catalog No.	GTIN	Biliary	Size (F)	Length (cm)	C-Code	Price
M001272600	27-260	08714729454953	Biliary RO	8	35 Reg	C1729	\$150
M001272610	27-261	08714729454960	Biliary RO	10	35 Reg	C1729	\$150
M001272620	27-262	08714729454977	Biliary RO	12	35 Reg	C1729	\$150
M001272630	27-263	08714729454984	Biliary RO	14	35 Reg	C1729	\$150
M001272640	27-264	08714729454991	Biliary RO	8	35 Firm	C1729	\$150
M001272650	27-265	08714729455004	Biliary RO	10	35 Firm	C1729	\$150

All catheters packaged with a luer cap with septum as well as both metal and flexible stiffening cannulas.

Flexima™ Material (Nephrostomy Catheter)

Catheters

UPN / Order Code	Catalog No.	GTIN	Nephrostomy	Size (F)	Length (cm) K (Kit)	C-Code	Price
M001271790	27-179	08714729323303	Neph	8	25 Reg	C1729	\$145
M001271800	27-180	08714729323310	Neph	10	25 Reg	C1729	\$145
M001271810	27-181	08714729323327	Neph	12	25 Reg	C1729	\$145
M001271820	27-182	08714729323334	Neph	14	25 Reg	C1729	\$145
M001271770	27-177	08714729195429	Neph	8	25 Firm	C1729	\$145
M001271780	27-178	08714729323297	Neph	10	25 Firm	C1729	\$145

All catheters packaged with both metal and flexible stiffening cannula.

Flexima™ Material (Nephrostomy Catheter)

Kits

UPN / Order Code	Catalog No.	GTIN	Nephrostomy	Size (F)	Length (cm) K (Kit)	C-Code	Price
M001271860	27-186	08714729323372	Neph	8	25 Reg - K	C1729	\$260
M001271870	27-187	08714729323389	Neph	10	25 Reg - K	C1729	\$260
M001271880	27-188	08714729323396	Neph	12	25 Reg - K	C1729	\$260
M001271900	27-190	08714729323419	Neph	8	25 Firm - K	C1729	\$260

Kits contain Accustick II Introducer Set with 0.038" and 0.018" wires, Percufix™ catheter cuff kit, dilators and connecting tube with stop cock

Flexima™ Material (Ureteral Stents)

Stents

UPN / Order Code	Catalog No.	GTIN	Ureteral Stent	Size (F)	Length (cm) Between Pigtails - K (Kit)	C-Code	Price
M001274000	27-400	08714729323754	US	8	20	C2617	\$195
M001274010	27-401	08714729323761	US	8	22	C2617	\$195
M001274020	27-402	08714729323778	US	8	24	C2617	\$195
M001274030	27-403	08714729323785	US	8	26	C2617	\$195
M001274040	27-404	08714729323792	US	8	28	C2617	\$195
M001274050	27-405	08714729323808	US	10	20	C2617	\$195
M001274060	27-406	08714729323815	US	10	22	C2617	\$195
M001274070	27-407	08714729323822	US	10	24	C2617	\$195
M001274080	27-408	08714729323839	US	10	26	C2617	\$195
M001274090	27-409	08714729323846	US	10	28	C2617	\$195

All stents packaged with flexible stiffening cannula and stabilizer

Kits

UPN / Order Code	Catalog No.	GTIN	Ureteral Stent	Size (F)	Length (cm) Between Pigtails - K (Kit)	C-Code	Price
M001274110	27-411	08714729323853	US	6	22-K	C2617	\$315
M001274120	27-412	08714729323860	US	6	24-K	C2617	\$315
M001274130	27-413	08714729323877	US	6	26-K	C2617	\$315
M001274140	27-414	08714729323884	US	6	28-K	C2617	\$315
M001274150	27-415	08714729323891	US	8	20-K	C2617	\$315
M001274160	27-416	08714729323907	US	8	22-K	C2617	\$315
M001274170	27-417	08714729323914	US	8	24-K	C2617	\$315
M001274180	27-418	08714729323921	US	8	26-K	C2617	\$315
M001274190	27-419	08714729323938	US	8	28-K	C2617	\$315

Kits contain an 0.038" guidewire, stabilizer, and flexible stiffening cannula

Percuflex™ Material (VTC Nephrostomy Catheter)

Non-Coated Catheter

UPN / Order Code	Catalog No.	GTIN	VTC	Catheter Size (F)	Length (cm) K (Kit)	C-Code	Price
M001225140	22-514	08714729031178	VTC Firm	8	25	C1729	\$130

All catheters packaged with both metal and flexible stiffening cannula.

Glidex™ Coated Catheters

UPN / Order Code	Catalog No.	GTIN	VTC	Catheter Size (F)	Length (cm) K (Kit)	C-Code	Price
M001245260	24-526	08714729032052	VTCTG	8	25	C1729	\$155
M001245270	24-527	08714729032069	VTCTG	10	25	C1729	\$155
M001245280	24-528	08714729048558	VTCTG	12	25	C1729	\$155
M001245290	24-529	08714729032076	VTCTG	14	25	C1729	\$155

All catheters packaged with both metal and flexible stiffening cannula.

Glidex™ Coated Kits

UPN / Order Code	Catalog No.	GTIN	VTC	Catheter Size (F)	Length (cm) K (Kit)	C-Code	Price
M001245300	24-530	08714729048565	VTCTG	8	25-K	C1729	\$270
M001245310	24-531	08714729048572	VTCTG	10	25-K	C1729	\$270
M001245320	24-532	08714729048589	VTCTG	12	25-K	C1729	\$270

All catheters packaged with both metal and flexible stiffening cannula.

Kits contain Accustick II Introducer Set with 0.038" and 0.018" wires, Percuflex™ catheter cuff kit, dilators and connecting tube with stop cock

Percuflex™ Material (Amplatz Anchor™ Catheters with Locking Malecot)

Non-Coated Catheter

UPN / Order Code	Catalog No.	GTIN	Amplatz Anchor	Size (F)	Length (cm)	C-Code	Price
M001221000	22-100	08714729030638	AA	8	30	C1729	\$145
M001221010	22-101	08714729030645	AA	10	30	C1729	\$145
M001221020	22-102	08714729030652	AA	12	30	C1729	\$145

All catheters packaged with both metal and flexible stiffening cannula.

Percuflex™ Material

Non-Coated Stents (Nephroureteral Stents)

UPN / Order Code	Catalog No.	GTIN	Nephroureteral Stent	Size (F)	Length (cm) Between Pigtails - K (Kit)	C-Code	Price
M001221360	22-136	08714729114079	IE	8	22	C2617	\$180
M001221370	22-137	08714729114086	IE	8	24	C2617	\$180
M001221380	22-138	08714729114093	IE	8	26	C2617	\$180
M001221390	22-139	08714729114109	IE	8	28	C2617	\$180
M001221400	22-140	08714729114116	IE	10	22	C2617	\$180
M001221410	22-141	08714729114123	IE	10	24	C2617	\$180
M001221420	22-142	08714729114130	IE	10	26	C2617	\$180
M001221430	22-143	08714729114147	IE	10	28	C2617	\$180

All stents packaged with metal stiffening cannula

Glidex™ Coated Stents (Ureteral Stents)

UPN / Order Code	Catalog No.	GTIN	Ureteral Stent	Size (F)	Length (cm) Between Pigtails - K (Kit)	C-Code	Price
M001245490	24-549	08714729048602	USG	8	20	C2617	\$240
M001245500	24-550	08714729178613	USG	8	22	C2617	\$240
M001245510	24-551	08714729178620	USG	8	24	C2617	\$240
M001245520	24-552	08714729170556	USG	8	26	C2617	\$240
M001245530	24-553	08714729048633	USG	8	28	C2617	\$240
M001245570	24-557	08714729069508	USG	10	20	C2617	\$240
M001245580	24-558	08714729069515	USG	10	22	C2617	\$240
M001245590	24-559	08714729069522	USG	10	24	C2617	\$240
M001245600	24-560	08714729069539	USG	10	26	C2617	\$240
M001245610	24-561	08714729069546	USG	10	28	C2617	\$240

All stents packaged with flexible stiffening cannula and stabilizer

Percuflex™ Material
Glidex™ Coated Kits (Ureteral Stents)

UPN / Order Code	Catalog No.	GTIN	Ureteral Stent	Size (F)	Length (cm) Between Pigtails - K (Kit)	C-Code	Price
M001245720	24-572	08714729048640	USGK	6	20-K	C2617	\$360
M001245730	24-573	08714729048657	USGK	6	22-K	C2617	\$360
M001245740	24-574	08714729048664	USGK	6	24-K	C2617	\$360
M001245750	24-575	08714729048671	USGK	6	26-K	C2617	\$360
M001245760	24-576	08714729048688	USGK	6	28-K	C2617	\$360
M001245620	24-562	08714729032090	USGK	8	20-K	C2617	\$360
M001245630	24-563	08714729032106	USGK	8	22-K	C2617	\$360
M001245640	24-564	08714729032113	USGK	8	24-K	C2617	\$360
M001245650	24-565	08714729032120	USGK	8	26-K	C2617	\$360
M001245660	24-566	08714729032137	USGK	8	28-K	C2617	\$360

Kits contain an 0.038" guidewire, stabilizer, and flexible stiffening cannula.

Stents

Eluvia™ Drug-Eluting Vascular Stent System
130 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code*	Price
H74939294600410	N/A	8714729876571	6	40	6	130	1,2	C1874	\$4,400
H74939294600610	N/A	8714729876588	6	60	6	130	1,2	C1874	\$4,400
H74939294600810	N/A	8714729876595	6	80	6	130	1,2	C1874	\$4,730
H74939294601010	N/A	8714729876601	6	100	6	130	1,2	C1874	\$4,730
H74939294601210	N/A	8714729876618	6	120	6	130	1,2	C1874	\$4,730
H74939294601510	N/A	8714729876625	6	150	6	130	1,2	C1874	\$5,060
H74939294700410	N/A	8714729876694	7	40	6	130	1,2	C1874	\$4,400
H74939294700610	N/A	8714729876700	7	60	6	130	1,2	C1874	\$4,400
H74939294700810	N/A	8714729876717	7	80	6	130	1,2	C1874	\$4,730
H74939294701010	N/A	8714729876724	7	100	6	130	1,2	C1874	\$4,730
H74939294701210	N/A	8714729876731	7	120	6	130	1,2	C1874	\$4,730
H74939294701510	N/A	8714729876748	7	150	6	130	1,2	C1874	\$5,060

Indications: 1 = Superficial Femoral Artery , 2 = Proximal Popliteal Artery

*As Eluvia is a new and novel technology, Boston Scientific believes that currently available C-Codes do not accurately describe it. However, the closest applicable C-Code to report the use of Eluvia is C1874, defined as "Stent, coated/covered, with delivery system."

Balloons

Ranger™ Drug-Coated Balloon
Low profile .018"

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939419400410	N/A	08714729975762	4	40	135	5	14	C2623	\$2,270
H74939419400480	N/A	08714729975779	4	40	80	5	14	C2623	\$2,270
H74939419400610	N/A	08714729975786	4	60	135	5	14	C2623	\$2,270
H74939419400680	N/A	08714729975793	4	60	80	5	14	C2623	\$2,270
H74939419400810	N/A	08714729975809	4	80	135	5	14	C2623	\$2,270
H74939419400880	N/A	08714729975816	4	80	80	5	14	C2623	\$2,270
H74939419401010	N/A	08714729975823	4	100	135	5	14	C2623	\$2,270
H74939419401080	N/A	08714729975830	4	100	80	5	14	C2623	\$2,270
H74939419401210	N/A	08714729975847	4	120	150	5	14	C2623	\$2,520
H74939419401510	N/A	08714729975861	4	150	150	5	14	C2623	\$2,775
H74939419401590	N/A	08714729975878	4	150	90	5	14	C2623	\$2,775
H74939419402010	N/A	08714729975885	4	200	150	5	14	C2623	\$3,155
H74939419402090	N/A	08714729975892	4	200	90	5	14	C2623	\$3,155
H74939419500410	N/A	08714729975922	5	40	135	5	14	C2623	\$2,270
H74939419500480	N/A	08714729975939	5	40	80	5	14	C2623	\$2,270
H74939419500610	N/A	08714729975946	5	60	135	5	14	C2623	\$2,270
H74939419500680	N/A	08714729975953	5	60	80	5	14	C2623	\$2,270
H74939419500810	N/A	08714729975960	5	80	135	5	14	C2623	\$2,270
H74939419500880	N/A	08714729975977	5	80	80	5	14	C2623	\$2,270
H74939419501010	N/A	08714729975984	5	100	135	5	14	C2623	\$2,270
H74939419501080	N/A	08714729975991	5	100	80	5	14	C2623	\$2,270
H74939419501210	N/A	08714729976004	5	120	150	5	14	C2623	\$2,520
H74939419501510	N/A	08714729976028	5	150	150	5	14	C2623	\$2,775
H74939419501590	N/A	08714729976035	5	150	90	5	14	C2623	\$2,775
H74939419502010	N/A	08714729976042	5	200	150	5	14	C2623	\$3,155
H74939419502090	N/A	08714729976059	5	200	90	5	14	C2623	\$3,155
H74939419600410	N/A	08714729976080	6	40	135	5	14	C2623	\$2,270
H74939419600480	N/A	08714729976097	6	40	80	5	14	C2623	\$2,270
H74939419600610	N/A	08714729976103	6	60	135	5	14	C2623	\$2,270
H74939419600680	N/A	08714729976110	6	60	80	5	14	C2623	\$2,270
H74939419600810	N/A	08714729976127	6	80	135	5	14	C2623	\$2,270
H74939419600880	N/A	08714729976134	6	80	80	5	14	C2623	\$2,270

Balloons

Ranger™ Drug-Coated Balloon (cont.)
Low profile .018"

UPN / Order Code	Catalog No.	GTIN	Balloon Diameter (mm)	Balloon Length (mm)	Shaft Length (cm)	Introducer Sheath (F)	Rated Burst Pressure (ATM)	C-Code	Price
H74939419601010	N/A	08714729976141	6	100	135	5	14	C2623	\$2,270
H74939419601080	N/A	08714729976158	6	100	80	5	14	C2623	\$2,270
H74939419601210	N/A	08714729976165	6	120	150	5	14	C2623	\$2,520
H74939419601510	N/A	08714729976189	6	150	150	5	14	C2623	\$2,775
H74939419601590	N/A	08714729976196	6	150	90	5	14	C2623	\$2,775
H74939419602010	N/A	08714729976202	6	200	150	5	14	C2623	\$3,155
H74939419602090	N/A	08714729976219	6	200	90	5	14	C2623	\$3,155
H74939419700410	N/A	08714729976240	7	40	135	6	14	C2623	\$2,270
H74939419700480	N/A	08714729976257	7	40	80	6	14	C2623	\$2,270
H74939419700610	N/A	08714729976264	7	60	135	6	14	C2623	\$2,270
H74939419700680	N/A	08714729976271	7	60	80	6	14	C2623	\$2,270
H74939419700810	N/A	08714729976288	7	80	135	6	14	C2623	\$2,270
H74939419700880	N/A	08714729976295	7	80	80	6	14	C2623	\$2,270
H74939419701010	N/A	08714729976301	7	100	135	6	14	C2623	\$2,270
H74939419701080	N/A	08714729976318	7	100	80	6	14	C2623	\$2,270
H74939419701210	N/A	08714729976325	7	120	150	6	14	C2623	\$2,520
H74939419701510	N/A	08714729976349	7	150	150	6	14	C2623	\$2,775
H74939419702010	N/A	08714729976363	7	200	150	6	14	C2623	\$3,155

EMBOLIC BEADS & PARTICLES

MICROSPHERES

LC Bead M1™

UPN / Order Code	GTIN	Particle Size (Microns)	Box (Units)	Packaging	C-Code	Price
VE020GS	05060116920222	70-150	1	Vial	C1889	\$3,300

LC Bead LUMI™ M0

UPN / Order Code	GTIN	Particle Size (Microns)	Box (Units)	Packaging	C-Code	Price
RO2L000	05060116921229	40-90	1	Vial	C1889	\$3,850

LC Bead LUMI™ M1

UPN / Order Code	GTIN	Particle Size (Microns)	Box (Units)	Packaging	C-Code	Price
RO2L001	05060116921151	70-150	1	Vial	C1889	\$3,850

LC Bead™

LC Bead™ 100-300

UPN / Order Code	GTIN	Particle Size (Microns)	Box (Units)	Packaging	C-Code	Price
VE220GS	05060116920185	100-300	1	Vial	C1889	\$2,310

LC Bead™ 300-500

UPN / Order Code	GTIN	Particle Size (Microns)	Box (Units)	Packaging	C-Code	Price
VE420GS	05060116920192	300-500	1	Vial	C1889	\$2,310

Particles

Contour™ - PVA Particles

Embolization Particles

UPN / Order Code	GTIN	Particle Size (Microns)	Box (Units)	Packaging	C-Code	Price
M0017600121	08714729159667	45 - 150	2	Vial	C1889	\$455
M0017600151	08714729159674	45 - 150	5	Vial	C1889	\$1,090
M0017600221	08714729159681	150 - 250	2	Vial	C1889	\$455
M0017600251	08714729159698	150 - 250	5	Vial	C1889	\$1,090
M0017600321	08714729159704	250 - 355	2	Vial	C1889	\$455
M0017600351	08714729159711	250 - 355	5	Vial	C1889	\$1,090
M0017600421	08714729159728	355 - 500	2	Vial	C1889	\$455
M0017600451	08714729159735	355 - 500	5	Vial	C1889	\$1,090
M0017600621	08714729159742	500 - 710	2	Vial	C1889	\$455
M0017600651	08714729159759	500 - 710	5	Vial	C1889	\$1,090
M0017600821	08714729159766	710 - 1000	2	Vial	C1889	\$455
M0017600851	08714729159773	710 - 1000	5	Vial	C1889	\$1,090
M0017601151	08714729159797	1000 - 1180	5	Vial	C1889	\$1,090

Mach1™ Peripheral Guide Catheter**Mach1™ Peripheral Guide Catheter 6 F**

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Length (cm)	I.D. (in)	Tip Shape	C-Code	Price
M001196000	19-600	08714729403142	6	55	0.070	HS	C1887	\$185
M001196040	19-604	08714729403173	6	55	0.070	RE S	C1887	\$185
M001196080	19-608	08714729403203	6	55	0.070	RE SS	C1887	\$185
M001196120	19-612	08714729403234	6	55	0.070	RE L	C1887	\$185
M001196160	19-616	08714729403265	6	55	0.070	C1	C1887	\$185
M001196200	19-620	08714729403296	6	55	0.070	C2	C1887	\$185
M001196230	19-623	08714729403326	6	55	0.070	MP	C1887	\$185
M001196260	19-626	08714729403357	6	90	0.070	MP	C1887	\$185
M001196300	19-630	08714729403388	6	55	0.070	MP2	C1887	\$185
M001196330	19-633	08714729403418	6	90	0.070	MP2	C1887	\$185
M001196380	19-638	08714729403456	6	90	0.070	40°	C1887	\$185
M001196450	19-645	08714729403500	6	55	0.070	ST	C1887	\$185
M001196480	19-648	08714729403531	6	90	0.070	ST	C1887	\$185
M001196530	19-653	08714729403562	6	55	0.070	CROSSOVER 1	C1887	\$185
M001196560	19-656	08714729403593	6	55	0.070	CROSSOVER 2	C1887	\$185
M001196620	19-662	08714729403623	6	55	0.070	LIMA	C1887	\$185
M001196660	19-666	08714729403654	6	55	0.070	RDC	C1887	\$185
M001196700	19-670	08714729403685	6	55	0.070	BATES 1*	C1887	\$185
M001196800	19-680	08714729403715	6	55	0.070	BATES 2*	C1887	\$185
M001196960	19-696	08714729403746	6	55	0.070	RDC 1	C1887	\$185

*Indicates two sideholes for this device

Mach1™ Peripheral Guide Catheter 7 F

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Length (cm)	I.D. (in)	Tip Shape	C-Code	Price
M001197000	19-700	08714729403753	7	55	0.081	HS	C1887	\$185
M001197040	19-704	08714729403760	7	55	0.081	RE S	C1887	\$185
M001197080	19-708	08714729403777	7	55	0.081	RE SS	C1887	\$185
M001197120	19-712	08714729403784	7	55	0.081	RE L	C1887	\$185
M001197160	19-716	08714729403791	7	55	0.081	C1	C1887	\$185
M001197200	19-720	08714729403807	7	55	0.081	C2	C1887	\$185
M001197230	19-723	08714729403814	7	55	0.081	MP	C1887	\$185
M001197260	19-726	08714729403821	7	90	0.081	MP	C1887	\$185
M001197300	19-730	08714729403838	7	55	0.081	MP2	C1887	\$185
M001197330	19-733	08714729403845	7	90	0.081	MP2	C1887	\$185
M001197380	19-738	08714729481584	7	90	0.081	40°	C1887	\$185
M001197450	19-745	08714729403852	7	55	0.081	ST	C1887	\$185
M001197480	19-748	08714729403869	7	90	0.081	ST	C1887	\$185
M001197530	19-753	08714729403876	7	55	0.081	CROSSOVER 1	C1887	\$185
M001197560	19-756	08714729403883	7	55	0.081	CROSSOVER 2	C1887	\$185
M001197620	19-762	08714729403890	7	55	0.081	LIMA	C1887	\$185
M001197660	19-766	08714729403906	7	55	0.081	RDC	C1887	\$185
M001197700	19-770	08714729403913	7	55	0.081	BATES 1*	C1887	\$185
M001197800	19-780	08714729403920	7	55	0.081	BATES 2*	C1887	\$185
M001197960	19-796	08714729403951	7	55	0.081	RDC 1	C1887	\$185

*Indicates two sideholes for this device

GUIDE CATHETERS

Guide Catheters

Mach1™ Peripheral Guide Catheter (cont.)

Mach1™ Peripheral Guide Catheter 8 F

UPN / Order Code	Catalog No.	GTIN	French Size (F)	Length (cm)	I.D. (in)	Tip Shape	C-Code	Price
M001198000	19-800	08714729403968	8	55	0.091	HS	C1887	\$185
M001198040	19-804	08714729403975	8	55	0.091	RE S	C1887	\$185
M001198080	19-808	08714729403982	8	55	0.091	RE SS	C1887	\$185
M001198120	19-812	08714729403999	8	55	0.091	RE L	C1887	\$185
M001198160	19-816	08714729404002	8	55	0.091	C1	C1887	\$185
M001198200	19-820	08714729404019	8	55	0.091	C2	C1887	\$185
M001198230	19-823	08714729404026	8	55	0.091	MP	C1887	\$185
M001198260	19-826	08714729404033	8	90	0.091	MP	C1887	\$185
M001198300	19-830	08714729404040	8	55	0.091	MP2	C1887	\$185
M001198330	19-833	08714729404057	8	90	0.091	MP2	C1887	\$185
M001198380	19-838	08714729481591	8	90	0.091	40°	C1887	\$185
M001198450	19-845	08714729404064	8	55	0.091	ST	C1887	\$185
M001198480	19-848	08714729404071	8	90	0.091	ST	C1887	\$185
M001198530	19-853	08714729404088	8	55	0.091	CROSSOVER 1	C1887	\$185
M001198560	19-856	08714729404095	8	55	0.091	CROSSOVER 2	C1887	\$185
M001198620	19-862	08714729404101	8	55	0.091	LIMA	C1887	\$185
M001198660	19-866	08714729404118	8	55	0.091	RDC	C1887	\$185
M001198700	19-870	08714729404125	8	55	0.091	BATES 1*	C1887	\$185
M001198960	19-896	08714729404163	8	55	0.091	RDC 1	C1887	\$185

*Indicates two sideholes for this device

Support Catheters

Rubicon™ Support Catheter

Rubicon™ Support Catheter

UPN / Order Code	Catalog No.	GTIN	Sheath Compatibility (F)	Length (cm)	Guidewire Compatibility (in)	Tip Shape	C-Code	Price
H74939240035061	39240-03506	08714729832065	5	65	0.035	Straight	C1887	\$1,190
H74939240035091	39240-03509	08714729832041	5	90	0.035	Straight	C1887	\$1,190
H74939240035131	39240-03513	08714729832072	5	135	0.035	Straight	C1887	\$1,190
H74939240035151	39240-03515	08714729832058	5	150	0.035	Straight	C1887	\$1,190
H74939239018091	39239-01809	08714729829232	4	90	0.018	Straight	C1887	\$1,190
H74939239018131	39239-01813	08714729829249	4	135	0.018	Straight	C1887	\$1,190
H74939239018151	39239-01815	08714729829256	4	150	0.018	Straight	C1887	\$1,190
H74939212014131	39212-01413	08714729799078	4	135	0.014	Straight	C1887	\$1,190
H74939212014151	39212-01415	08714729799085	4	150	0.014	Straight	C1887	\$1,190

Packaged 5 per box

Rubicon™ Control Support Catheter

Rubicon™ Control Support Catheter

UPN / Order Code	Catalog No.	GTIN	Sheath Compatibility	Length (cm)	Catheter Size	Tip Shape	Units Per Box	C-Code	List Price
H749394323506A1		191506019716	5F	65	35	Steep Angle	5	C-1887	\$2,346
H749394323506A2		191506019754	5F	90	35	Steep Angle	5	C-1887	\$2,346
H749394323513A1		191506019792	5F	135	35	Steep Angle	5	C-1887	\$2,346
H749394323515A1		191506019853	5F	150	35	Steep Angle	5	C-1887	\$2,346
H749394321809A1		191506019631	4F	90	18	Steep Angle	5	C-1887	\$2,346
H749394321813A1		191506019655	4F	135	18	Steep Angle	5	C-1887	\$2,346
H749394321815A1		191506019679	4F	150	18	Steep Angle	5	C-1887	\$2,346
H749394321409A1		191506019570	4F	90	14	Steep Angle	5	C-1887	\$2,346
H749394321413A1		191506019594	4F	135	14	Steep Angle	5	C-1887	\$2,346
H749394321415A1		191506019617	4F	150	14	Steep Angle	5	C-1887	\$2,346

GUIDE CATHETERS

Support Catheters

Rubicon™ Control Support Catheter (cont.)

Rubicon™ Control Support Catheter

UPN / Order Code	Catalog No.	GTIN	Sheath Compatibility	Length (cm)	Catheter Size	Tip Shape	Units Per Box	C-Code	List Price
H749394323506B1		191506019730	5F	65	35	Shallow Angle	5	C-1887	\$2,346
H749394323509B1		191506019778	5F	90	35	Shallow Angle	5	C-1887	\$2,346
H749394323513B1		191506019815	5F	135	35	Shallow Angle	5	C-1887	\$2,346
H749394323515B1		191506019877	5F	150	35	Shallow Angle	5	C-1887	\$2,346
H749394323513S1		191506019839	5F	135	35	Straight Tip	5	C-1887	\$2,346
H749394323515S1		191506019891	5F	150	35	Straight Tip	5	C-1887	\$2,346
H749394321815S1		191506019693	4F	150	18	Straight Tip	5	C-1887	\$2,346
H749394323506A0		191506019709	5F	65	35	Steep Angle	Single	C-1887	\$507
H749394323509A0		191506019747	5F	90	35	Steep Angle	Single	C-1887	\$507
H749394323513A0		191506019785	5F	135	35	Steep Angle	Single	C-1887	\$507
H749394323515A0		191506019846	5F	150	35	Steep Angle	Single	C-1887	\$507
H749394321809A0		191506019624	4F	90	18	Steep Angle	Single	C-1887	\$507
H749394321813A0		191506019648	4F	135	18	Steep Angle	Single	C-1887	\$507
H749394321815A0		191506019662	4F	150	18	Steep Angle	Single	C-1887	\$507
H749394321409A0		191506019563	4F	90	14	Steep Angle	Single	C-1887	\$507
H749394321413A0		191506019587	4F	135	14	Steep Angle	Single	C-1887	\$507
H749394321415A0		191506019600	4F	150	14	Steep Angle	Single	C-1887	\$507
H749394323506B0		191506019723	5F	65	35	Shallow Angle	Single	C-1887	\$507
H749394323509B0		191506019761	5F	90	35	Shallow Angle	Single	C-1887	\$507
H749394323513B0		191506019808	5F	135	35	Shallow Angle	Single	C-1887	\$507
H749394323515B0		191506019860	5F	150	35	Shallow Angle	Single	C-1887	\$507
H749394323513S0		191506019822	5F	135	35	Straight Tip	Single	C-1887	\$507
H749394323515S0		191506019884	5F	150	35	Straight Tip	Single	C-1887	\$507
H749394321815S0		191506019686	4F	150	18	Straight Tip	Single	C-1887	\$507

Amplatz Super Stiff™ Guidewires
3 mm J

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper Length (cm)	Tip Shape	Quantity	C-Code	Price
M001465001	46-500	08714729115847	0.035	145	7	3 mm J	5	C1769	\$415
M001465011	46-501	08714729115854	0.035	180	7	3 mm J	5	C1769	\$415
M001465021	46-502	08714729115861	0.035	260	7	3 mm J	5	C1769	\$415
M001465031	46-503	08714729115915	0.038	145	6	3 mm J	5	C1769	\$415
M001465041	46-504	08714729115922	0.038	180	6	3 mm J	5	C1769	\$415

Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper Length (cm)	Tip Shape	Quantity	C-Code	Price
M001465171	46-517	08714729115885	0.038	145	6	Straight	5	C1769	\$415
M001465191	46-519	08714729115892	0.038	180	6	Straight	5	C1769	\$415
M001465201	46-520	08714729115908	0.038	260	6	Straight	5	C1769	\$415
M001465231	46-523	08714729115809	0.035	145	7	Straight	5	C1769	\$415
M001465251	46-525	08714729115823	0.035	180	7	Straight	5	C1769	\$415
M001465261	46-526	08714729115830	0.035	260	7	Straight	5	C1769	\$415
M001465631	46-563	08714729114604	0.035	75	7	Straight	5	C1769	\$415
M001465641	46-564	08714729115878	0.038	75	6	Straight	5	C1769	\$415

Short Taper

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Length (cm)	Tip Shape	Quantity	C-Code	Price
M001465091	46-509	08714729115953	0.035	260	1	Straight	5	C1769	\$415
M001465101	46-510	08714729115960	0.038	260	1	Straight	5	C1769	\$415
M001465181	46-518	08714729115939	0.038	145	4	Straight	5	C1769	\$415
M001465241	46-524	08714729115946	0.035	145	4	Straight	5	C1769	\$415

Meier Guidewire

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Length (cm)	Tip Shape	Quantity	C-Code	Price
H965SCH306001	SCH-30600	08714729227885	0.035	185	15	J	5	C1769	\$532
H965SCH306011	SCH-30601	08714729227908	0.035	300	10	C	5	C1769	\$532
M001306021	30-602	08714729455790	0.035	260	10	C	5	C1769	\$532

Platinum Plus™ Guidewires
Platinum Plus™ Guidewires with Glidex™ Hydrophilic Coating

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001467120	46-712	08714729486886	0.014	300	Short	Straight	1	C1769	\$195
M001467300	46-730	08714729326199	0.018	60	Short	Straight	1	C1769	\$195
M001467310	46-731	08714729326205	0.018	180	Short	Straight	1	C1769	\$195
M001467320	46-732	08714729326212	0.018	260	Short	Straight	1	C1769	\$195
M001467400	46-740	08714729005049	0.014	180	Short	Straight	1	C1769	\$195

Platinum Plus™ Guidewires with Mediglide™ Silicone Coating

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Quantity	C-Code	Price
M001466000	46-600	08714729035947	0.018	145	Long	1	C1769	\$195
M001466010	46-601	08714729035954	0.018	180	Long	1	C1769	\$195
M001466020	46-602	08714729035961	0.018	260	Long	1	C1769	\$195
M001466040	46-604	08714729114567	0.018	180	Short	1	C1769	\$195
M001466050	46-605	08714729114574	0.018	260	Short	1	C1769	\$195
M001466100	46-610	08714729114581	0.025	180	Short	1	C1769	\$195
M001466110	46-611	08714729114598	0.025	260	Short	1	C1769	\$195

SAVION™ 0.014 Guidewires
SAVION FLX™ Guidewires

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Load (grams)	Tip Shape	Quantity	C-Code	Price
H7493942318510	N/A	08714729978831	0.014	185	2.9	Straight	1	C1769	\$490
H749394231851J0	N/A	08714729978848	0.014	185	2.9	Angled	1	C1769	\$490
H7493942230010	N/A	08714729978817	0.014	300	2.9	Straight	1	C1769	\$490
H749394223001J0	N/A	08714729978824	0.014	300	2.9	Angled	1	C1769	\$490

SAVION DLVR™ Guidewires

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Load (grams)	Tip Shape	Quantity	C-Code	Price
H7493942118210	N/A	08714729978879	0.014	182	0.8	Straight	1	C1769	\$425
H749394211821J0	N/A	08714729978886	0.014	182	0.8	Angled	1	C1769	\$425
H7493942030010	N/A	08714729978855	0.014	300	0.8	Straight	1	C1769	\$425
H749394203001J0	N/A	08714729978862	0.014	300	0.8	Angled	1	C1769	\$425

Thruway™ 0.014 Guidewires
130 cm Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492891	49-289	08714729717102	0.014	130	Long	Straight	1	C1769	\$195

190 cm J-Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492931	49-293	08714729717140	0.014	190	Short	J-Tip	1	C1769	\$200

190 cm Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492921	49-292	08714729717133	0.014	190	Short	Straight	1	C1769	\$200
M001492941	49-294	08714729717157	0.014	190	Long	Straight	1	C1769	\$200
M001492961	49-296	08714729717171	0.014	190	Extra Long	Straight	1	C1769	\$200

300 cm Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492971	49-297	08714729717188	0.014	300	Short	Straight	1	C1769	\$210
M001492981	49-298	08714729717195	0.014	300	Long	Straight	1	C1769	\$210
M001492991	49-299	08714729717201	0.014	300	Extra Long	Straight	1	C1769	\$210

Thruway™ 0.018 Guidewires
130 cm Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492751	49-275	08714729405542	0.018	130	Short	Straight	1	C1769	\$160
M001492771	49-277	08714729405566	0.018	130	Long	Straight	1	C1769	\$160

190 cm J-Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492801	49-280	08714729405597	0.018	190	Short	J-Tip	1	C1769	\$170
M001492821	49-282	08714729405610	0.018	190	Long	J-Tip	1	C1769	\$170

190 cm Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492791	49-279	08714729405580	0.018	190	Short	Straight	1	C1769	\$170
M001492811	49-281	08714729405603	0.018	190	Long	Straight	1	C1769	\$170

Thruway™ 0.018 Guidewires
300 cm J-Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492841	49-284	08714729405634	0.018	300	Short	J-Tip	1	C1769	\$185
M001492861	49-286	08714729405658	0.018	300	Long	J-Tip	1	C1769	\$185

300 cm Straight

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Taper	Tip Shape	Quantity	C-Code	Price
M001492831	49-283	08714729405627	0.018	300	Short	Straight	1	C1769	\$185
M001492851	49-285	08714729405641	0.018	300	Long	Straight	1	C1769	\$185

V-14™ ControlWire™ Guidewire

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper Length (cm)	Tip Shape	Quantity	C-Code	Price
H74939216718210	39216-71821	08714729799207	0.014	182	8	Straight	1	C1769	\$315
H74939216718220	39216-71822	08714729799214	0.014	182	8	Angled	1	C1769	\$315
H74939216730010	39216-73001	08714729799221	0.014	300	8	Straight	1	C1769	\$315
H74939216730020	39216-73002	08714729799238	0.014	300	8	Angled	1	C1769	\$315
H74939216918210	39216-91821	08714729799245	0.014	182	11	Straight	1	C1769	\$315
H74939216918220	39216-91822	08714729799252	0.014	182	11	Angled	1	C1769	\$315
H74939216930010	39216-93001	08714729799269	0.014	300	11	Straight	1	C1769	\$315
H74939216930020	39216-93002	08714729799276	0.014	300	11	Angled	1	C1769	\$315

V-18™ ControlWire™ Guidewire

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper Length (cm)	Tip Shape	Unit of Measure	C-Code	Price
M001468480	46-848	08714729715566	0.018	110	8	Straight	1	C1769	\$180
M001468500	46-850	08714729070184	0.018	150	8	Straight	1	C1769	\$180
M001468520	46-852	08714729070191	0.018	200	8	Straight	1	C1769	\$180
M001468540	46-854	08714729070207	0.018	300	8	Straight	1	C1769	\$180
M001468560	46-856	08714729070214	0.018	150	12	Straight	1	C1769	\$180
M001468580	46-858	08714729070221	0.018	200	12	Straight	1	C1769	\$180
M001468600	46-860	08714729070238	0.018	300	12	Straight	1	C1769	\$180

Victory™ 14 Guidewire

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Load (Grams)	Tip Shape	Quantity	C-Code	Price
H74939470195120	39470-19512	00816349012433	0.014	195	12	Straight	1	C1769	\$325
H74939470300120	39470-30012	00816349012372	0.014	300	12	Straight	1	C1769	\$325
H74939470195180	39470-19518	00816349012389	0.014	195	18	Straight	1	C1769	\$325
H74939470300180	39470-30018	00816349012440	0.014	300	18	Straight	1	C1769	\$325
H74939471195250	39471-19525	00816349012501	0.014	195	25	Straight	1	C1769	\$390
H74939471300250	39471-30025	00816349012464	0.014	300	25	Straight	1	C1769	\$390
H74939471195300	39471-19530	00816349012419	0.014	195	30	Straight	1	C1769	\$390
H74939471300300	39471-30030	00816349012402	0.014	300	30	Straight	1	C1769	\$390

Victory™ 18 Guidewire

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Tip Load (Grams)	Tip Shape	Quantity	C-Code	Price
H74939472195120	39472-19512	00816349012358	0.018	195	12	Straight	1	C1769	\$285
H74939472300120	39472-30012	00816349012457	0.018	300	12	Straight	1	C1769	\$285
H74939472195180	39472-19518	00816349012471	0.018	195	18	Straight	1	C1769	\$285
H74939472300180	39472-30018	00816349012365	0.018	300	18	Straight	1	C1769	\$285
H74939473195250	39473-19525	00816349012495	0.018	195	25	Straight	1	C1769	\$360
H74939473300250	39473-30025	00816349012396	0.018	300	25	Straight	1	C1769	\$360
H74939473195300	39473-19530	00816349012488	0.018	195	30	Straight	1	C1769	\$360
H74939473300300	39473-30030	00816349012426	0.018	300	30	Straight	1	C1769	\$360

ZIPwire™ Guidewire - J-Tip
J-Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146350B1	46-350B	08714729755050	0.035	150	3	J-Tip	5	C1769	\$515

ZIPwire™ Guidewire - Long Taper Standard Shaft
Angled Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146325B1	46-325B	08714729754930	0.035	150	8	Angled	5	C1769	\$515

ZIPwire™ Guidewire - Long Taper Stiff Shaft
Angled Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146318B1	46-318B	08714729754954	0.035	150	5	Angled	5	C1769	\$515
M00146319B1	46-319B	08714729754961	0.035	260	5	Angled	5	C1769	\$515

Straight Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146320B1	46-320B	08714729761686	0.035	150	5	Straight	5	C1769	\$515

ZIPwire™ Guidewire - Small Vessel
Angled Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146125B1	46-125B	08714729754992	0.025	150	3	Angled	5	C1769	\$595
M00146231B1	46-231B	08714729770268	0.018	150	3	Angled	5	C1769	\$555
M00146232B1	46-232B	08714729770282	0.018	180	3	Angled	5	C1769	\$555
M00146240B1	46-240B	08714729755012	0.025	180	3	Angled	5	C1769	\$595
M00146254B1	46-254B	08714729755036	0.025	260	3	Angled	5	C1769	\$595

Straight Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146126B1	46-126B	08714729755005	0.025	150	3	Straight	5	C1769	\$595
M00146235B1	46-235B	08714729770299	0.018	150	3	Straight	5	C1769	\$555
M00146236B1	46-236B	08714729770312	0.018	180	3	Straight	5	C1769	\$555
M00146249B1	46-249B	08714729755029	0.025	180	3	Straight	5	C1769	\$595
M00146255B1	46-255B	08714729755043	0.025	260	3	Straight	5	C1769	\$595

Starter™ Guidewires (PTFE Coated)

Straight Fixed Core

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491001	49-100	08714729326458	0.025	150	3	Straight	10	C1769	\$175
M001491011	49-101	08714729201281	0.035	150	4	Straight	10	C1769	\$175
M001491021	49-102	08714729326465	0.038	150	4	Straight	10	C1769	\$175
M001491031	49-103	08714729326472	0.035	260	4	Straight	5	C1769	\$140
M001492311	49-231	08714729327172	0.035	80	4	Straight	10	C1769	\$175

Starter™ Guidewires (PTFE Coated)
Straight Moveable Core

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491101	49-110	08714729201298	0.035	150	0	Straight	10	C1769	\$225
M001492351	49-235	08714729327189	0.035	150	12	Straight	10	C1769	\$225

J Curved Fixed Core

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491151	49-115	08714729326496	0.035	150	4	1.5 mm J	10	C1769	\$205
M001491161	49-116	08714729326502	0.025	150	3	3 mm J	10	C1769	\$205
M001491171	49-117	08714729326519	0.035	125	10	3 mm J FS†	10	C1769	\$205
M001491181	49-118	08714729201304	0.035	150	10	3 mm J FS†	10	C1769	\$205
M001491191	49-119	08714729326526	0.038	150	11	3 mm J FS†	10	C1769	\$205
M001491201	49-120	08714729326533	0.035	180	10	3 mm J FS†	5	C1769	\$205
M001491211	49-121	08714729326540	0.035	260	10	3 mm J FS†	5	C1769	\$205
M001491221	49-122	08714729326557	0.035	150	8	15 mm J FS†	10	C1769	\$205
M001491231	49-123	08714729326564	0.038	150	8	15 mm J FS†	10	C1769	\$205
M001492201	49-220	08714729327066	0.035	150	4	1.5 mm J FS†	10	C1769	\$205
M001492241	49-224	08714729327103	0.035	150	4	3 mm J FS†	10	C1769	\$205
M001492291	49-229	08714729327158	0.035	150	9	7.5 mm J FS†	10	C1769	\$205
M001492301	49-230	08714729327165	0.035	80	10	3 mm J FS†	10	C1769	\$205

†Finger Straightenable

J Curved Moveable Core

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491301	49-130	08714729326571	0.035	150	0	1.5 mm J	10	C1769	\$240
M001491321	49-132	08714729201311	0.035	150	0	3 mm J FS†	10	C1769	\$240
M001491331	49-133	08714729326595	0.038	150	0	3 mm J FS†	10	C1769	\$240
M001491341	49-134	08714729326601	0.035	260	0	3 mm J FS†	5	C1769	\$150
M001491381	49-138	08714729326632	0.035	150	12	3 mm J Tapered	10	C1769	\$240
M001491391	49-139	08714729326649	0.038	150	12	3 mm J Tapered	10	C1769	\$240

†Finger Straightenable

Bentson

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491451	49-145	08714729201328	0.035	150	15	Straight	10	C1769	\$270
M001491461	49-146	08714729326656	0.038	150	15	Straight	10	C1769	\$270
M001491471	49-147	08714729326663	0.035	180	15	Straight	5	C1769	\$180
M001491481	49-148	08714729326670	0.035	260	15	Straight	5	C1769	\$180

Newton

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491501	49-150	08714729326687	0.035	150	10	Straight	10	C1769	\$225
M001491511	49-151	08714729326694	0.035	150	15	Straight	10	C1769	\$225
M001491521	49-152	08714729326700	0.035	150	15	3 mm J FS†	10	C1769	\$225
M001491531	49-153	08714729326717	0.035	150	15	15 mm J FS†	10	C1769	\$225
M001492361	49-236	08714729327196	0.018	150	14	Straight	10	C1769	\$225

†Finger Straightenable

Rosen

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491551	49-155	08714729326724	0.035	150	4	1.5 mm J	10	C1769	\$190
M001491561	49-156	08714729326731	0.035	180	4	1.5 mm J	5	C1769	\$150
M001491571	49-157	08714729326748	0.035	260	4	1.5 mm J	5	C1769	\$150

Starter™ Guidewires (PTFE/Heparin Coated)

Straight Fixed Core

UPN / Order Code	Catalog No.		Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491601	49-160	08714729326755	0.035	150	4	Straight	10	C1769	\$175
M001491611	49-161	08714729326762	0.035	260	4	Straight	5	C1769	\$140

J Curved Fixed Core

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491651	49-165	08714729326779	0.035	150	10	3 mm JFS†	10	C1769	\$205
M001491661	49-166	08714729326786	0.038	150	11	3 mm JFS†	10	C1769	\$205
M001491671	49-167	08714729326793	0.035	180	10	3 mm JFS†	5	C1769	\$160
M001491681	49-168	08714729326809	0.035	260	10	3 mm JFS†	5	C1769	\$160
M001491691	49-169	08714729326816	0.035	150	8	15 mm JFS†	10	C1769	\$205

J Curve Moveable Core

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491701	49-170	08714729326823	0.035	150	0	3 mm JFS†	10	C1769	\$240

†Finger Straightenable

Bentson

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491751	49-175	08714729326830	0.035	150	15	Straight	10	C1769	\$270
M001491761	49-176	08714729326847	0.038	150	15	Straight	10	C1769	\$270
M001491771	49-177	08714729326854	0.035	180	15	Straight	5	C1769	\$180

Newton

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491801	49-180	08714729326861	0.035	150	10	Straight	10	C1769	\$225
M001491811	49-181	08714729326878	0.035	150	15	Straight	10	C1769	\$225
M001491821	49-182	08714729326885	0.035	150	10	3 mm J FS†	10	C1769	\$225
M001491831	49-183	08714729326892	0.038	150	10	3 mm J FS†	10	C1769	\$225

†Finger Straightenable

Rosen

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M001491851	49-185	08714729326908	0.035	180	4	1.5 mm J	5	C1769	\$150

†Finger Straightenable

ZIPwire™ Hydrophilic Guidewires

ZIPwire™ Guidewire - Standard Shaft

Angled Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146151B1	46-151B	08714729754695	0.035	150	3	Angled	5	C1769	\$465
M00146152B1	46-152B	08714729754718	0.035	180	3	Angled	5	C1769	\$465
M00146154B1	46-154B	08714729754732	0.035	260	3	Angled	5	C1769	\$465
M00146161B1	46-161B	08714729754756	0.038	150	3	Angled	5	C1769	\$465
M00146162B1	46-162B	08714729754770	0.038	180	3	Angled	5	C1769	\$465
M00146164B1	46-164B	08714729754794	0.038	260	3	Angled	5	C1769	\$465

ZIPwire™ Guidewire - Standard Shaft

Straight Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146155B1	46-155B	08714729754701	0.035	150	3	Straight	5	C1769	\$465
M00146156B1	46-156B	08714729754725	0.035	180	3	Straight	5	C1769	\$465
M00146159B1	46-159B	08714729754749	0.035	260	3	Straight	5	C1769	\$465
M00146165B1	46-165B	08714729754763	0.038	150	3	Straight	5	C1769	\$465
M00146166B1	46-166B	08714729754787	0.038	180	3	Straight	5	C1769	\$465

ZIPwire™ Guidewire - Stiff Shaft

Angled Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146300B1	46-300B	08714729754817	0.035	80	3	Angled	5	C1769	\$515
M00146301B1	46-301B	08714729754886	0.038	80	3	Angled	5	C1769	\$515
M00146305B1	46-305B	08714729754824	0.035	150	3	Angled	5	C1769	\$515
M00146308B1	46-308B	08714729754848	0.035	180	3	Angled	5	C1769	\$515
M00146311B1	46-311B	08714729754893	0.038	150	3	Angled	5	C1769	\$515
M00146314B1	46-314B	08714729754916	0.038	180	3	Angled	5	C1769	\$515
M00146316B1	46-316B	08714729754862	0.035	260	3	Angled	5	C1769	\$515

Straight Tip

UPN / Order Code	Catalog No.	GTIN	Diameter (in)	Length (cm)	Taper (cm)	Tip Shape	Quantity	C-Code	Price
M00146306B1	46-306B	08714729754831	0.035	150	3	Straight	5	C1769	\$515
M00146309B1	46-309B	08714729754855	0.035	180	3	Straight	5	C1769	\$515
M00146312B1	46-312B	08714729754909	0.038	150	3	Straight	5	C1769	\$515
M00146315B1	46-315B	08714729754923	0.038	180	3	Straight	5	C1769	\$515
M00146317B1	46-317B	08714729754879	0.035	260	3	Straight	5	C1769	\$515

Fathom™ Peripheral Embolization Guidewires

Fathom™ 14 Steerable Guidewires

UPN / Order Code	Catalog No.	GTIN	Total Length (cm)	Tip Length (cm)	Proximal/ Distal OD (in)	Tip Shape	C-Code	Price
M001508100	50-810	08714729801023	200	10	0.014	Straight	C1769	\$485
M001508110	50-811	08714729801030	200	10	0.014	Angled	C1769	\$505
M001508140	50-814	08714729801047	300	10	0.014	Straight	C1769	\$485
M001508150	50-815	08714729801054	300	10	0.014	Angled	C1769	\$505

Fathom™ 16 Steerable Guidewires

UPN / Order Code	Catalog No.	GTIN	Total Length (cm)	Tip Length (cm)	Proximal/ Distal OD (in)	Tip Shape	C-Code	Price
M001509000	50-900	08714729764410	140	10	0.016	Straight	C1769	\$485
M001509010	50-901	08714729764427	140	20	0.016	Straight	C1769	\$485
M001509100	50-910	08714729762553	180	10	0.016	Straight	C1769	\$485
M001509110	50-911	08714729762560	180	20	0.016	Straight	C1769	\$485
M001509120	N/A	08714729939658	180	10	0.016	Angled	C1769	\$505
M001509200	N/A	08714729961116	200	10	0.016	Straight	C1769	\$485
M001509210	N/A	08714729961123	200	10	0.016	Angled	C1769	\$505
M001509300	N/A	08714729961130	215	10	0.016	Straight	C1769	\$485
M001509310	N/A	08714729961147	215	10	0.016	Angled	C1769	\$505

Transend™ Peripheral Embolization Guidewires

Transend™ - 14 Steerable Guidewires

UPN / Order Code	Catalog No.	GTIN	Total Length (cm)	Shapeable Tip Length (cm)	Proximal/ Distal OD (in)	Tip Shape	C-Code	Price
M001468100	46-810	08714729326274	135	2	0.014	Straight	C1769	\$350
M001468110	46-811	08714729326281	165	2	0.014	Straight	C1769	\$350
M001468010	46-801	08714729326243	190	2	0.014	Straight	C1769	\$350

Transend™ - 18 Steerable Guidewires

UPN / Order Code	Catalog No.	GTIN	Total Length (cm)	Tip Length (cm)	Proximal/ Distal OD (in)	Tip Shape	C-Code	Price
M001468120	46-812	08714729326298	135	2	0.018	Straight	C1769	\$350
M001468130	46-813	08714729326304	165	2	0.018	Straight	C1769	\$350

Microcatheters

TruSelect™ Microcatheter

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Tip Shape	C-Code	Price
M001394101050	41010-50	0814729976608	105	2.8/2.0	0.021	Straight	C1887	\$935
M001394101300	41013-00	08714729976615	130	2.8/2.0	0.021	Straight	C1887	\$935
M001394101550	41015-50	08714729976622	155	2.8/2.0	0.021	Straight	C1887	\$935
M001394101750	41017-50	08714729976639	175	2.8/2.0	0.021	Straight	C1887	\$990
M001394111050	41110-50	08714729976646	105	2.8/2.0	0.021	Bern	C1887	\$935
M001394111300	41113-00	08714729976653	130	2.8/2.0	0.021	Bern	C1887	\$935
M001394111550	41115-50	08714729976660	155	2.8/2.0	0.021	Bern	C1887	\$935
M001394111750	41117-50	08714729976677	175	2.8/2.0	0.021	Bern	C1887	\$990

Microcatheters

Direxion™ Torqueable Microcatheters

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Tip Shape	RO Markers	C-Code	Price
M001195200	19-520	08714729839316	105	2.7/2.4	0.021	Straight	1	C1887	\$700
M001195210	19-521	08714729839323	130	2.7/2.4	0.021	Straight	1	C1887	\$700
M001195220	19-522	08714729839330	155	2.7/2.4	0.021	Straight	1	C1887	\$700
M001195230	19-523	08714729839347	105	2.7/2.4	0.021	Bern	1	C1887	\$760
M001195240	19-524	08714729839354	130	2.7/2.4	0.021	Bern	1	C1887	\$760
M001195250	19-525	08714729839361	155	2.7/2.4	0.021	Bern	1	C1887	\$760
M001195260	19-526	08714729839378	105	2.7/2.4	0.021	J	1	C1887	\$760
M001195270	19-527	08714729839385	130	2.7/2.4	0.021	J	1	C1887	\$760
M001195280	19-528	08714729839392	155	2.7/2.4	0.021	J	1	C1887	\$760
M001195290	19-529	08714729839408	105	2.7/2.4	0.021	Swan	1	C1887	\$760
M001195300	19-530	08714729839415	130	2.7/2.4	0.021	Swan	1	C1887	\$760
M001195310	19-531	08714729839422	155	2.7/2.4	0.021	Swan	1	C1887	\$760
M001195320	19-532	08714729839439	130	2.7/2.4	0.021	Straight	2	C1887	\$760
M001195330	19-533	08714729839446	155	2.7/2.4	0.021	Straight	2	C1887	\$760
M001195340	19-534	08714729839453	130	2.7/2.4	0.021	Bern	2	C1887	\$815
M001195350	19-535	08714729839460	155	2.7/2.4	0.021	Bern	2	C1887	\$815

Direxion™ Fathom™ - 16 System Pre-Loaded Torqueable Microcatheters

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Tip Shape	Fathom™ 16 Steerable Guidewire Overall Length (cm)	C-Code	Price
M001195610	19-561	08714729839606	130	2.7/2.4	0.021	Straight	180	C1887	\$1,065
M001195620	19-562	08714729839613	155	2.7/2.4	0.021	Straight	180	C1887	\$1,065
M001195640	19-564	08714729839637	130	2.7/2.4	0.021	Bern	180	C1887	\$1,115
M001195650	19-565	08714729839644	155	2.7/2.4	0.021	Bern	180	C1887	\$1,115
M001195660	19-566	08714729839651	130	2.7/2.4	0.021	Straight (2RO)	180	C1887	\$1,115
M001195670	19-567	08714729839668	130	2.7/2.4	0.021	Bern (2RO)	180	C1887	\$1,115

Direxion™ Transend™ - 14 System Pre-Loaded Torqueable Microcatheters

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Tip Shape	Transend™ 14 Steerable Guidewire Overall Length (cm)	C-Code	Price
M001195810	19-581	08714729839743	130	2.7/2.4	0.021	Straight	165	C1887	\$925
M001195840	19-584	08714729839774	130	2.7/2.4	0.021	Bern	165	C1887	\$980
M001195850	19-585	08714729839781	155	2.7/2.4	0.021	Bern	165	C1887	\$980

Direxion™ HI-FLO™ Torqueable Microcatheters

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Tip Shape	C-Code	Price
M001195400	19-540	08714729839477	105	3.0/2.8	0.027	Straight	C1887	\$700
M001195410	19-541	08714729839484	130	3.0/2.8	0.027	Straight	C1887	\$700
M001195420	19-542	08714729839491	155	3.0/2.8	0.027	Straight	C1887	\$700
M001195430	19-543	08714729839507	105	3.0/2.8	0.027	Bern	C1887	\$760
M001195440	19-544	08714729839514	130	3.0/2.8	0.027	Bern	C1887	\$760
M001195450	19-545	08714729839521	155	3.0/2.8	0.027	Bern	C1887	\$760
M001195470	19-547	08714729839545	130	3.0/2.8	0.027	J	C1887	\$760
M001195480	19-548	08714729839552	155	3.0/2.8	0.027	J	C1887	\$760
M001195500	19-550	08714729839576	130	3.0/2.8	0.027	Swan	C1887	\$760
M001195510	19-551	08714729839583	155	3.0/2.8	0.027	Swan	C1887	\$760

Direxion™ HI-FLO Fathom™ - 16 System Pre-Loaded Torqueable Microcatheters

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Tip Shape	Fathom™ 16 Steerable Guidewire Overall Length (cm)	C-Code	Price
M001195710	19-571	08714729839682	130	3.0/2.8	0.027	Straight	180	C1887	\$1,065
M001195720	19-572	08714729839699	155	3.0/2.8	0.027	Straight	180	C1887	\$1,065
M001195740	19-574	08714729839712	130	3.0/2.8	0.027	Bern	180	C1887	\$1,115
M001195750	19-575	08714729839729	155	3.0/2.8	0.027	Bern	180	C1887	\$1,115

Microcatheters

Direxion™ HI-FLO Transend™ - 18 System Pre-Loaded Torqueable Microcatheters

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Tip Shape	Transend™ 18 Steerable Guidewire Overall Length (cm)	C-Code	Price
M001195910	19-591	08714729839828	130	3.0/2.8	0.027	Straight	165	C1887	\$925
M001195940	19-594	08714729839842	130	3.0/2.8	0.027	Bern	165	C1887	\$980

Renegade™ Fiber Braided Microcatheters with Hydro Pass™ Hydrophilic Coating

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Distal Tip Length (cm)	C-Code	Price
M001182520	18-252	08714729235125	130	3.0 / 2.5	0.021	20	C1887	\$590
M001182510	18-251	08714729235118	150	3.0 / 2.5	0.021	10	C1887	\$590
M001182530	18-253	08714729322832	150	3.0 / 2.5	0.021	20	C1887	\$590

Renegade™ STC-18 Microcatheters

Angled

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Distal Tip Length (cm)	Tip Design	C-Code	Price
M001181270	18-127	08714729444404	105	3.0 / 2.4	0.021	20	Angled	C1887	\$655
M001181280	18-128	08714729444411	105	3.0 / 2.4	0.021	30	Angled	C1887	\$655
M001181330	18-133	08714729444466	130	3.0 / 2.4	0.021	20	Angled	C1887	\$655
M001181340	18-134	08714729444473	130	3.0 / 2.4	0.021	30	Angled	C1887	\$655
M001181390	18-139	08714729444527	150	3.0 / 2.4	0.021	20	Angled	C1887	\$655
M001181400	18-140	08714729444534	150	3.0 / 2.4	0.021	30	Angled	C1887	\$655

Straight

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Distal Tip Length (cm)	Tip Design	C-Code	Price
M001181250	18-125	08714729444381	105	3.0 / 2.4	0.021	20	Straight	C1887	\$635
M001181260	18-126	08714729444398	105	3.0 / 2.4	0.021	30	Straight	C1887	\$635
M001181310	18-131	08714729444442	130	3.0 / 2.4	0.021	20	Straight	C1887	\$635
M001181320	18-132	08714729444459	130	3.0 / 2.4	0.021	30	Straight	C1887	\$635
M001181370	18-137	08714729444503	150	3.0 / 2.4	0.021	20	Straight	C1887	\$635
M001181380	18-138	08714729444510	150	3.0 / 2.4	0.021	30	Straight	C1887	\$635

Renegade™ HI-FLO™ Microcatheters

UPN / Order Code	Catalog No.	GTIN	Usable Length (cm)	Proximal/ Distal OD (F)	Lumen Inner Diameter (in)	Distal Tip Length (cm)	Tip Design	C-Code	Price
M001182840	18-284	08714729280590	105	3.0 / 2.8	0.027	10	Straight	C1887	\$590
M001182850	18-285	08714729280606	105	3.0 / 2.8	0.027	20	Straight	C1887	\$590
M001182860	18-286	08714729280613	115	3.0 / 2.8	0.027	10	Straight	C1887	\$590
M001182870	18-287	08714729280620	115	3.0 / 2.8	0.027	20	Straight	C1887	\$590
M001182880	18-288	08714729280637	135	3.0 / 2.8	0.027	10	Straight	C1887	\$590
M001182890	18-289	08714729280644	135	3.0 / 2.8	0.027	20	Straight	C1887	\$590
M001182900	18-290	08714729280651	150	3.0 / 2.8	0.027	10	Straight	C1887	\$590
M001182910	18-291	08714729280668	150	3.0 / 2.8	0.027	20	Straight	C1887	\$590

Microcatheter Kits and Systems

Renegade™ HI-FLO™ Microcatheter Kits with Transend™-18 Guidewires

UPN / Order Code	Catalog No.	GTIN	Transend-18 Guidewire Overall Length (cm)	Transend-18 Guidewire Distal Tip Length (cm)	Renegade HI-FLO Microcatheter Usable Length (cm)	Renegade HI-FLO Microcatheter Distal Tip Length (cm)	C-Code	Price
M001182980	18-298	08714729280736	135	38	105	10	C1887	\$840
M001182990	18-299	08714729280743	135	38	105	20	C1887	\$840
M001183000	18-300	08714729280750	135	38	115	10	C1887	\$840
M001183010	18-301	08714729280767	135	38	115	20	C1887	\$840
M001183020	18-302	08714729280774	165	38	135	10	C1887	\$840
M001183030	18-303	08714729280781	165	38	135	20	C1887	\$840

Renegade™ HI-FLO™ Fathom™ Pre-Loaded Systems with Fathom™-16 Guidewires

UPN / Order Code	Catalog No.	GTIN	Fathom-16 Steerable Guidewire Overall Length (cm)	Fathom-16 Steerable Guidewire Distal Tip Length (cm)	Renegade HI-FLO Microcatheter Usable Length (cm)	Renegade HI-FLO Microcatheter Distal Tip Length (cm)	C-Code	Price
M001184500	18-450	08714729790815	140	10	105	10	C1887	\$955
M001184510	18-451	08714729790822	140	10	105	20	C1887	\$955
M001184520	18-452	08714729790839	140	10	115	10	C1887	\$955
M001184530	18-453	08714729790846	140	10	115	20	C1887	\$955
M001184540	18-454	08714729790853	180	10	135	10	C1887	\$955
M001184550	18-455	08714729790860	180	10	135	20	C1887	\$955
M001184560	18-456	08714729790877	180	10	150	10	C1887	\$955
M001184570	18-457	08714729790884	180	10	150	20	C1887	\$955

Occlusion Balloons

Standard Occlusion Balloon Catheters

UPN / Order Code	Catalog No.	GTIN	Size (F)	Lumens	Usable Length (cm)	Recommended Guidewire	Inflated Balloon Diameter (mm)	C-Code	Price
M001171020	17-102	08714729112389	7	2	65	0.038	11.5	C2628	\$365
M001171030	17-103	08714729112396	7	2	100	0.038	11.5	C2628	\$365

Berenstein Occlusion Balloon Catheters

UPN / Order Code	Catalog No.	GTIN	Size (F)	Lumens	Usable Length (cm)	Recommended Guidewire	Inflated Balloon Diameter (mm)	C-Code	Price
M001173010	17-301	08714729112495	6	2	80	0.038	8.5 / 11.5	C2628	\$395

Conformable Embolic

Obsidio™

UPN / Order Code	Catalog No.	GTIN	Description	Volume (mL)	Box (Units)	C-Code	Price
M0013972101010	10-101	00191506039332	Obsidio Syringe	1	1	n/a	\$6,370

EMBOLD™ Fibered Detachable Coil System

UPN / Order Code	GTIN	Description	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001394240020040	08714729983521	EMBOLD FIBERED 2X4	2	4	1	N/A	\$1,460
M001394240020060	08714729983538	EMBOLD FIBERED 2X6	2	6	1	N/A	\$1,460
M001394240030060	08714729983545	EMBOLD FIBERED 3X6	3	6	1	N/A	\$1,460
M001394240030120	08714729983552	EMBOLD FIBERED 3X12	3	12	1	N/A	\$1,620
M001394240040080	08714729983569	EMBOLD FIBERED 4X8	4	8	1	N/A	\$1,460
M001394240040150	08714729983576	EMBOLD FIBERED 4X15	4	15	1	N/A	\$1,620
M001394240050080	08714729983583	EMBOLD FIBERED 5X8	5	8	1	N/A	\$1,460
M001394240050150	08714729983590	EMBOLD FIBERED 5X15	5	15	1	N/A	\$1,620
M001394240060100	08714729983606	EMBOLD FIBERED 6X10	6	10	1	N/A	\$1,620
M001394240060200	08714729983613	EMBOLD FIBERED 6X20	6	20	1	N/A	\$1,620
M001394240080200	08714729983620	EMBOLD FIBERED 8X20	8	20	1	N/A	\$1,620
M001394240100200	08714729983637	EMBOLD FIBERED 10X20	10	20	1	N/A	\$1,620
M001394240100300	08714729983644	EMBOLD FIBERED 10X30	10	30	1	N/A	\$2,050
M001394240100500	08714729983651	EMBOLD FIBERED 10X50	10	50	1	N/A	\$2,050
M001394240120200	08714729983668	EMBOLD FIBERED 12X20	12	20	1	N/A	\$1,620
M001394240120300	08714729983675	EMBOLD FIBERED 12X30	12	30	1	N/A	\$2,050
M001394240120500	08714729983682	EMBOLD FIBERED 12X50	12	50	1	N/A	\$2,050
M001394240140200	08714729983699	EMBOLD FIBERED 14X20	14	20	1	N/A	\$1,620
M001394240140300	08714729983705	EMBOLD FIBERED 14X30	14	30	1	N/A	\$2,050
M001394240140500	08714729983712	EMBOLD FIBERED 14X50	14	50	1	N/A	\$2,050
M001394240180500	08714729983729	EMBOLD FIBERED 18X50	18	50	1	N/A	\$2,050
M001394240200500	08714729983736	EMBOLD FIBERED 20X50	20	50	1	N/A	\$2,050
M001394240220600	08714729983743	EMBOLD FIBERED 22X60	22	60	1	N/A	\$2,050
M001394240320600	0191506010485	EMBOLD FIBERED 32X60	32	60	1	N/A	\$2,050

EMBOLD™ Soft Detachable Coil System

UPN / Order Code	GTIN	Description	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001394250020020	00191506026363	EMBOLD SOFT 2X2	2	2	1	N/A	\$1,460
M001394250020040	00191506026370	EMBOLD SOFT 2X4	2	4	1	N/A	\$1,460
M001394250020060	00191506026387	EMBOLD SOFT 2X6	2	6	1	N/A	\$1,460
M001394250020080	00191506026394	EMBOLD SOFT 2X8	2	8	1	N/A	\$1,460
M001394250030040	00191506026400	EMBOLD SOFT 3X4	3	4	1	N/A	\$1,460
M001394250030080	00191506026417	EMBOLD SOFT 3X8	3	8	1	N/A	\$1,460
M001394250030120	00191506026424	EMBOLD SOFT 3X12	3	12	1	N/A	\$1,620
M001394250030150	00191506026431	EMBOLD SOFT 3X15	3	15	1	N/A	\$1,620
M001394250040080	00191506026448	EMBOLD SOFT 4X8	4	8	1	N/A	\$1,460
M001394250040100	00191506026455	EMBOLD SOFT 4X10	4	10	1	N/A	\$1,620
M001394250040150	00191506026462	EMBOLD SOFT 4X15	4	15	1	N/A	\$1,620
M001394250040300	00191506026479	EMBOLD SOFT 4X30	4	30	1	N/A	\$2,050
M001394250050100	00191506026486	EMBOLD SOFT 5X10	5	10	1	N/A	\$1,620
M001394250050150	00191506026493	EMBOLD SOFT 5X15	5	15	1	N/A	\$1,620
M001394250050200	00191506026509	EMBOLD SOFT 5X20	5	20	1	N/A	\$1,620
M001394250050300	00191506026516	EMBOLD SOFT 5X30	5	30	1	N/A	\$2,050
M001394250060150	00191506026523	EMBOLD SOFT 6X15	6	15	1	N/A	\$1,620
M001394250060200	00191506026530	EMBOLD SOFT 6X20	6	20	1	N/A	\$1,620
M001394250060300	00191506026547	EMBOLD SOFT 6X30	6	30	1	N/A	\$2,050
M001394250060600	00191506026554	EMBOLD SOFT 6X60	6	60	1	N/A	\$2,050
M001394250080300	00191506026561	EMBOLD SOFT 8X30	8	30	1	N/A	\$2,050
M001394250080600	00191506026578	EMBOLD SOFT 8X60	8	60	1	N/A	\$2,050

EMBOLD™ Packing Detachable Coil System

UPN / Order Code	GTIN	Description	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001394260000150	00191506026585	EMBOLD PACKING US 15	N/A	15	1	N/A	\$1,775
M001394260000300	00191506026592	EMBOLD PACKING US 30	N/A	30	1	N/A	\$2,205
M001394260000450	00191506026608	EMBOLD PACKING US 45	N/A	45	1	N/A	\$2,205
M001394260000600	00191506026615	EMBOLD PACKING US 60	N/A	60	1	N/A	\$2,205

Detachable Coils

Interlock™ - 35 Fibered IDC Occlusion System

2D Helical

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001363500	36-350	08714729793151	3	4	1	N/A	\$870
M001363520	36-352	08714729792987	4	10	1	N/A	\$870
M001363540	36-354	08714729793007	6	10	1	N/A	\$870
M001363550	36-355	08714729793021	6	20	1	N/A	\$1,225
M001363570	36-357	08714729797425	8	10	1	N/A	\$870
M001363580	36-358	08714729793045	8	20	1	N/A	\$1,225
M001363590	36-359	08714729793014	8	40	1	N/A	\$1,595
M001363600	36-360	08714729793137	10	20	1	N/A	\$1,225
M001363610	36-361	08714729793168	10	40	1	N/A	\$1,595
M001363620	36-362	08714729793175	12	20	1	N/A	\$1,225
M001363630	36-363	08714729792970	12	40	1	N/A	\$1,595
M001363640	36-364	08714729793120	15	20	1	N/A	\$1,225
M001363650	36-365	08714729793144	15	40	1	N/A	\$1,595
M001363660	36-366	08714729797432	18	20	1	N/A	\$1,225
M001363670	36-367	08714729795414	18	40	1	N/A	\$1,595

Cube

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001363700	36-370	08714729792994	4	6	1	N/A	\$870
M001363720	36-372	08714729793069	6	10	1	N/A	\$870
M001363730	36-373	08714729793090	6	20	1	N/A	\$1,225
M001363760	36-376	08714729793052	8	20	1	N/A	\$1,225
M001363790	36-379	08714729793106	10	25	1	N/A	\$1,225
M001363800	36-380	08714729795513	10	40	1	N/A	\$1,595
M001363810	36-381	08714729793113	15	25	1	N/A	\$1,225
M001363820	36-382	08714729795438	15	40	1	N/A	\$1,595
M001363830	36-383	08714729795445	20	40	1	N/A	\$1,595

Diamond

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001363910	36-391	08714729795483	4	4.5	1	N/A	\$870
M001363920	36-392	08714729795520	6	9	1	N/A	\$870
M001363930	36-393	08714729796114	8	14	1	N/A	\$870

Interlock™ - 18 Fibered IDC Occlusion System

2D

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001361480	36-148	08714729845263	2	4	1	N/A	\$795
M001361490	36-149	08714729845270	2	6	1	N/A	\$795
M001361500	36-150	08714729764939	3	6	1	N/A	\$795
M001361510	36-151	08714729764946	3	12	1	N/A	\$795
M001361520	36-152	08714729764953	4	8	1	N/A	\$795
M001361530	36-153	08714729764960	4	15	1	N/A	\$795
M001361540	36-154	08714729764977	5	8	1	N/A	\$795
M001361550	36-155	08714729764984	5	15	1	N/A	\$795
M001361560	36-156	08714729764991	6	10	1	N/A	\$795
M001361570	36-157	08714729765004	6	20	1	N/A	\$860
M001361580	36-158	08714729765011	8	20	1	N/A	\$860
M001361590	36-159	08714729765028	10	20	1	N/A	\$860
M001361600	36-160	08714729765035	10	30	1	N/A	\$860
M001361920	36-192	08714729845287	10	50	1	N/A	\$1,180
M001361610	36-161	08714729765042	12	20	1	N/A	\$860
M001361620	36-162	08714729765059	12	30	1	N/A	\$860
M001361630	36-163	08714729765066	14	20	1	N/A	\$860
M001361640	36-164	08714729765073	14	30	1	N/A	\$860
M001361930	36-193	08714729845294	14	50	1	N/A	\$1,180
M001361940	36-194	08714729845300	18	50	1	N/A	\$1,180
M001361950	36-195	08714729845317	20	50	1	N/A	\$1,180
M001361960	36-196	08714729845324	22	60	1	N/A	\$1,180

VortX™ Diamond

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M001361740	36-174	08714729765080	3	2.3	1	N/A	\$795
M001361750	36-175	08714729765097	4	4.1	1	N/A	\$795
M001361760	36-176	08714729765103	5	5.8	1	N/A	\$795
M001361770	36-177	08714729765110	6	8	1	N/A	\$795

IDC™ Interlocking Detachable Coil

Soft

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Restrained Length (cm)	Box (Units)	C-Code	Price
M0033612020	361202	08714729158684	2	2	1	N/A	\$900
M0013612040	361204	08714729158691	2	4	1	N/A	\$900
M0013613060	361306	08714729158714	3	6	1	N/A	\$900
M0013613100	361310	08714729158721	3	10	1	N/A	\$900
M0013614040	361404	08714729158738	4	4	1	N/A	\$900
M0013614080	361408	08714729158745	4	8	1	N/A	\$900
M0013614120	361412	08714729158752	4	12	1	N/A	\$900
M0013615080	361508	08714729158769	5	8	1	N/A	\$900
M0013615120	361512	08714729158776	5	12	1	N/A	\$900

Pushable Coils

VortX™ Diamond - 18

Fibered Platinum Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013822031	382203	08714729158912	3	3.3	23	5	N/A	\$920
M0013822041	382204	08714729158929	4	3.7	41	5	N/A	\$920
M0013822051	382205	08714729158936	5	5.5	58	5	N/A	\$920
M0013822061	382206	08714729158943	6	6.7	80	5	N/A	\$920

VortX™ - 18 Coil

Fibred Platinum Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013812031	381203	08714729158875	3	2.5	22	5	N/A	\$920
M0013812041	381204	08714729158882	4	4	42	5	N/A	\$920
M0013812051	381205	08714729158899	5	5.5	60	5	N/A	\$920
M0013812061	381206	08714729158905	6	6.5	85	5	N/A	\$920

Straight-18 Coil

Fibred Platinum Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013120021	312002	08714729155980	Straight	2	2	5	N/A	\$700
M0013120051	312005	08714729155997	Straight	5	5	5	N/A	\$700

Figure 8 - 18 Coil

Fibred Platinum Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013120211	312021	08714729156000	2	5	10	5	N/A	\$700

Multi-Loop - 18 Coil

Fibred Platinum Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013120431	312043	08714729156048	4	7	30	5	N/A	\$700

Complex Helical - 18 Coil

Fibred Platinum Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013120221	312022	08714729156017	4	4	20	5	N/A	\$700
M0013120331	312033	08714729156031	6	6	30	5	N/A	\$700
M0013120441	312044	08714729156055	7	10	40	5	N/A	\$700
M0013120551	312055	08714729156062	8	12	50	5	N/A	\$700
M0013120661	312066	08714729156079	10	14	60	5	N/A	\$700
M0013120771	312077	08714729156086	11	17	70	5	N/A	\$700

Coil Pusher - 16 Coil

Coil Pusher-16

UPN / Order Code	Catalog No.	GTIN	Total Length (cm)	Outer Diameter (in)	Minimum Catheter Lumen Diameter (in)	Box (Units)	C-Code	Price
M0014012160	401216	08714729158967	177	0.016	0.018	1	N/A	\$265

2D Helical - 35

Fibred Platinum Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013723021	372302	08714729192701	3	2.6	20	5	N/A	\$720
M0013723041	372304	08714729158806	3	5.2	40	5	N/A	\$720
M0013724031	372403	08714729158813	4	2.9	30	5	N/A	\$720
M0013725031	372503	08714729192718	5	2.4	30	5	N/A	\$720
M0013725051	372505	08714729158837	5	4	50	5	N/A	\$720
M0013726041	372604	08714729158844	6	2.6	40	5	N/A	\$720
M0013727041	372704	08714729158851	7	2.3	40	5	N/A	\$720
M0013729061	372906	08714729158868	9	2.7	60	5	N/A	\$720

VortX™ - 35 Coil
Vascular Occlusion Coils

UPN / Order Code	Catalog No.	GTIN	Maximum Diameter (mm)	Unrestrained Coil Length (mm)	Length in Introducer (mm)	Box (Units)	C-Code	Price
M0013732041	373204	08714729162728	4	4	30	5	N/A	\$730
M0013733051	373305	08714729162735	5	4.5	35	5	N/A	\$730
M0013733061	373306	08714729162742	6	5	53	5	N/A	\$730
M0013733071	373307	08714729162759	7	5.5	67	5	N/A	\$730

BSC Part Number	GTIN	Description	C-Code
M00139552US0300	05060116921373	Y-90 TheraSphere (3 GBq) US Comm	C2616
M00139552US0350	05060116921380	Y-90 TheraSphere (3.5 GBq) US Comm	C2616
M00139552US0400	05060116921397	Y-90 TheraSphere (4 GBq) US Comm	C2616
M00139552US0450	05060116921403	Y-90 TheraSphere (4.5 GBq) US Comm	C2616
M00139552US0500	05060116921410	Y-90 TheraSphere (5 GBq) US Comm	C2616
M00139552US0550	05060116921427	Y-90 TheraSphere (5.5 GBq) US Comm	C2616
M00139552US0600	05060116921434	Y-90 TheraSphere (6 GBq) US Comm	C2616
M00139552US0650	05060116921441	Y-90 TheraSphere (6.5 GBq) US Comm	C2616
M00139552US0700	05060116921458	Y-90 TheraSphere (7 GBq) US Comm	C2616
M00139552US0750	05060116921465	Y-90 TheraSphere (7.5 GBq) US Comm	C2616
M00139552US0800	05060116921472	Y-90 TheraSphere (8 GBq) US Comm	C2616
M00139552US0850	05060116921489	Y-90 TheraSphere (8.5 GBq) US Comm	C2616
M00139552US0900	05060116921496	Y-90 TheraSphere (9 GBq) US Comm	C2616
M00139552US0950	05060116921502	Y-90 TheraSphere (9.5 GBq) US Comm	C2616
M00139552US1000	05060116921519	Y-90 TheraSphere (10 GBq) US Comm	C2616
M00139552US1050	05060116921526	Y-90 TheraSphere (10.5 GBq) US Comm	C2616
M00139552US1100	05060116921533	Y-90 TheraSphere (11 GBq) US Comm	C2616
M00139552US1150	05060116921540	Y-90 TheraSphere (11.5 GBq) US Comm	C2616
M00139552US1200	05060116921557	Y-90 TheraSphere (12 GBq) US Comm	C2616
M00139552US1250	05060116921564	Y-90 TheraSphere (12.5 GBq) US Comm	C2616
M00139552US1300	05060116921571	Y-90 TheraSphere (13 GBq) US Comm	C2616
M00139552US1350	05060116921588	Y-90 TheraSphere (13.5 GBq) US Comm	C2616
M00139552US1400	05060116921595	Y-90 TheraSphere (14 GBq) US Comm	C2616
M00139552US1450	05060116921601	Y-90 TheraSphere (14.5 GBq) US Comm	C2616
M00139552US1500	05060116921618	Y-90 TheraSphere (15 GBq) US Comm	C2616
M00139552US1550	05060116921625	Y-90 TheraSphere (15.5 GBq) US Comm	C2616
M00139552US1600	05060116921632	Y-90 TheraSphere (16 GBq) US Comm	C2616
M00139552US1650	05060116921649	Y-90 TheraSphere (16.5 GBq) US Comm	C2616
M00139552US1700	05060116921656	Y-90 TheraSphere (17 GBq) US Comm	C2616
M00139552US1750	05060116921663	Y-90 TheraSphere (17.5 GBq) US Comm	C2616
M00139552US1800	05060116921670	Y-90 TheraSphere (18 GBq) US Comm	C2616
M00139552US1850	05060116921687	Y-90 TheraSphere (18.5 GBq) US Comm	C2616
M00139552US1900	05060116921694	Y-90 TheraSphere (19 GBq) US Comm	C2616
M00139552US1950	05060116921700	Y-90 TheraSphere (19.5 GBq) US Comm	C2616
M00139552US2000	05060116921717	Y-90 TheraSphere (20 GBq) US Comm	C2616

Customer Ordering Information

Customer Service – TheraSphere

In the US or Canada, for orders or any supply-related questions about TheraSphere please contact:

➤ Phone: 1.866.363.3330 (toll free)

➤ Fax: 1.800.268.5299 (toll free)

➤ Fax: 1.613.482.4843

➤ Email: theraspherecustomersupport@bsci.com

➤ Service Hours of Operation: Monday through Friday, 8:00 a.m. to 8:00 p.m. ET

Order Online (US Sales only): www.TheraSphereNow.com

Simplicit90Y™ Personalized Dosimetry

UPN / Order Code	Description
S90Y-1YR-AS-FXD	Annual Subscription fixed
S90Y-1YR-AS-FLT	Annual Subscription floating
S90Y-1YR-FXD-RNW	Annual Subscription fixed renewal
S90Y-1YR-FLT-RNW	Annual Subscription floating renewal
S90Y-2YR-FXD	2YR Subscription fixed
S90Y-2YR-FLT	2YR Subscription floating
S90Y-2YR-FXD-RNW	2YR Subscription fixed renewal
S90Y-2YR-FLT-RNW	2YR Subscription floating renewal
S90Y-3YR-FXD	3YR Subscription fixed
S90Y-3YR-FLT	3YR Subscription floating
S90Y-3YR-FXD-RNW	3YR Subscription fixed renewal
S90Y-3YR-FLT-RNW	3YR Subscription floating renewal
S90Y-1YR-AS-FXD-A	Annual Subscription fixed - Additional
S90Y-1YR-AS-FLT-A	Annual Subscription floating- Additional
S90Y-FXD-ADD-RNW	Fixed licence additional renewal
S90Y-FLT-ADD-RNW	Floating licence additional renewal
S90Y-PERP-FXD	Perpetual fixed
S90Y-PERP-FXD-A	Perpetual fixed - Additional
S90Y-PERP-FLT	Perpetual floating
S90Y-PERP-FLT-A	Perpetual floating - Additional
S90Y-UP-PERP-FXD	Upgrade-Perpetual fixed
S90Y-UP-PERP-FXD-A	Upgrade-Perpetual fixed - Additional
S90Y-UP-PERP-FLT	Upgrade-Perpetual floating
S90Y-UP-FLT-A	Upgrade-Perpetual floating - Additional
S90Y-SLA	Service Level Agreement
S90Y-RESEARCH	Research License
S90Y-TRIAL	Trial License

Simplicit90Y is developed by Mirada Medical Ltd. and used under license by Biocompatibles UK Ltd. Boston Scientific is the sales agent for Simplit90Y.

Customer Ordering Information

Customer Service – Simplit90Y
 For questions please contact our dedicated TheraSphere Solution Specialists:
 ➤ Email: SimplicitySupport@bsci.com

Express™ LD Biliary Over-the-Wire Premounted Stent System
75 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
H74938046520750	38046-52075	08714729484912	5	17	20	6	75	C1876	\$2,345
H74938046530750	38046-53075	08714729484929	5	27	30	6	75	C1876	\$2,470
H74938046540750	38046-54075	08714729484936	5	37	40	6	75	C1876	\$2,470
H74938046560750	38046-56075	08714729484943	5	57	60	6	75	C1876	\$2,595
H74938046620750	38046-62075	08714729392132	6	17	20	6	75	C1876	\$2,345
H74938046630750	38046-63075	08714729579656	6	27	30	6	75	C1876	\$2,470
H74938046640750	38046-64075	08714729392156	6	37	40	6	75	C1876	\$2,470
H74938046660750	38046-66075	08714729392170	6	57	60	6	75	C1876	\$2,595
H74938046720750	38046-72075	08714729392194	7	17	20	6	75	C1876	\$2,345
H74938046730750	38046-73075	08714729579663	7	27	30	6	75	C1876	\$2,470
H74938046740750	38046-74075	08714729392217	7	37	40	6	75	C1876	\$2,470
H74938046760750	38046-76075	08714729392231	7	57	60	6	75	C1876	\$2,595
H74938046820750	38046-82075	08714729392255	8	17	20	6	75	C1876	\$2,345
H74938046830750	38046-83075	08714729579670	8	27	30	6	75	C1876	\$2,470
H74938046840750	38046-84075	08714729392279	8	37	40	6	75	C1876	\$2,470
H74938046860750	38046-86075	08714729392293	8	57	60	7	75	C1876	\$2,595
H74938046920750	38046-92075	08714729392316	9	25	30	7	75	C1876	\$2,345
H74938046940750	38046-94075	08714729392330	9	37	40	7	75	C1876	\$2,470
H74938046960750	38046-96075	08714729392354	9	57	60	7	75	C1876	\$2,595
H74938046102070	38046-10207	08714729392071	10	25	30	7	75	C1876	\$2,345
H74938046104070	38046-10407	08714729392095	10	37	40	7	75	C1876	\$2,470
H74938046106070	38046-10607	08714729392118	10	57	60	7	75	C1876	\$2,595

135 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Balloon Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
H74938047520130	38047-52013	08714729484950	5	17	20	6	135	C1876	\$2,345
H74938047530130	38047-53013	08714729484967	5	27	30	6	135	C1876	\$2,470
H74938047540130	38047-54013	08714729484974	5	37	40	6	135	C1876	\$2,470
H74938047560130	38047-56013	08714729484981	5	57	60	6	135	C1876	\$2,595
H74938047620130	38047-62013	08714729392408	6	17	20	6	135	C1876	\$2,345
H74938047630130	38047-63013	08714729579687	6	27	30	6	135	C1876	\$2,470
H74938047640130	38047-64013	08714729392415	6	37	40	6	135	C1876	\$2,470
H74938047660130	38047-66013	08714729392422	6	57	60	6	135	C1876	\$2,595
H74938047720130	38047-72013	08714729392439	7	17	20	6	135	C1876	\$2,345
H74938047730130	38047-73013	08714729579694	7	27	30	6	135	C1876	\$2,470
H74938047740130	38047-74013	08714729392446	7	37	40	6	135	C1876	\$2,470
H74938047760130	38047-76013	08714729392453	7	57	60	6	135	C1876	\$2,595
H74938047820130	38047-82013	08714729392460	8	17	20	6	135	C1876	\$2,345
H74938047830130	38047-83013	08714729579700	8	27	30	6	135	C1876	\$2,470
H74938047840130	38047-84013	08714729392477	8	37	40	6	135	C1876	\$2,470
H74938047860130	38047-86013	08714729392484	8	57	60	7	135	C1876	\$2,595
H74938047920130	38047-92013	08714729392491	9	25	30	7	135	C1876	\$2,345
H74938047940130	38047-94013	08714729392507	9	37	40	7	135	C1876	\$2,470
H74938047960130	38047-96013	08714729392514	9	57	60	7	135	C1876	\$2,595
H74938047120130	38047-12013	08714729392378	10	25	30	7	135	C1876	\$2,345
H74938047140130	38047-14013	08714729392385	10	37	40	7	135	C1876	\$2,470
H74938047160130	38047-16013	08714729392392	10	57	60	7	135	C1876	\$2,595

Express™ SD Renal Monorail™ Premounted Stent System
90 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Sheath Compatibility (F)	Guide Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
H74937911415900	37911-41590	08714729390381	4	15	5	6	90	C1876	\$2,345
H74937912419900	37912-41990	08714729390466	4	19	5	6	90	C1876	\$2,345
H74937911515900	37911-51590	08714729390404	5	15	5	6	90	C1876	\$2,345
H74937912519900	37912-51990	08714729390480	5	19	5	6	90	C1876	\$2,345
H74937911614900	37911-61490	08714729484837	6	14	5	6	90	C1876	\$2,345
H74937912618900	37912-61890	08714729484851	6	18	5	6	90	C1876	\$2,345
H74937911715900	37911-71590	08714729390442	7	15	6	7	90	C1876	\$2,345
H74937912719900	37912-71990	08714729390527	7	19	6	7	90	C1876	\$2,345

Express™ SD Renal Monorail™ Premounted Stent System
150 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Sheath Compatibility (F)	Guide Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
H74937911415150	37911-41515	08714729390374	4	15	5	6	150	C1876	\$2,345
H74937912419150	37912-41915	08714729390459	4	19	5	6	150	C1876	\$2,345
H74937911515150	37911-51515	08714729390398	5	15	5	6	150	C1876	\$2,345
H74937912519150	37912-51915	08714729390473	5	19	5	6	150	C1876	\$2,345
H74937911614150	37911-61415	08714729484820	6	14	5	6	150	C1876	\$2,345
H74937912618150	37912-61815	08714729484844	6	18	5	6	150	C1876	\$2,345
H74937911715150	37911-71515	08714729390435	7	15	6	7	150	C1876	\$2,345
H74937912719150	37912-71915	08714729390510	7	19	6	7	150	C1876	\$2,345

WALLSTENT™ RP Endoprosthesis
75 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
M001711060	71-106	08714729406174	5	80	6	100	1	C1876	\$4,835
M001711080	71-108	08714729406198	6	20	6	100	1	C1876	\$4,835
M001711100	71-110	08714729406211	6	45	6	100	1	C1876	\$4,835
M001711240	71-124	08714729406358	8	20	6	100	1, 3	C1876	\$4,835
M001711260	71-126	08714729406372	8	40	6	100	1, 3	C1876	\$4,835
M001711280	71-128	08714729406396	8	60	6	100	1, 3	C1876	\$4,835
M001711300	71-130	08714729406419	8	80	6	100	1, 3	C1876	\$4,835

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous

Fully open dimensions: to determine implanted length, refer to Box Label

135 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
M001711030	71-103	08714729406143	5	40	6	160	1	C1876	\$4,835
M001711050	71-105	08714729406167	5	55	6	160	1	C1876	\$4,835
M001711070	71-107	08714729406181	5	80	6	160	1	C1876	\$4,835
M001711230	71-123	08714729406341	7	90	6	160	1	C1876	\$4,835
M001711290	71-129	08714729406402	8	60	6	160	1	C1876	\$4,835
M001711310	71-131	08714729406426	8	80	6	160	1	C1876	\$4,835
M001711350	71-135	08714729406464	10	42	7	160	1	C1876	\$4,835
M001711370	71-137	08714729406488	10	68	7	160	1	C1876	\$4,835
M001711390	71-139	08714729406501	10	94	7	160	1	C1876	\$4,835

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous

Fully open dimensions: to determine implanted length, refer to Box Label

WALLSTENT™ Endoprosthesis with Unistep™ Plus Delivery System
75 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
H965404500	40450	08714729204121	22	35	11	100	1	C1876	\$5,715
H965404510	40451	08714729204138	22	45	11	100	1	C1876	\$5,715
H965404520	40452	08714729204145	22	70	11	100	1	C1876	\$5,715
H965405100	40510	08714729204152	24	35	12	100	1	C1876	\$5,715
H965405110	40511	08714729204169	24	45	12	100	1	C1876	\$5,715
H965405120	40512	08714729204176	24	70	12	100	1	C1876	\$5,715

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous

Fully open dimensions: to determine implanted length, refer to Box Label

135 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
H965412000	41-200	08714729204398	12	20	9	160	1	C1876	\$4,835

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous

Fully open dimensions: to determine implanted length, refer to Box Label

WallFlex™ Biliary Transhepatic Stent System

Fully Covered Stents with Permalume™ Covering

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
M00574750	57475	08714729820536	8	60	9	75 cm	C1874	\$4,395
M00574760	57476	08714729820543	8	80	9	75 cm	C1874	\$4,395
M00574770	57477	08714729820505	10	40	9	75 cm	C1874	\$4,395
M00574780	57478	08714729820512	10	60	9	75 cm	C1874	\$4,395
M00574790	57479	08714729820529	10	80	9	75 cm	C1874	\$4,395

Partially Covered Stents with Permalume™ Covering

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
M00574700	57470	08714729799382	8	60	9	75 cm	C1874	\$4,205
M00574710	57471	08714729799399	8	80	9	75 cm	C1874	\$4,205
M00574720	57472	08714729799368	10	40	9	75 cm	C1874	\$4,205
M00574730	57473	08714729799344	10	60	9	75 cm	C1874	\$4,205
M00574740	57474	08714729799375	10	80	9	75 cm	C1874	\$4,205

Uncovered Stents

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
M00574620	57462	08714729820581	8	40	9	75 cm	C1876	\$3,125
M00574630	57463	08714729799412	8	60	9	75 cm	C1876	\$3,125
M00574640	57464	08714729799429	8	80	9	75 cm	C1876	\$3,125
M00574650	57465	08714729799436	8	100	9	75 cm	C1876	\$3,125
M00574660	57466	08714729820550	10	40	9	75 cm	C1876	\$3,125
M00574670	57467	08714729820567	10	60	9	75 cm	C1876	\$3,125
M00574680	57468	08714729820574	10	80	9	75 cm	C1876	\$3,125
M00574690	57469	08714729799351	10	100	9	75 cm	C1876	\$3,125

Carotid WALLSTENT™ Monorail™ Endoprosthesis

UPN / Order Code	Catalog No.	GTIN	Unconstrained Diameter* (mm)	Unconstrained Length* (mm)	Guide Sheath Compatibility (F)	Guide Catheter Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
M001719000	71-900	08714729781165	6	22	5	7	135	C1876	\$4,175
M001719010	71-901	08714729781189	8	21	5	7	135	C1876	\$4,175
M001719020	71-902	08714729781202	8	29	5	7	135	C1876	\$4,175
M001719030	71-903	08714729781226	8	36	5	7	135	C1876	\$4,175
M001719040	71-904	08714729781103	10	24	6	8	135	C1876	\$4,175
M001719050	71-905	08714729781127	10	31	6	8	135	C1876	\$4,175
M001719060	71-906	08714729781141	10	37	6	8	135	C1876	\$4,175

*Unconstrained diameter selected should be 1-2 mm larger than largest vessel to be stented.

Innova™ Vascular Self-Expanding Stent System

75 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
H74939293050270	39293-05027	08714729873532	5	20	6	105	1,2	C1876	\$2,975
H74939293054070	39293-05407	08714729873549	5	40	6	105	1,2	C1876	\$2,975
H74939293056070	39293-05607	08714729873556	5	60	6	105	1,2	C1876	\$2,975
H74939293058070	39293-05807	08714729873563	5	80	6	105	1,2	C1876	\$2,975
H74939293051070	39293-05107	08714729873570	5	100	6	105	1,2	C1876	\$2,975
H74939293051270	39293-05127	08714729873587	5	120	6	105	1,2	C1876	\$2,975
H74939293051570	39293-05157	08714729873594	5	150	6	105	1,2	C1876	\$3,570
H74939293052070	39293-05207	08714729873617	5	200	6	105	1,2	C1876	\$4,565
H74939293060270	39293-06027	08714729873624	6	20	6	105	1,2	C1876	\$2,975
H74939293064070	39293-06407	08714729873631	6	40	6	105	1,2	C1876	\$2,975
H74939293066070	39293-06607	08714729873648	6	60	6	105	1,2	C1876	\$2,975
H74939293068070	39293-06807	08714729873655	6	80	6	105	1,2	C1876	\$2,975
H74939293061070	39293-06107	08714729873662	6	100	6	105	1,2	C1876	\$2,975
H74939293061270	39293-06127	08714729873679	6	120	6	105	1,2	C1876	\$2,975
H74939293061570	39293-06157	08714729873686	6	150	6	105	1,2	C1876	\$3,570
H74939293062030	39293-06203	08714729873709	6	200	6	105	1,2	C1876	\$4,565
H74939293070270	39293-07027	08714729873716	7	20	6	105	1,2	C1876	\$2,975
H74939293074070	39293-07407	08714729873723	7	40	6	105	1,2	C1876	\$2,975
H74939293076070	39293-07607	08714729873730	7	60	6	105	1,2	C1876	\$2,975
H74939293078070	39293-07807	08714729873747	7	80	6	105	1,2	C1876	\$2,975
H74939293071070	39293-07107	08714729873754	7	100	6	105	1,2	C1876	\$2,975
H74939293071270	39293-07127	08714729873761	7	120	6	105	1,2	C1876	\$2,975
H74939293071570	39293-07157	08714729873778	7	150	6	105	1,2	C1876	\$3,570
H74939293072070	39293-07207	08714729873792	7	200	6	105	1,2	C1876	\$4,565
H74939293080270	39293-08027	08714729873808	8	20	6	105	1,2	C1876	\$2,975
H74939293084070	39293-08407	08714729873815	8	40	6	105	1,2	C1876	\$2,975
H74939293086070	39293-08607	08714729873822	8	60	6	105	1,2	C1876	\$2,975
H74939293088070	39293-08807	08714729873839	8	80	6	105	1,2	C1876	\$2,975
H74939293081070	39293-08107	08714729873846	8	100	6	105	1,2	C1876	\$2,975
H74939293081270	39293-08127	08714729873853	8	120	6	105	1,2	C1876	\$2,975
H74939293081570	39293-08157	08714729873860	8	150	6	105	1,2	C1876	\$3,570
H74939293082070	39293-08207	08714729873884	8	200	6	105	1,2	C1876	\$4,565

130 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
H74939293050230	39293-05023	08714729873891	5	20	6	160	1,2	C1876	\$2,975
H74939293054030	39293-05403	08714729873907	5	40	6	160	1,2	C1876	\$2,975
H74939293056030	39293-05603	08714729873914	5	60	6	160	1,2	C1876	\$2,975
H74939293058030	39293-05803	08714729873921	5	80	6	160	1,2	C1876	\$2,975
H74939293051030	39293-05103	08714729873938	5	100	6	160	1,2	C1876	\$2,975
H74939293051230	39293-05123	08714729873945	5	120	6	160	1,2	C1876	\$2,975
H74939293051530	39293-05153	08714729873952	5	150	6	160	1,2	C1876	\$3,570
H74939293052030	39293-05203	08714729873976	5	200	6	160	1,2	C1876	\$4,565
H74939293060230	39293-06023	08714729873983	6	20	6	160	1,2	C1876	\$2,975
H74939293064030	39293-06403	08714729873990	6	40	6	160	1,2	C1876	\$2,975
H74939293066030	39293-06603	08714729874003	6	60	6	160	1,2	C1876	\$2,975
H74939293068030	39293-06803	08714729874010	6	80	6	160	1,2	C1876	\$2,975

Metallic Stents

Innova™ Vascular Self-Expanding Stent System (cont.)
130 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
H74939293061030	39293-06103	08714729874027	6	100	6	160	1,2	C1876	\$2,975
H74939293061230	39293-06123	08714729874034	6	120	6	160	1,2	C1876	\$2,975
H74939293061530	39293-06153	08714729874041	6	150	6	160	1,2	C1876	\$3,570
H74939293062030	39293-062030	08714729874065	6	200	6	160	1,2	C1876	\$4,565
H74939293070230	39293-07023	08714729874072	7	20	6	160	1,2	C1876	\$2,975
H74939293074030	39293-07403	08714729874089	7	40	6	160	1,2	C1876	\$2,975
H74939293076030	39293-07603	08714729874096	7	60	6	160	1,2	C1876	\$2,975
H74939293078030	39293-07803	08714729874102	7	80	6	160	1,2	C1876	\$2,975
H74939293071030	39293-07103	08714729874119	7	100	6	160	1,2	C1876	\$2,975
H74939293071230	39293-07123	08714729874126	7	120	6	160	1,2	C1876	\$2,975
H74939293071530	39293-07153	08714729874133	7	150	6	160	1,2	C1876	\$3,570
H74939293072070	39293-072070	08714729874157	7	200	6	160	1,2	C1876	\$4,565
H74939293080230	39293-08023	08714729874164	8	20	6	160	1,2	C1876	\$2,975
H74939293084030	39293-08403	08714729874171	8	40	6	160	1,2	C1876	\$2,975
H74939293086030	39293-08603	08714729874188	8	60	6	160	1,2	C1876	\$2,975
H74939293088030	39293-08803	08714729874195	8	80	6	160	1,2	C1876	\$2,975
H74939293081030	39293-08103	08714729874201	8	100	6	160	1,2	C1876	\$2,975
H74939293081230	39293-08123	08714729874218	8	120	6	160	1,2	C1876	\$2,975
H74939293081530	39293-08153	08714729874225	8	150	6	160	1,2	C1876	\$3,570
H74939293082030	39293-082030	8714729874249	8	200	6	160	1,2	C1876	\$4,565

Indications: 1 = Superficial Femoral Artery, 2 = Proximal Popliteal Artery

Express™ LD Iliac Over-the-Wire Premounted Stent System
75 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Balloon Length (mm)	Minimum Sheath Size/I.D. (in)	Catheter Total Length (cm)	C-Code	Price
H74938046620750	38046-62075	08714729392132	6	17	20	6	75	C1876	\$2,345
H74938046630750	38046-63075	08714729579656	6	27	30	6	75	C1876	\$2,470
H74938046640750	38046-64075	08714729392156	6	37	40	6	75	C1876	\$2,470
H74938046660750	38046-66075	08714729392170	6	57	60	6	75	C1876	\$2,595
H74938046720750	38046-72075	08714729392194	7	17	20	6	75	C1876	\$2,345
H74938046730750	38046-73075	08714729579663	7	27	30	6	75	C1876	\$2,470
H74938046740750	38046-74075	08714729392217	7	37	40	6	75	C1876	\$2,470
H74938046760750	38046-76075	08714729392231	7	57	60	6	75	C1876	\$2,595
H74938046820750	38046-82075	08714729392255	8	17	20	6	75	C1876	\$2,345
H74938046830750	38046-83075	08714729579670	8	27	30	6	75	C1876	\$2,470
H74938046840750	38046-84075	08714729392279	8	37	40	6	75	C1876	\$2,470
H74938046860750	38046-86075	08714729392293	8	57	60	7	75	C1876	\$2,595
H74938046920750	38046-92075	08714729392316	9	25	30	7	75	C1876	\$2,345
H74938046940750	38046-94075	08714729392330	9	37	40	7	75	C1876	\$2,470
H74938046960750	38046-96075	08714729392354	9	57	60	7	75	C1876	\$2,595
H74938046102070	38046-10207	08714729392071	10	25	30	7	75	C1876	\$2,345
H74938046104070	38046-10407	08714729392095	10	37	40	7	75	C1876	\$2,470
H74938046106070	38046-10607	08714729392118	10	57	60	7	75	C1876	\$2,595

135 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Balloon Length (mm)	Minimum Sheath Size/I.D. (in)	Catheter Total Length (cm)	C-Code	Price
H74938047620130	38047-62013	08714729392408	6	17	20	6	135	C1876	\$2,345
H74938047630130	38047-63013	08714729579687	6	27	30	6	135	C1876	\$2,470
H74938047640130	38047-64013	08714729392415	6	37	40	6	135	C1876	\$2,470
H74938047660130	38047-66013	08714729392422	6	57	60	6	135	C1876	\$2,595
H74938047720130	38047-72013	08714729392439	7	17	20	6	135	C1876	\$2,345
H74938047730130	38047-73013	08714729579694	7	27	30	6	135	C1876	\$2,470
H74938047740130	38047-74013	08714729392446	7	37	40	6	135	C1876	\$2,470
H74938047760130	38047-76013	08714729392453	7	57	60	6	135	C1876	\$2,595
H74938047820130	38047-82013	08714729392460	8	17	20	6	135	C1876	\$2,345
H74938047830130	38047-83013	08714729579700	8	27	30	6	135	C1876	\$2,470
H74938047840130	38047-84013	08714729392477	8	37	40	6	135	C1876	\$2,470
H74938047860130	38047-86013	08714729392484	8	57	60	7	135	C1876	\$2,595
H74938047920130	38047-92013	08714729392491	9	25	30	7	135	C1876	\$2,345
H74938047940130	38047-94013	08714729392507	9	37	40	7	135	C1876	\$2,470
H74938047960130	38047-96013	08714729392514	9	57	60	7	135	C1876	\$2,595
H74938047120130	38047-12013	08714729392378	10	25	30	7	135	C1876	\$2,345
H74938047140130	38047-14013	08714729392385	10	37	40	7	135	C1876	\$2,470
H74938047160130	38047-16013	08714729392392	10	57	60	7	135	C1876	\$2,595

Express™ SD Renal Monorail™ Premounted Stent System

90 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Sheath Compatibility (F)	Guide Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
H74937911415900	37911-41590	08714729390381	4	15	5	6	90	C1876	\$2,345
H74937912419900	37912-41990	08714729390466	4	19	5	6	90	C1876	\$2,345
H74937911515900	37911-51590	08714729390404	5	15	5	6	90	C1876	\$2,345
H74937912519900	37912-51990	08714729390480	5	19	5	6	90	C1876	\$2,345
H74937911614900	37911-61490	08714729484837	6	14	5	6	90	C1876	\$2,345
H74937912618900	37912-61890	08714729484851	6	18	5	6	90	C1876	\$2,345
H74937911715900	37911-71590	08714729390442	7	15	6	7	90	C1876	\$2,345
H74937912719900	37912-71990	08714729390527	7	19	6	7	90	C1876	\$2,345

150 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Stent Diameter (mm)	Stent Length (mm)	Sheath Compatibility (F)	Guide Compatibility (F)	Catheter Total Length (cm)	C-Code	Price
H74937911415150	37911-41515	08714729390374	4	15	5	6	150	C1876	\$2,345
H74937912419150	37912-41915	08714729390459	4	19	5	6	150	C1876	\$2,345
H74937911515150	37911-51515	08714729390398	5	15	5	6	150	C1876	\$2,345
H74937912519150	37912-51915	08714729390473	5	19	5	6	150	C1876	\$2,345
H74937911614150	37911-61415	08714729484820	6	14	5	6	150	C1876	\$2,345
H74937912618150	37912-61815	08714729484844	6	18	5	6	150	C1876	\$2,345
H74937911715150	37911-71515	08714729390435	7	15	6	7	150	C1876	\$2,345
H74937912719150	37912-71915	08714729390510	7	19	6	7	150	C1876	\$2,345

WALLSTENT™ RP Endoprosthesis

75 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
M001712000	71-200	08714729345589	6	24	6	100	4	C1876	\$4,835
M001712020	71-202	08714729345602	6	36	6	100	4	C1876	\$4,835
M001712040	71-204	08714729345626	6	46	6	100	4	C1876	\$4,835
M001712060	71-206	08714729345640	6	59	6	100	4	C1876	\$4,835
M001712080	71-208	08714729338130	7	23	6	100	4	C1876	\$4,835
M001712100	71-210	08714729345671	7	34	6	100	4	C1876	\$4,835
M001712120	71-212	08714729345695	7	55	6	100	4	C1876	\$4,835
M001712140	71-214	08714729345701	7	67	6	100	4	C1876	\$4,835
M001712160	71-216	08714729345725	8	20	6	100	4	C1876	\$4,835
M001712180	71-218	08714729345749	8	38	6	100	4	C1876	\$4,835
M001712200	71-220	08714729320067	8	47	6	100	4	C1876	\$4,835
M001712220	71-222	08714729345770	8	66	6	100	4	C1876	\$4,835
M001712260	71-226	08714729345817	9	35	6	100	4	C1876	\$4,835
M001712280	71-228	08714729320074	9	52	6	100	4	C1876	\$4,835
M001712300	71-230	08714729345848	9	61	6	100	4	C1876	\$4,835
M001712320	71-232	08714729345862	10	20	6	100	4	C1876	\$4,835
M001712340	71-234	08714729345886	10	39	6	100	4	C1876	\$4,835
M001712360	71-236	08714729345909	10	49	6	100	4	C1876	\$4,835
M001712380	71-238	08714729345923	10	69	6	100	4	C1876	\$4,835

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous
Fully open dimensions: to determine implanted length, refer to Box Label

135 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
M001712070	71-207	08714729345657	6	59	6	160	4	C1876	\$4,835

WALLSTENT™ RP Endoprosthesis (cont.)
135 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
M001712090	71-209	08714729345664	7	23	6	160	4	C1876	\$4,835
M001712110	71-211	08714729345688	7	34	6	160	4	C1876	\$4,395
M001712130	71-213	08714729338147	7	55	6	160	4	C1876	\$4,395
M001712190	71-219	08714729345756	8	38	6	160	4	C1876	\$4,395
M001712210	71-221	08714729345763	8	47	6	160	4	C1876	\$4,395
M001712230	71-223	08714729345787	8	66	6	160	4	C1876	\$4,395
M001712270	71-227	08714729345824	9	35	6	160	4	C1876	\$4,395
M001712350	71-235	08714729345893	10	39	6	160	4	C1876	\$4,395
M001712370	71-237	08714729345916	10	49	6	160	4	C1876	\$4,395
M001712390	71-239	08714729345930	10	69	6	160	4	C1876	\$4,395

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous
 Fully open dimensions: to determine implanted length, refer to Box Label

Epic™ Self-Expanding Vascular Stent System
75 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Working Length (cm)	Indications	C-Code	Price
H74939200062070	39200-06207	08714729805250	6	20	6	75	4	C1876	\$2,625
H74939200063070	39200-06307	08714729805274	6	30	6	75	4	C1876	\$2,625
H74939200064070	39200-06407	08714729805298	6	40	6	75	4	C1876	\$2,625
H74939200066070	39200-06607	08714729805335	6	60	6	75	4	C1876	\$2,625
H74939200068070	39200-06807	08714729805373	6	80	6	75	4	C1876	\$2,625
H74939200061070	39200-06107	08714729805212	6	100	6	75	4	C1876	\$2,625
H74939200061270	39200-06127	08714729805236	6	120	6	75	4	C1876	\$2,875
H74939200072070	39200-07207	08714729805410	7	20	6	75	4	C1876	\$2,625
H74939200073070	39200-07307	08714729805434	7	30	6	75	4	C1876	\$2,625
H74939200074070	39200-07407	08714729805458	7	40	6	75	4	C1876	\$2,625
H74939200076070	39200-07607	08714729805496	7	60	6	75	4	C1876	\$2,625
H74939200078070	39200-07807	08714729805533	7	80	6	75	4	C1876	\$2,625
H74939200071070	39200-07107	08714729805557	7	100	6	75	4	C1876	\$2,625
H74939200071270	39200-07127	08714729805397	7	120	6	75	4	C1876	\$2,875
H74939200082070	39200-08207	08714729805618	8	20	6	75	4	C1876	\$2,625
H74939200083070	39200-08307	08714729805632	8	30	6	75	4	C1876	\$2,625
H74939200084070	39200-08407	08714729805656	8	40	6	75	4	C1876	\$2,625
H74939200086070	39200-08607	08714729805694	8	60	6	75	4	C1876	\$2,625
H74939200088070	39200-08807	08714729805731	8	80	6	75	4	C1876	\$2,625
H74939200081070	39200-08107	08714729805571	8	100	6	75	4	C1876	\$2,625
H74939200081270	39200-08127	08714729805595	8	120	6	75	4	C1876	\$2,875
H74939200092070	39200-09207	08714729805755	9	20	6	75	4	C1876	\$2,625
H74939200093070	39200-09307	08714729805779	9	30	6	75	4	C1876	\$2,625
H74939200094070	39200-09407	08714729805793	9	40	6	75	4	C1876	\$2,625
H74939200096070	39200-09607	08714729805830	9	60	6	75	4	C1876	\$2,625
H74939200098070	39200-09807	08714729805878	9	80	6	75	4	C1876	\$2,625
H74939200091070	39200-09107	08714729805892	9	100	6	75	4	C1876	\$2,625
H74939200102070	39200-10207	08714729804970	10	20	6	75	4	C1876	\$2,625
H74939200103070	39200-10307	08714729804994	10	30	6	75	4	C1876	\$2,625
H74939200104070	39200-10407	08714729805014	10	40	6	75	4	C1876	\$2,625
H74939200106070	39200-10607	08714729805052	10	60	6	75	4	C1876	\$2,625
H74939200108070	39200-10807	08714729805090	10	80	6	75	4	C1876	\$2,625
H74939200101070	39200-10107	08714729805113	10	100	6	75	4	C1876	\$2,625
H74939200123070	39200-12307	08714729805137	12	30	6	75	4	C1876	\$2,625
H74939200124070	39200-12407	08714729805151	12	40	6	75	4	C1876	\$2,625
H74939200126070	39200-12607	08714729805199	12	60	6	75	4	C1876	\$2,625

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous
 Fully open dimensions: to determine implanted length, refer to Box Label

Epic™ Self-Expanding Vascular Stent System (cont.)
120 cm Catheter Working Length

UPN / Order Code	Catalog No.	GTIN	Unconstrained Stent Diameter (mm)	Unconstrained Stent Length (mm)	Sheath Compatibility (F)	Catheter Working Length (cm)	Indications	C-Code	Price
H74939200062020	39200-06202	08714729805243	6	20	6	120	4	C1876	\$2,625
H74939200063020	39200-06302	08714729805267	6	30	6	120	4	C1876	\$2,625
H74939200064020	39200-06402	08714729805281	6	40	6	120	4	C1876	\$2,625
H74939200066020	39200-06602	08714729805328	6	60	6	120	4	C1876	\$2,625
H74939200068020	39200-06802	08714729805366	6	80	6	120	4	C1876	\$2,625
H74939200061020	39200-06102	08714729805205	6	100	6	120	4	C1876	\$2,875
H74939200061220	39200-06122	08714729805229	6	120	6	120	4	C1876	\$2,625
H74939200072020	39200-07202	08714729805403	7	20	6	120	4	C1876	\$2,625
H74939200073020	39200-07302	08714729805427	7	30	6	120	4	C1876	\$2,625
H74939200074020	39200-07402	08714729805441	7	40	6	120	4	C1876	\$2,625
H74939200076020	39200-07602	08714729805489	7	60	6	120	4	C1876	\$2,625
H74939200078020	39200-07802	08714729805526	7	80	6	120	4	C1876	\$2,625
H74939200071020	39200-07102	08714729805540	7	100	6	120	4	C1876	\$2,625
H74939200071220	39200-07122	08714729805380	7	120	6	120	4	C1876	\$2,875
H74939200082020	39200-08202	08714729805601	8	20	6	120	4	C1876	\$2,625
H74939200083020	39200-08302	08714729805625	8	30	6	120	4	C1876	\$2,625
H74939200084020	39200-08402	08714729805649	8	40	6	120	4	C1876	\$2,625
H74939200086020	39200-08602	08714729805687	8	60	6	120	4	C1876	\$2,625
H74939200088020	39200-08802	08714729805724	8	80	6	120	4	C1876	\$2,625
H74939200081020	39200-08102	08714729805564	8	100	6	120	4	C1876	\$2,625
H74939200081220	39200-08122	08714729805588	8	120	6	120	4	C1876	\$2,875
H74939200092020	39200-09202	08714729805748	9	20	6	120	4	C1876	\$2,625
H74939200093020	39200-09302	08714729805762	9	30	6	120	4	C1876	\$2,625
H74939200094020	39200-09402	08714729805786	9	40	6	120	4	C1876	\$2,625
H74939200096020	39200-09602	08714729805823	9	60	6	120	4	C1876	\$2,625
H74939200098020	39200-09802	08714729805861	9	80	6	120	4	C1876	\$2,625
H74939200091020	39200-09102	08714729805885	9	100	6	120	4	C1876	\$2,625
H74939200102020	39200-10202	08714729804963	10	20	6	120	4	C1876	\$2,625
H74939200103020	39200-10302	08714729804987	10	30	6	120	4	C1876	\$2,625
H74939200104020	39200-10402	08714729805007	10	40	6	120	4	C1876	\$2,625
H74939200106020	39200-10602	08714729805045	10	60	6	120	4	C1876	\$2,625
H74939200108020	39200-10802	08714729805083	10	80	6	120	4	C1876	\$2,625
H74939200101020	39200-10102	08714729805106	10	100	6	120	4	C1876	\$2,625
H74939200123020	39200-12302	08714729805120	12	30	6	120	4	C1876	\$2,625
H74939200124020	39200-12402	08714729805144	12	40	6	120	4	C1876	\$2,625
H74939200126020	39200-12602	08714729805182	12	60	6	120	4	C1876	\$2,625

Indications: 1 = Tracheobronchial, 2 = Transjugular Intrahepatic Portosystemic Shunt (TIPS), 3 = Transhepatic Biliary, 4 = Iliac, 5 = Venous
 Fully open dimensions: to determine implanted length, refer to Box Label

Metallic Stents

VENOUS WALLSTENT™

75cm Working Catheter Length

UPN / Order Code	GTIN	Unconstrained Diameter* (mm)	Unconstrained Length (mm)	Guide Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
H74912044102070	8714729979401	10	20	6	75	Venous	C1876	\$4,835
H74912044104270	8714729979418	10	42	7	75	Venous	C1876	\$4,835
H74912044106870	8714729979425	10	68	7	75	Venous	C1876	\$4,835
H74912044109470	8714729979432	10	94	7	75	Venous	C1876	\$4,835
H74912044122070	8714729979449	12	20	9	75	Venous	C1876	\$4,835
H74912044124070	8714729979456	12	40	9	75	Venous	C1876	\$4,835
H74912044126070	8714729979470	12	60	9	75	Venous	C1876	\$4,835
H74912044129070	8714729979494	12	90	9	75	Venous	C1876	\$4,835
H74912044142070	8714729979517	14	20	10	75	Venous	C1876	\$4,835
H74912044144070	8714729979524	14	40	10	75	Venous	C1876	\$4,835
H74912044146070	8714729979531	14	60	10	75	Venous	C1876	\$4,835
H74912044149070	8714729979548	14	90	10	75	Venous	C1876	\$4,835
H74912044162070	8714729979555	16	20	10	75	Venous	C1876	\$4,835
H74912044164070	8714729979562	16	40	10	75	Venous	C1876	\$4,835
H74912044166070	8714729979579	16	60	10	75	Venous	C1876	\$4,835
H74912044169070	8714729979586	16	90	10	75	Venous	C1876	\$4,835
H74912044184070	8714729979593	18	40	11	75	Venous	C1876	\$4,835
H74912044186070	8714729979609	18	60	11	75	Venous	C1876	\$4,835
H74912044189070	8714729979616	18	90	11	75	Venous	C1876	\$4,835
H74912044204070	8714729979623	20	40	11	75	Venous	C1876	\$4,835
H74912044205570	8714729979630	20	55	11	75	Venous	C1876	\$4,835
H74912044208070	8714729979647	20	80	11	75	Venous	C1876	\$4,835

135cm Working Catheter Length

UPN / Order Code	GTIN	Unconstrained Diameter* (mm)	Unconstrained Length (mm)	Guide Sheath Compatibility (F)	Catheter Total Length (cm)	Indications	C-Code	Price
H74912044124010	8714729979463	12	40	9	135	Venous	C1876	\$4,835
H74912044126010	8714729979487	12	60	9	135	Venous	C1876	\$4,835
H74912044129010	8714729979500	12	90	9	135	Venous	C1876	\$4,835

THROMBECTOMY

Thrombectomy System

AngioJet® Ultra Thrombectomy Console

UPN / Order Code	Catalog No.	GTIN	Equipment Class	Power Requirements	Processor Memory Backup	Enclosure Protection	Foot Switch Protection	Size (D x W x H) (in)	Weight (lbs)	C-Code	Price
105650-001	105650-001	08714729890584	IEC Class 1 Type CF Device	110/130 VAC 220/240 VAC 50/60 Hz 10 Amps Max	60 Seconds	IPX1	IPX8	23 x 15 x 52.5	150	C1757	\$49,995

Thrombectomy Catheters

SmartClaw® Thrombectomy Catheter

UPN / Order Code	Catalog No.	GTIN	Delivery Platform	Minimum Vessel Diameter (mm)	Catheter Working Length (cm)	Catheter Diameter (F)	GuideWire (in)	Introducer Sheath (F)	Basket Size	C-Code	Price
FD0660-01	FD0660-01	191506050016	OTW	6	145	4.95	0.035	8 minimum	20mm	C1757	\$5,995
FD0660-02	FD0660-02	191506050023	OTW	6	145	4.95	0.035	8 minimum	32mm	C1757	\$5,995

ZelanteDVT™ Thrombectomy Set Ultra

UPN / Order Code	Catalog No.	GTIN	Delivery Platform	Minimum Vessel Diameter (mm)	Catheter Length (cm)	Catheter Diameter (F)	GuideWire (in)	Introducer Sheath (F)	Power Pulse® Delivery Enabled	C-Code	Price
114610-001	114610-001	08714729904724	OTW	6	105	8	0.035	8	Yes	C1757	\$6,615

Solent® Proxi

UPN / Order Code	Catalog No.	GTIN	Delivery Platform	Minimum Vessel Diameter (mm)	Catheter Length (cm)	Catheter Diameter (F)	GuideWire (in)	Introducer Sheath (F)	Power Pulse® Delivery Enabled	C-Code	Price
109676-001	109676-001	08714729889595	OTW	3	90	6	0.035	6	Yes	C1757	\$4,405

Solent® Omni

UPN / Order Code	Catalog No.	GTIN	Delivery Platform	Minimum Vessel Diameter (mm)	Catheter Length (cm)	Catheter Diameter (F)	GuideWire (in)	Introducer Sheath (F)	Power Pulse® Delivery Enabled	C-Code	Price
109681-001	109681-001	08714729889663	OTW	3	120	6	0.035	6	Yes	C1757	\$4,405

Solent® Dista

UPN / Order Code	Catalog No.	GTIN	Delivery Platform	Minimum Vessel Diameter (mm)	Catheter Length (cm)	Catheter Diameter (F)	GuideWire (in)	Introducer Sheath (F)	Power Pulse® Delivery Enabled	C-Code	Price
111303-001	111303-001	08714729889731	OTW	1.5	145	4 / 3	0.014	4	Yes	C1757	\$4,520

AVX®

UPN / Order Code	Catalog No.	GTIN	Delivery Platform	Minimum Vessel Diameter (mm)	Catheter Length (cm)	Catheter Diameter (F)	GuideWire (in)	Introducer Sheath (F)	Power Pulse® Delivery Enabled	C-Code	Price
105039-001	105039-001	08714729889045	OTW	3	50	6	0.035	6	No	C1757	\$1,525

Accessories

Ultra Power Pulse® Delivery Kit

UPN / Order Code	Catalog No.	GTIN	C-Code	Price
104834-0021	104834-0021	0871472988970	NA	\$390

For use only with the AngioJet Ultra Console
Packaged 5 per box

EKOS™ Endovascular System Control Unit 4.0

UPN/Order Code	GTIN	Equipment Class	Power Requirements	Processor Memory Back up	Size	Weight	C-Code	Price
600-40500	858593006462	Class II (product code KRA)	250VAC, 50/60Hz	Yes	26.9 x 22.7 x 11.4 cm	12.5 lb	N/A	\$49,995

EKOS™ Endovascular System

UPN/Order Code	GTIN	Product Description	Infusion Catheter Diameter (F)	Catheter Length (cm)	Treatment Zone Length (cm)	Guidewire (in)	C-Code	Price
500-55106	858593006028	EKOS™ 106CM 6 CM	5.4	106	6	0.035	C1889, C1887	\$6,015
500-55112	858593006134	EKOS™ 106CM 12 CM	5.4	106	12	0.035	C1889, C1887	\$6,015
500-55118	858593006141	EKOS™ 106CM 18 CM	5.4	106	18	0.035	C1889, C1887	\$6,015
500-55124	858593006158	EKOS™ 106CM 24 CM	5.4	106	24	0.035	C1889, C1887	\$6,015
500-55130	858593006165	EKOS™ 106CM 30 CM	5.4	106	30	0.035	C1889, C1887	\$6,015
500-55140	858593006172	EKOS™ 106CM 40 CM	5.4	106	40	0.035	C1889, C1887	\$6,015
500-55150	858593006189	EKOS™ 106CM 50 CM	5.4	106	50	0.035	C1889, C1887	\$6,015
500-56112	858593006264	EKOS™ 135CM 12 CM	5.4	135	12	0.035	C1889, C1887	\$6,015
500-56130	858593006295	EKOS™ 135CM 30 CM	5.4	135	30	0.035	C1889, C1887	\$6,015
500-56140	858593006301	EKOS™ 135CM 40 CM	5.4	135	40	0.035	C1889, C1887	\$6,015
500-56150	858593006318	EKOS™ 135CM 50 CM	5.4	135	50	0.035	C1889, C1887	\$6,015

EKOS+™ Endovascular System

UPN/Order Code	GTIN	Product Description	Infusion Catheter Diameter (F)	Catheter Length (cm)	Treatment Zone Length (cm)	Guidewire (in)	C-Code	Price
H74939605106080	00191506015459	EKOS+™ 106CM 8 CM	7.8	106	8	0.035	C1889, C1887	\$6,785
H74939605106120	00191506015466	EKOS+™ 106CM 12 CM	7.8	106	12	0.035	C1889, C1887	\$6,785
H74939605106160	00191506015473	EKOS+™ 106CM 16 CM	7.8	106	16	0.035	C1889, C1887	\$6,785
H74939605106200	00191506015480	EKOS+™ 106CM 20 CM	7.8	106	20	0.035	C1889, C1887	\$6,785
H74939605135080	00191506015497	EKOS+™ 135CM 8 CM	7.8	135	8	0.035	C1889, C1887	\$6,785
H74939605135120	00191506015503	EKOS+™ 135CM 12 CM	7.8	135	12	0.035	C1889, C1887	\$6,785
H74939605135160	00191506015510	EKOS+™ 135CM 16 CM	7.8	135	16	0.035	C1889, C1887	\$6,785
H74939605135200	00191506015527	EKOS+™ 135CM 20 CM	7.8	135	20	0.035	C1889, C1887	\$6,785

Intravascular Ultrasound (IVUS) - Catheters

Peripheral Imaging Catheter

OPTICROSS™ 18 Peripheral Imaging Catheter

UPN / Order Code	Catalog No.	GTIN	French Size	Guide Wire Size (in)	Usable Length (cm)	Transducer Frequency (MHz)	C-Code	Price
H7493932800180	N/A	08714729904366	4	0.018	135	30	C1753	\$1,375
*These two C-Codes are billed when IVUS is used in conjunction with atherectomy or stenting in fem-pop revascularization procedures. For more information, please contact your PI HEMA Field Reimbursement representative.								

OPTICROSS™ 35 Peripheral Imaging Catheter

UPN / Order Code	Catalog No.	GTIN	French Size	Guide Wire Size (in)	Usable Length (cm)	Transducer Frequency (MHz)	C-Code	Price
H7493932800350	N/A	8714729984542	8	0.035	115	15	C1753	\$1,480

Intravascular Ultrasound (IVUS) - Hardware

Avvigo™+ Multi-Modality Guidance System (Integrated Version)

UPN / Order Code	Catalog No.	GTIN	Dimensions	Frequency of Operation	Power Requirements	Weight	C-Code	Price
H749 249312010	N/A	00191506033163	N/A	9 - 60 MHz	100 - 240 VAC, 50/60 Hz	N/A	N/A	\$120,000
*List price does not include installation costs - contact your Boston Scientific Representative for details								

Avvigo™+ Multi-Modality Guidance System (Cart Version)

UPN / Order Code	Catalog No.	GTIN	Dimensions	Frequency of Operation	Power Requirements	Weight	C-Code	Price
H749 2493120C0	N/A	00191506033125	68" Height 28" Width 28" Depth	9 - 60 MHz	100 - 240 VAC, 50/60 Hz	118 lbs.	N/A	\$120,000

Intravascular Ultrasound Accessories

Avvigo™+ Guidance System Accessories

UPN / Order Code	Catalog No.	GTIN	Item	C-Code	Price
H749MDU5PLUS0	MDU5PLUS	08714729842279	Motor Drive	N/A	\$18,288
H749A70200	A7020	08714729228684	Sterile Pullback Sled	N/A	\$40
H749MDU5PLUSBAG0	MDU5PLUSBAG	08714729842255	MDU5 PLUS Sterile Bag	N/A	\$15
H749393160100	39316010	08714729896838	Permanent Sled	N/A	\$1,599
H749393150100	39315010	08714729896845	Permanent Sled Sterile Bag	N/A	\$150

Varithena

UPN / Order Code	GTIN	Description	CPT Code
FP00019	00360635118014	Varithena Polidocanol Injectable Foam 1% (18 mL)	36465 / 36466
FP00024	n/a	Varithena Universal Admin Pack	N/A

Customer Ordering Information

Customer Care – Varithena

In the US for orders or any supply-related questions about Varithena please contact:

- Phone: 1.855.971.8346 (toll free)
- Email: VarithenaCustomerCare@bsci.com
- Service Hours of Operation: Monday through Friday, 8:30 a.m. to 6:00 p.m. ET

Order Online (US Sales only): www.VarithenaProfessional.com