

ORISE™ Tissue Retractor System

The ORISE Tissue Retractor System is a medical device that provides and facilitates tissue retraction. The system is comprised of the ORISE Tissue Retractor and the ORISE Instrument Guides.


Designed to:

- Stabilize intraluminal space
- Help visualize the dissection plane
- Enable tissue manipulation
- Provide tissue retraction and counter-traction
- Allow independent movement of scope and tissue graspers




CAGE LENGTH

8.0 cm ± 0.5 cm




CAGE WIDTH

≥4.5 cm



CAGE HEIGHT

6.0 cm +1.0 cm, -0.5 cm



LENGTH

60 ± 2 cm



2.6 cm



ORISE Tissue Retractor

Order Number	Anatomy	Reference	Diameter (mm)	Cage Length (cm)	Working Length (cm)	# of Channels	Scope Compatibility (mm)
M00519510	Left Colon	OTR-60	26	8	60	3	≤11.7mm outer diameter

ORISE Instrument Guides

Order Number	Anatomy	Reference	Working Length (cm)	Tip Angles (degrees)	Tip Color	Instrument Compatibility (mm)
M00519520	Left Colon	OIG-90B	81	90	Blue	2.8
M00519530	Left Colon	OIG-90G	81	90	Gray	2.8
M00519540	Left Colon	OIG-60B	81	60	Blue	2.8
M00519550	Left Colon	OIG-45B	81	45	Blue	2.8

All photographs taken by Boston Scientific.

All trademarks are the property of their respective owners.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

Rx Only.

Caution: The law restricts this device to sale by or on the order of a physician. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.